

Parts to fit A-Dec Utility Module and Filter Regulators

UTILITY MODULE
OEM Part #24-0352-00
Exploded View

*It's clear!
The RPI Filter Housing
is clear so you can see
the filter inside!*

WATER SIDE

AIR SIDE

WHITE VALVE BODY REPAIR KIT
RPI Part #ADK153

BLACK VALVE BODY REPAIR KIT
RPI Part #ADK173

FILTER REGULATORS

GOOD NEWS!
These parts also fit the Regulators:

Filter Kit
(RPI Part #ADK201)

Black Valve Body Repair Kit
(RPI Part #ADK173)

White Valve Body Repair Kit
(RPI Part #ADK153)



replacement parts industries, inc.

20338 Corisco Street • Chatsworth, California • 91311
Call (800) 221-9723 • Fax (818) 882-7028 • Email www.order@rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit A-Dec equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. 07-10