

Model: Cascade 1040; Decade 1010, 1011, 1020 & 1021; Performer I, II & III; Positioner 1030; Priority 1005



Introducing New Switches to fit A-dec Dental Chairs!



RPI Part #ADS211
New OEM Part #044.184.01
Old OEM Part #044.049.01


LIMIT SWITCH (MODIFIED)

- SPDT
- 15.1A @ 125/250 VAC
- .187" tabs

Fits: Various positions & chairs

Models:

Cascade 1040
Decade 1010, 1011, 1020 & 1021
Performer III
Positioner 1030
Priority 1005



RPI Part #ADS212
New OEM Part #044.184.00
Old OEM Part #044.049.00


MICRO-SWITCH

- SPDT
- 15.1A @ 125/250 VAC
- .187" Tabs

Fits: Various positions and chairs

Models:

Decade 1010, 1011, 1020 & 1021
Performer I, II & III



RPI Part #ADS223
OEM Part #041.119.00


SWITCH (MINIATURE)

- 10.1A @ 125/250VAC

Fits: 8 Function Footswitch (ADEC part number 61.3043.80) 100VAC/120VAC & 240 VAC

Models:

Priority 1005 (Serial # M1-0001 to Present)



RPI Part #ADS224
OEM Part #041.059.00


SWITCH (PUSHBUTTON)

- 5A @ 120VAC, 2.5A @ 220VAC
- SPST; Momentary (N.O.)
- Black Button
- Panel mount with mounting nut

Fits: 8 Function Footswitch, Multi-Function Footswitch, Auto Return and Pre-position Assemblies

Models:

Cascade 1040
Decade 1021
Priority 1005 (Serial # D3-0135 to present)



RPI also carries Cylinders & Hydraulic Fluid to fit A-dec Dental Chairs!

RPI Part #	OEM Part #	Description	Cascade	Decade				Priority	
			1040	1010	1011	1015	1020	1021	1005
ADC174	61-0860-00	TILT CYLINDER KIT							•
ADC175	61-0850-00	LIFT CYLINDER KIT							•
ADC176	61-1287-00	LIFT CYLINDER KIT	•	•	•	•	•	•	
ADC177	61-2050-00	TILT CYLINDER	•						
ADC178	61-1267-00	TILT CYLINDER		•	•	•	•	•	
RPF449	No OEM # available	HYDRAULIC FLUID - 16 oz.	•	•	•	•	•	•	•
RPF453	No OEM # available	HYDRAULIC FLUID, CASE - 12 (16 oz.)	•	•	•	•	•	•	•

Watch for more parts to fit A-dec Dental Chairs!

Don't forget all of our other parts to fit A-dec Century II and Century Plus Dental Delivery Systems including the Control Block Service Kit, Handpiece Control Block Kit, Air Valve Cover Kit, Water Cover Valve Kit, Master Control Block Kit, 3 Block Tie Bolt Kit, Black and Gray Knob Sets, and more!



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit A-dec equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Website Fax 09/10