

SAFETY VALVE & KITS TO FIT TUTTNAUER E, M, EK & MK




RPI Kits now include both the ...
ELBOW FITTING and THREADED ADAPTER
For easy plumbing!

Now Take Your Pick ...
37 PSI or 40 PSI?
RPI Offers Both!

THE RPI ADVANTAGE!
Safety Valve & Kits


- Safety Valve now available for 37 PSI and 40 PSI!
- Updated Safety Valve Kits include both Elbow Fitting and Threaded Adapter for easy plumbing!



DESCRIPTION & PART #	FITS E & M SERIES		FITS EK & MK SERIES	
	37 PSI	40 PSI	37 PSI	40 PSI
AIR JET VALVE	TUJ034 (BLACK TOP)	TUJ034 (BLACK TOP)	TUJ033 (RED TOP)	TUJ033 (RED TOP)
SAFETY VALVE HOLDER KIT <i>Includes Air Jet Valve, Safety Valve & Safety Valve Holder (see below)</i>	TUK054	TUK078	TUK053	TUK077
SAFETY VALVE HOLDER <i>Includes Safety Valve Holder, Threaded Adapter, Elbow Fitting & Hardware</i>	TUH032	TUH032	TUH031	TUH031
SAFETY VALVE	TUV011	TUV065	TUV011	TUV065
THREADED ADAPTER	TUA060	TUA060	TUA060	TUA060
ELBOW FITTING	TUF079	TUF079	TUF079	TUF079
MOUNTING HARDWARE	TUK055	TUK055	TUK055	TUK055




GASKET FITS 3850/3870 E/EA/M



RPI Part #TUG074
OEM Part #2610019
DOOR GASKET
• Red Silicone; 15.945" OD
Models: 3850/3870 E, EA & M Series



HEPA FILTER & COVER FITS EA, EKA & EZ



RPI Part #RPF644
OEM Part #03140036
HEPA FILTER
• ~ 2.5" lg.
Models: EA, EKA & EZ Series



RPI Part #TUC067
OEM Part #02550019
WATER RESERVOIR/FILTER COVER
• Can also be used as Reservoir Cover
Models: EA, EKA & EZ Series



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@riparts.com • WEBSITE www.riparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Tuttnauer equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. web 08/07

GREAT PRICES AND ALL PARTS ARE GUARANTEED TO MEET YOUR 100% COMPLETE SATISFACTION!



RPI PART #	DESCRIPTION
RPA369	THREADLOCKER 545
RPF644	HEPA FILTER
RPT579	1/4" TEFLON THREAD SEALING TAPE
RPT580	1/2" TEFLON THREAD SEALING TAPE
TUA060	THREADED ADAPTER
TUC067	WATER RESERVOIR/FILTER COVER
TUF079	ELBOW FITTING
TUG074	DOOR GASKET
TUH031	SAFETY VALVE HOLDER
TUH032	SAFETY VALVE HOLDER
TUJ033	AIR JET VALVE (RED TOP)
TUJ034	AIR JET VALVE (BLACK TOP)
TUK053	SAFETY VALVE HOLDER KIT (37PSI)
TUK054	SAFETY VALVE HOLDER KIT (37PSI)
TUK055	MOUNTING HARDWARE (SAFETY VALVE HOLDER)
TUK077	SAFETY VALVE HOLDER KIT (40PSI)
TUK078	SAFETY VALVE HOLDER KIT (40PSI)
TUV011	SAFETY VALVE (37PSI)
TUV065	SAFETY VALVE (40PSI)

SERVICE TIP!
 When servicing the Safety Valve & Holder, use Threadlocker or Sealing Tape for a secure fit.

Threadlocker 545
 RPI Part #RPA369

Teflon Thread Sealing Tape
 1/4" wide - RPI Part #RPT579


Teflon Thread Sealing Tape
 1/2" wide - RPI Part #RPT580

————— An ISO 9001 Certified Company —————

• Address Service Requested •

PERMITTED #2399
 VAN NUYS, CA
 U.S. POSTAGE PAID
 PRESORTED STANDARD

id!
 replacement parts industries, inc.
 P.O. Box 5019 • Chatsworth, California • 91313-5019



NEW & IMPROVED PARTS TO FIT TUTTNAUER AUTOCLAVES!

NEW

DOOR GASKET
 Fits 3850/3870 E, EA & M Series



IMPROVED

SAFETY VALVE HOLDER KITS
 Fits E, M, EK & MK Series

- Kits now include both the Threaded Adapter and the Elbow Fitting!
- Safety Valve now available for 37 PSI and 40 PSI!



NEW

HEPA FILTER & WATER RESERVOIR/FILTER COVER
 Fits EA, EKA & EZ Series

