




Models: 100, 200, 400/500, 700, 800, 3522, 3523, 3525, 3533, 3622, 3623, 3633, 4200MPS, 6900AMPS & MC3600

# New Parts to fit Castle/Getinge® Bulk Sterilizers

**CASTLE/GETINGE**




**RPI Part #CSK048**  
**OEM Part #68164**

**SOLENOID VALVE REPAIR KIT**

- Includes all parts as shown

**Fits:** Solenoid Valve  
(OEM Part #517338)

**Models:** 3522, 3525 & 3533




**RPI Part #CSK050**  
**OEM Part #54187**

**SOLENOID VALVE REPAIR KIT**

- Includes all parts as shown

**Fits:** Solenoid Valve  
(OEM Part #517343)

**Models:** 3525, 4200MPS & 6900AMPS




**RPI Part #CSK065**  
**OEM Part #513700**

**CHECK VALVE REBUILD KIT**

- For use with a 3/4" or 1" NPT valve
- Includes all parts as shown

**Fits:** 3/4" or 1" NPT In-Line Check Valve

**Models:** 3522, 3523, 3525, 3533, 3622 & 3633



**RPI Part #CSK066**  
**OEM Part #513369**

**CHECK VALVE REBUILD KIT**

- For use with a 1/2", 3/8" or 3/4" valve
- Includes all parts as shown

**Fits:** Check Valve (OEM Part #'s 513361, 513362 & 513363)

**Models:** 100, 200, 400/500, 700, 800, 3522, 3523, 3525, 3533, 3622, 3633, 4200MPS & 6900AMPS

\* Fits 1/2" Check Valve (OEM Part #513361)  
\*\* Fits 3/8" Check Valve (OEM Part #513362) and 3/4" Check Valve (OEM Part #513363)



**Look Inside! →**

**RPI has them ...**

**PM and Spare Parts Packages to fit  
Castle/Getinge 3500 & 3600 Series  
Bulk Sterilizers!**



**replacement parts industries, inc.**
















CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@rpiparts.com](mailto:order@rpiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Castle/Getinge equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Web 09/09

SECTION  
**10**  
PAGE  
**254**  
UPDATED  
09/09



## New PM and Spare Parts Packages

RPI PART #	OEM PART #	DESCRIPTION	RPI PART #'S INCLUDED IN KIT
 CSP036	61301602085	<b>PM PKG (SOLENOID VALVES)</b> Models: 3522, 3523 & 3533	AMK174, (2) AMK177, AMK183 & (2) CSK048
 CSP037	61301602090	<b>PM PKG (SOLENOID VALVES)</b> Models: 3525	AMK174, (2) AMK177, AMK183, CSD051, (2) CSK048 & CSK050
 CSP038	61301602091	<b>PM PKG (FILTER &amp; SCREEN)</b> Models: 3522, 3523, 3525 & 3533	CSF052 & (2) CSS053
 CSP039	61301602087	<b>PM PKG (STM TRAP &amp; CHK VLV)</b> Models: 3525	(3) CSK066 & (2) CST059
 CSP040	61301602092	<b>PM PKG (STM TRAP &amp; CHK VLV)</b> Models: 3525	CSK066, CSV056, CSV057 & CSV058
 CSP068	61301602079	<b>PM PKG (STM TRAP &amp; CHK VLV)</b> Models: 3622 & 3633	CSK060, CSK065, (3) CSK066 & CST059
 CSP069	No OEM Part # Available	<b>PM PKG (STM TRAP &amp; CHK VLV)</b> Models: 3522, 3523, 3525 & 3533	(3) CSK066 & CSK067
 CSP070	No OEM Part # Available	<b>PM PKG (STM TRAP &amp; CHK VLV)</b> Models: 3622 & 3633	CSK060, CSK065, (3) CSK066 & CSK067
 CSP041	61301602510	<b>SPARE PARTS PACKAGE</b> Models: 3525	AMK174, (2) AMK177, AMK183, CSD051, CSF052, CSG002, (3) CSK048, CSK050, (3) CSK066, CSR028, CSR029, CSS054, CSV056, CSV057 & CSV058
 CSP042	61301602498	<b>SPARE PARTS PACKAGE</b> Models: 3522	AMK174, (2) AMK177, AMK183, CSG002, (2) CSK048, (3) CSK066, CSR028, CSR029 & (2) CST059
 CSP043	61301602499	<b>SPARE PARTS PACKAGE</b> Models: 3523	AMK174, (2) AMK177, AMK183, CSF052, CSG002, (3) CSK048, (3) CSK066, CSR028, CSR029 & (3) CST059
 CSP044	61301602500	<b>SPARE PARTS PACKAGE</b> Models: 3533	AMK174, (2) AMK177, AMK183, CSF052, CSG002, CSK018, (2) CSK048, (3) CSK066, CSR028, CSR029 & (2) CST059
 CSP045	61301602501	<b>SPARE PARTS PACKAGE</b> Models: 3522	AMK174, (2) AMK177, AMK183, CSG001, (2) CSK048, (3) CSK066, CSR028, CSR029 & (2) CST059
 CSP046	61301602502	<b>SPARE PARTS PACKAGE</b> Models: 3523	AMK174, (2) AMK177, AMK183, CSF052, CSG001, (3) CSK048, (3) CSK066, CSR028, CSR029 & (3) CST059
 CSP047	61301602503	<b>SPARE PARTS PACKAGE</b> Models: 3533	AMK174, (2) AMK177, AMK183, CSF052, CSG001, CSK018, (2) CSK048, (3) CSK066, CSR028, CSR029 & (2) CST059

RPI Part #AMK174  
OEM Part #56161

### SOLENOID VALVE REPAIR KIT


• Includes all parts as shown

**Included in PM Pkg:**

RPI Part #'s CSP036 & CSP037

**Included in Spare Parts Pkg:**

RPI Part #'s CSP041, CSP042, CSP043, CSP044, CSP045, CSP046 & CSP047



RPI Part #AMK177  
OEM Part #61301601739

### SOLENOID VALVE REPAIR KIT


• Includes all parts as shown

**Included in PM Pkg:**

RPI Part #'s CSP036 & CSP037

**Included in Spare Parts Pkg:**

RPI Part #'s CSP041, CSP042, CSP043, CSP044, CSP045, CSP046 & CSP047



RPI Part #AMK183  
OEM Part #54190

### SOLENOID VALVE REPAIR KIT


• Includes all parts as shown

**Included in PM Pkg:**

RPI Part #'s CSP036 & CSP037

**Included in Spare Parts Pkg:**

RPI Part #'s CSP041, CSP042, CSP043, CSP044, CSP045, CSP046 & CSP047



RPI Part #CSD051  
OEM Part #517153


### DIAPHRAGM REPLACEMENT 1.0

**Included in PM Pkg:**

RPI Part #'s CSP037

**Included in Spare Parts Pkg:**

RPI Part #'s CSP041



RPI Part #CSF052  
OEM Part #61301600303


### FILTER, AIR-IN

**Included in PM Pkg:**

RPI Part #'s CSP038

**Included in Spare Parts Pkg:**

RPI Part #'s CSP041, CSP043, CSP044, CSP046 & CSP047



**RPI Part #CSG001**  
**OEM Part #61301605016**  
 Old OEM Part #58398




**DOOR GASKET (20" x 20")**  
**Included in Spare Parts Pkg:**  
 RPI Part #'s CSP045, CSP046  
 & CSP047

**RPI Part #CSG002**  
**OEM Part #61301605017**  
 Old OEM Part #58397




**DOOR GASKET (16" x 16")**  
**Included in Spare Parts Pkg:**  
 RPI Part #'s CSP041, CSP042,  
 CSP043 & CSP044

**RPI Part #CSK018**  
**OEM Part #517295**




**SOLENOID VALVE REPAIR KIT**  
 • Includes all parts as shown  
**Included in Spare Parts Pkg:**  
 RPI Part #'s CSP044 & CSP047

**RPI Part #CSK048**  
**OEM Part #68164**




**SOLENOID VALVE REPAIR KIT**  
 • Includes all parts as shown  
**Included in PM Pkg:**  
 RPI Part #'s CSP036 & CSP037  
**Included in Spare Parts Pkg:**  
 RPI Part #'s CSP041, CSP042,  
 CSP043, CSP044, CSP045,  
 CSP046 & CSP047

**RPI Part #CSK050**  
**OEM Part #54187**



**SOLENOID VALVE REPAIR KIT**  
 • Includes all parts as shown  
**Included in PM Pkg:**  
 RPI Part #'s CSP037  
**Included in Spare Parts Pkg:**  
 RPI Part #'s CSP041


**RPI Part #CSK060**  
**OEM Part #27557**



**STEAM TRAP REBUILD KIT**  
 • Includes all parts as shown  
**Included in PM Pkg:**  
 RPI Part #'s CSP068 & CSP070




**RPI Part #CSK065**  
**OEM Part #513700**



**CHECK VALVE REBUILD KIT**  
 • Includes all parts as shown  
**Included in PM Pkg:**  
 RPI Part #'s CSP068 & CSP070




**RPI Part #CSK066**  
**OEM Part #513369**



See Page 1 for more details!


**CHECK VALVE REBUILD KIT**  
 • Includes all parts as shown  
**Included in PM Pkg:**  
 RPI Part #'s CSP039, CSP040, CSP068,  
 CSP069 & CSP070  
**Included in Spare Parts Pkg:** RPI Part #'s  
 CSP041, CSP042, CSP043, CSP044,  
 CSP045, CSP046 & CSP047

**RPI Part #CSK067**  
**OEM Part # (No OEM Part # Available)**




**STEAM TRAP REBUILD KIT**  
 • Includes all parts as shown  
 • Also available, complete  
 Steam Trap (RPI Part #CST059)  
**Included in PM Pkg:**  
 RPI Part #'s CSP069 & CSP070

**RPI Part #CSR028**  
**OEM Part #86407**




**SOLID STATE RELAY**  
**Included in Spare Parts Pkg:**  
 RPI Part #'s CSP041, CSP042,  
 CSP043, CSP044, CSP045,  
 CSP046 & CSP047

**RPI Part #CSR029**  
**OEM Part #86408**



**SOLID STATE RELAY**  
**Included in Spare Parts Pkg:**  
 RPI Part #'s CSP041, CSP042,  
 CSP043, CSP044, CSP045,  
 CSP046 & CSP047

**RPI Part #CSS053**  
**OEM Part #43907**




**SCREEN FILTER**  
**Included in PM Pkg:**  
 RPI Part #'s CSP038

**RPI Part #CSS054**  
**OEM Part #45690**



**SCREEN, DETERGENT**  
**Included in Spare Parts Pkg:**  
 RPI Part #'s CSP041

**RPI Part #CST059**  
**OEM Part #61301602492**




**STEAM TRAP**  
 • Includes all parts as shown  
**Included in PM Pkg:**  
 RPI Part #'s CSP039 & CSP068  
**Included in Spare Parts Pkg:**  
 RPI Part #'s CSP042, CSP043,  
 CSP044, CSP045, CSP046 &  
 CSP047

**RPI Part #CSV056**  
**OEM Part #14010**




**CHECK VALVE – 1/4"**  
**Included in PM Pkg:**  
 RPI Part #'s CSP040  
**Included in Spare Parts Pkg:**  
 RPI Part #'s CSP041

**RPI Part #CSV057**  
**OEM Part #19114**



**CHECK VALVE – 1/2"**  
**Included in PM Pkg:**  
 RPI Part #'s CSP040  
**Included in Spare Parts Pkg:**  
 RPI Part #'s CSP041

**RPI Part #CSV058**  
**OEM Part #08045**



**CHECK VALVE – 3/4"**  
**Included in PM Pkg:**  
 RPI Part #'s CSP040  
**Included in Spare Parts Pkg:**  
 RPI Part #'s CSP041



## LIST OF RPI PARTS TO FIT CASTLE/GETINGE® BULK STERILIZERS

RPI PART #	OEM PART #	DESCRIPTION	MODELS
<b>AMK174</b>	56161	Solenoid Valve Repair Kit	3522, 3523, 3525, 3533
<b>AMK177</b>	61301601739	Solenoid Valve Repair Kit	3522, 3523, 3525, 3533, 3622, 3623, 3633
<b>AMK183</b>	54190	Solenoid Valve Repair Kit	3522, 3523, 3525, 3533
CSC063	(No OEM Part # Available)	Coil	3522, 3533, 3622, 3623, 3633
CSC064	(No OEM Part # Available)	Coil	3523
<b>CSD051</b>	517153	Diaphragm Replacement 1.0	3522, 3523, 3525, 3533
<b>CSF052</b>	61301600303	Filter, Air-In	100, 200, 400/500, 3522, 3523, 3525, 3533, 3622, 3633, 4200MPS, 6900AMPS
<b>CSG001</b>	61301605016 58398 (old)	Door Gasket (20" X 20")	100, 800, 3522, 3523, 3533
<b>CSG002</b>	61301605017 58397 (old)	Door Gasket (16" X 16")	100, 3522, 3523, 3525, 3533
CSG003	40441	Door Gasket (Green)	6300, 6400, 6900
<b>CSK018</b>	517295	Solenoid Valve Repair Kit	3523, 3533
CSK019	61301601738 & 518018 (old)	Dual Solenoid Valve Repair Kit	3522, 3533, 3622, 3623, 3633
CSK020	517298	Solenoid Valve Repair Kit	3523
CSK022	517297	Solenoid Valve Repair Kit	3523
CSK023	61301600442	Solenoid Valve Repair Kit	3522, 3523, 3525, 3533, 4200MPS, 6900AMPS
CSK025	54186	Solenoid Valve Repair Kit	3622, 3623, 3633, 4200MPS, 6900AMPS
<b>CSK048</b>	68164	Solenoid Valve Repair Kit	3522, 3525, 3533
<b>CSK050</b>	54187	Solenoid Valve Repair Kit	3525, 4200MPS, 6900AMPS
<b>CSK060</b>	27557	Steam Trap Rebuild Kit	200, 3622, 3633, 4200MPS, 6900AMPS
CSK061	(No OEM Part # Available)	Solenoid Valve Repair Kit	3522, 3533, 3622, 3623, 3633
CSK062	(No OEM Part # Available)	Solenoid Valve Repair Kit	3523
<b>CSK065</b>	513700	Check Valve Rebuild Kit	3522, 3523, 3525, 3533, 3622, 3633
<b>CSK066</b>	513369	Check Valve Rebuild Kit	100, 200, 400/500, 700, 800, 3522, 3523, 3525, 3533, 3622, 3633, 4200MPS, 6900AMPS
<b>CSK067</b>	(No OEM Part # Available)	Steam Trap Rebuild Kit	100, 200, 400/500, 700, 800, 3522, 3523, 3525, 3533, 3622, 3633
<b>CSP036</b>	61301602085	PM Pkg (Solenoid Valves) <i>See page 2-3 for details.</i>	3522, 3523, 3533
<b>CSP037</b>	61301602090	PM Pkg (Solenoid Valves) <i>See page 2-3 for details.</i>	3525
<b>CSP038</b>	61301602091	PM Pkg (Filter & Screen) <i>See page 2-3 for details.</i>	3522, 3523, 3525, 3533
<b>CSP039</b>	61301602087	PM Pkg (Stm Trap & Chk Vlv) <i>See page 2-3 for details.</i>	3525
<b>CSP040</b>	61301602092	PM Pkg (Stm Trap & Chk Vlv) <i>See page 2-3 for details.</i>	3525



RPI PART #	OEM PART #	DESCRIPTION	MODELS
<b>CSP041</b>	61301602510	Spare Parts Package <i>See page 2-3 for details.</i>	3525
<b>CSP042</b>	61301602498	Spare Parts Package <i>See page 2-3 for details.</i>	3522
<b>CSP043</b>	61301602499	Spare Parts Package <i>See page 2-3 for details.</i>	3523
<b>CSP044</b>	61301602500	Spare Parts Package <i>See page 2-3 for details.</i>	3533
<b>CSP045</b>	61301602501	Spare Parts Package <i>See page 2-3 for details.</i>	3522
<b>CSP046</b>	61301602502	Spare Parts Package <i>See page 2-3 for details.</i>	3523
<b>CSP047</b>	61301602503	Spare Parts Package <i>See page 2-3 for details.</i>	3533
<b>CSP068</b>	61301602079	PM Pkg (Stm Trap & Chk Vlv) <i>See page 2-3 for details.</i>	3622, 3633
<b>CSP069</b>	(No OEM Part # Available)	PM Pkg (Stm Trap & Chk Vlv) <i>See page 2-3 for details.</i>	3522, 3523, 3525, 3533
<b>CSP070</b>	(No OEM Part # Available)	PM Pkg (Stm Trap & Chk Vlv) <i>See page 2-3 for details.</i>	3622, 3633
<b>CSR028</b>	86407	Solid State Relay	3522, 3523, 3525, 3533, 3622, 3623, 3633
<b>CSR029</b>	86408	Solid State Relay	3522, 3523, 3525, 3533, 3622, 3623, 3633
<b>CSS053</b>	43907	Screen Filter	3522, 3523, 3525, 3533, 3622, 3633
<b>CSS054</b>	45690	Screen, Detergent	3522, 3523, 3525, 3533, 3622, 3633
<b>CST059</b>	61301602492	Steam Trap	100, 200, 400/500, 700, 800, 3522, 3523, 3525, 3533, 3622, 3633
CSV004	517339	Solenoid Valve Assembly	3522, 3523, 3525, 3533, 3622, 3623, 3633
CSV005	517345	Solenoid Valve Assembly	3523
CSV007	517347	Solenoid Valve Assembly	3523, 3533
CSV008	517356	Dual Coil Solenoid Valve	3522, 3533, 3622, 3623, 3633
CSV009	517354	Solenoid Valve Assembly	3523
CSV010	517349	Solenoid Valve Assembly	3522, 3525
CSV011	517353	Solenoid Valve Assembly	3523
CSV012	541582	Solenoid Valve Assembly	3522, 3523, 3525, 3533, 4200MPS, 6900AMPS
CSV014	517342	Solenoid Valve Assembly	3622, 3623, 3633, 4200MPS, 6900AMPS
CSV033	56913	Safety Relief Valve	3622, 3633
CSV034	61301604180	Safety Relief Valve	3522, 3523, 3525, 3533
CSV035	532554	Relief Valve	3622, 3633
<b>CSV056</b>	14010	Check Valve	3522, 3523, 3525, 3533
<b>CSV057</b>	19114	Check Valve	3522, 3523, 3525, 3533
<b>CSV058</b>	08045	Check Valve	3522, 3523, 3525, 3533

RPI PM Packages and Spare Parts Packages shown in **BOLD** are composed of RPI Parts that are also sold separately. Please see pages 2-3 for more information. RPI Part #'s shown in **BOLD ITALICS** are also available in one or more PM Package or Spare Parts Package. Please see pages 2-3 for more information.