



# New Parts to fit CASTLE/GETINGE

**CASTLE/GETINGE**

## Other Solenoid Valve Assemblies & Kits Available

### SOL-1 STEAM TO CHAMBER

RPI Part #CSV004 OEM Part #517339  
 Models: 3522, 3523, 3525, 3533, 3622, 3623 & 3633  
 See also: Valve Repair Kit (RPI Part #AMK177)

### SOL-1 STEAM TO CHAMBER

RPI Part #CSV005 OEM Part #517345  
 Model: 3523  
 See also: Valve Repair Kit (RPI Part #AMK183)

### SOL-3 AIR-IN

RPI Part #CSV007 OEM Part #517347  
 Models: 3523 & 3533  
 See also: Valve Repair Kit (RPI Part #CSK018)

### SOL-4 MULTI BYPASS

RPI Part #CSV004 OEM Part #517339  
 Models: 3523, 3622, 3623 & 3633  
 See also: Valve Repair Kit (RPI Part #AMK177)

### SOL-5 SLOW EXHAUST

RPI Part #CSV009 OEM Part #517354  
 Model: 3523  
 See also: Valve Repair Kit (RPI Part #CSK020)

### SOL-5 & 6 SLOW EXHAUST & VAC-WATER

DUAL SOLENOID VALVE REPAIR KIT  
 RPI Part #CSK019 OEM Part #61301601738  
 Models: 3522, 3533, 3622, 3623 & 3633

### SOL-6 VAC-WATER

RPI Part #CSV010 OEM Part #517349  
 Models: 3522 & 3525  
 See also: Valve Repair Kit (RPI Part #AMK174)

### SOL-7 MULTI BYPASS SOLENOID VALVE REPAIR KIT

RPI Part #CSK022 OEM Part #517297  
 Model: 3523

### SOL-7 MULTI STEAM

RPI Part #CSV004 OEM Part #517339  
 Models: 3622, 3623 & 3633  
 See also: Valve Repair Kit (RPI Part #AMK177)

### SOL-8 JACKET


RPI Part #CSV004 OEM Part #517339  
 Models: 3522, 3523, 3533, 3622, 3623 & 3633  
 See also: Valve Repair Kit (RPI Part #AMK177)

### SOL-8 JACKET

RPI Part #CSV014 OEM Part #517342  
 Models: 3622, 3623, 3633, 4200 & 6900  
 See also: Valve Repair Kit (RPI Part #CSK025)

### SOL-10 & 11 DOOR

RPI Part #CSV012 OEM Part #541582  
 Models: 3522, 3523, 3525, 3533, 4200 & 6900  
 See also: Valve Repair Kit (RPI Part #CSK023)




### SOL-5 & 6 Slow Exhaust & Vac-Water

RPI Part #CSV008  
 OEM Part #517356

### DUAL COIL SOLENOID VALVE

- This is **not** an ASCO® Valve Assembly
- Includes all parts as shown

Models: 3522, 3533, 3622, 3623 & 3633



RPI Part #CSK061  
 OEM Part #


(No OEM Part # Available)

### SOLENOID VALVE REPAIR KIT

- Fits RPI Dual Coil Solenoid Valve (RPI Part #CSV008) **only**; For kit to repair OEM Dual Coil Solenoid Valve NO. 517356, see RPI Part #CSK019
- Includes all parts as shown

Models: 3522, 3533, 3622, 3623 & 3633

Also available, **DUAL SOLENOID VALVE REPAIR KIT** (RPI Part #CSK019) to fit Castle/Getinge valve #517356 **only**.




RPI Part #CSC063  
 OEM Part # (No OEM Part # Available)

### COIL

- Fits Dual Coil Solenoid Valve (RPI Part #CSV008) **only**

Models: 3522, 3533, 3622, 3623 & 3633




### SOL-7 Multitherm Steam

RPI Part #CSV011  
 OEM Part #517353

### SOLENOID VALVE ASSEMBLY

- This is **not** an ASCO® Valve Assembly
- Includes all parts as shown

Model: 3523 (Sol-7 Multitherm Steam)




RPI Part #CSK062  
 OEM Part #

(No OEM Part # Available)

### SOLENOID VALVE REPAIR KIT

- Fits Solenoid Valve Assembly (RPI Part #CSV011) **only**
- Includes all parts as shown

Model: 3523 (Sol-7 Multitherm Steam)



RPI Part #CSC064  
 OEM Part # (No OEM Part # Available)

### COIL

- Fits Solenoid Valve Assembly (RPI Part #CSV011) **only**

Model: 3523 (Sol-7 Multitherm Steam)



**replacement parts industries, inc.**

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@riparts.com](mailto:order@riparts.com) • WEBSITE [www.riparts.com](http://www.riparts.com)

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Castle/Getinge equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. WEB 03/09