



... Proud to be a **Healthcare Components Group** company ...



**HEALTHCARE
Components Group**



AUTOCLAVE LOGGER THING

The easy to use, easy to set up, wireless, automatic autoclave cycle logger.

Automatic • Wireless • Paperless

The Autoclave Logger Thing securely stores autoclave data to the cloud and on a removable SD card and eliminates the need for printing.

Currently available to fit 1st, 2nd and 3rd generation Midmark® M9 & M11 Autoclaves and Tuttnauer E Series Autoclaves.

Features:

- 15 minute set up
- 1 wire hookup
- Secure Wi-Fi
- Password protected
- Unique naming of each autoclave
- Removeable micro SD card
- FCC Certified
- 5 year time and date battery
- Comes with everything to get set up and running out of the box.

Works with Dropbox™ Cloud Storage  **Dropbox**

Learn More:

- www.autoclaveloggerthing.com
- Setup Video
- User Manual
- Automation and Integration

fits Midmark® • Ritter® M9 & M11



RPI Part #MIL337
OEM Part #9A682001

DATA LOGGER (M19)

- Autoclave Logger Thing M19

Fits: MIDMARK® - RITTER

Models: M9-001 thru -019 & M11-001 thru -019

RPI Part #MIL338
OEM Part #9A682001

DATA LOGGER (M22)

- Autoclave Logger Thing M22

Fits: MIDMARK® - RITTER

Models: M9-020 thru -022 & M11-020 thru -022

RPI Part #MIL340
OEM Part #9A682001


DATA LOGGER (M42)

- Autoclave Logger Thing M42

Fits: MIDMARK® - RITTER

Models: M9-040 thru -043 & M11-040 thru -043

fits Tuttnauer E Series



RPI Part #TUL149
OEM Part # (No OEM Part # Available)

DATA LOGGER (T01)

- The Autoclave Logger Thing T01

Fits: Autoclave

Models: E, EA, EK, EKA, EZ9, EZ10, EZ10K

LEARN MORE

