




**New Additions to our Line of Handpiece Hangers
To Fit PELTON & CRANE SPIRIT® DELIVERY UNITS**

DENTAL DELIVERY SYSTEMS/UNITS – PELTON & CRANE




RPI Part #PCK759
OEM Part #15 32 535

HANDPIECE HANGER ASSEMBLY

- White Nylon
- Photo-interrupter sensor with connector
- Includes (2) anti-rotational inserts, and (2) mounting screws

Fits: Delivery Unit Handpiece Rail
Models: Spirit EC (After Oct, 1993)




(4) #8-32 x 1/2" Socket Head Cap Screws

RPI Part #PCK776
New OEM #33 22 315 (Assy.), 15 35 546 (Top) & 15 35 827 (Bottom)

HANDPIECE HANGER ASSEMBLY (WHITE)

- White Nylon
- Photo-interrupter sensor with connector
- Includes all parts as shown

Fits: Delivery Units
Models: Spirit IIE & Spirit EC (prior to October, 1993)




(4) #8-32 x 1/2" Socket Head Cap Screws

RPI Part #PCK732
OEM Part #021814 (Assy.), 020282 (Top) & 020185 (Bottom)

HANDPIECE HANGER ASSEMBLY (BLACK)

- Black Nylon
- On-off switch with LED
- Photo-interrupter sensor with connector
- Includes all parts as shown

Fits: Delivery Units
Models: Spirit (1983-1990) & Spirit II (1990-1995)




PCV787

RPI Part #PCK783
OEM Part #33 39 970

HANDPIECE HANGER ASSEMBLY

- White Nylon
- Includes Autoholder Valve (RPI Part #PCV787), (2) anti-rotational inserts, and (2) mounting screws

Fits: Delivery Unit Handpiece Rail
Models: Spirit AC




(4) #8-32 x 1/2" Socket Head Cap Screws

RPI Part #PCK778
OEM Part #33 39 889 (Assy.), 15 24 813 (Top) & 15 35 827 (Bottom)

SYRINGE HANGER ASSEMBLY (WHITE)

- White Nylon
- Includes all parts as shown

Fits: Delivery Units
Models: Spirit IIE & Spirit EC (prior to October, 1993)




(4) #8-32 x 1/2" Socket Head Cap Screws

RPI Part #PCK741
OEM Part #021813 (Assy.), 020547 (Top) & 020185 (Bottom)

SYRINGE HANGER ASSEMBLY (BLACK)

- Black Nylon
- Includes all parts as shown

Fits: Delivery Units
Models: Spirit (1983-1990) & Spirit II (1990-1995)




RPI Part #PCV787
OEM Part #33 07 431

AUTOHOLDER VALVE

- Nickel plated brass
- Normally closed
- 1/16" barb

Fits: Handpiece Hanger Assembly
Model: Spirit AC



RPI Part #PCI753
OEM Part #020892

INSERT (SYRINGE)

- Black Nylon
- Fits:** Assistant's Holders (both manual and automatic activation units)


Models: Spirit, Spirit II, Spirit IIE, Spirit AC & Spirit EC

Visit the RPI website:
www.rpiparts.com

Order parts online:
order@rpiparts.com

FREE Tech Support
(800) 221-9723


More parts to fit Pelton & Crane Spirit® Delivery Units!



RPI Part #PCC777
OEM Part #
(No OEM Part # Available)

WATER SOLENOID COIL

- Coil Rating: 12VDC; ~44Ω; 3.28W; Coil: .990" lg. x .620" OD with .363" dia. through hole
 - Wire #: 24 AWG; 5-3/8" lg with JST connector
- Fits:** Drive Solenoid & Plug Assembly
Models: Spirit, Spirit II, Spirit IIE & Spirit EC



RPI Part #PCK754
OEM Part # *(No OEM Part # Available)*

WATER SOLENOID REPAIR KIT

- Includes all parts as shown
- Fits:** Drive Solenoid & Plug Assembly (OEM Part #15 35 371)
Models: Spirit, Spirit II, Spirit IIE & Spirit EC



RPI Part #PCC790
OEM Part #33 02 762

HANDPIECE BLOCK CONNECTOR

- 8 per package
 - .365" OD x .365" lg.; Material: Brass
 - Includes (2) O-rings (RPI Part #RPO305) per connector (installed)
- Fits:** High and Low Speed Handpiece Blocks
Models: Spirit, Spirit II, Spirit IIE & Spirit EC

RPI Part #PCS788
OEM Part #021833

AIR VALVE STEM

- 2.56" lg.; Stainless steel
 - Includes O-ring (RPI Part #RPO312)
 - Uses 3/32" hex key
- Fits:** High and Low Speed Handpiece Blocks
Models: Spirit, Spirit II, Spirit IIE & Spirit EC

RPI Part #PCS789
OEM Part #020177




WATER VALVE STEM

- 3.12" lg.; Stainless steel
 - Includes O-ring (RPI Part #RPO312)
- Fits:** High and Low Speed Handpiece Blocks
Models: Spirit, Spirit II, Spirit IIE & Spirit EC

RPI Part #PCS791
OEM Part #94 47 897

WATER VALVE STEM

- 5.60" lg.; Stainless Steel
 - Includes O-ring (RPI Part #RPO312)
 - Uses 1/8" hex key
- Fits:** Scaler or Electronic Motor Assy.
Models: Spirit IIE & EC



Listing of parts to fit Pelton & Crane Spirit® Delivery Units

RPI PART #	OEM PART #	DESCRIPTION	SPIRIT	SPIRIT II	SPIRIT IIE	SPIRIT AC	SPIRIT EC
PCA752	33 03 042/51 74 052	MICRO-SWITCH & DIAPHRAGM ASSEMBLY	•	•	•		•
PCC756	94 33 525	IN-LINE COUPLING (1/8")	•	•	•		•
PCC757	94 33 517	MALE QUICK CONNECT (1/8")	•	•	•		•
PCC772	94 33 541	FEMALE QUICK CONNECT (1/16")				•	
PCC773	94 33 533	MALE QUICK CONNECT (1/16")				•	
PCC777	<i>(No OEM Part # Available)</i>	WATER SOLENOID COIL	•	•	•		•
PCC790	33 02 762	HANDPIECE BLOCK CONNECTOR (8/pkg)	•	•	•		•
PCD748	33 02 978	REGULATOR DIAPHRAGM	•	•	•		•
PCD749	33 02 853	FOOT SWITCH DIAPHRAGM	•	•	•		•
PCG760	30 04 484	PRESSURE GAUGE				•	
PCI753	020892	INSERT (SYRINGE)	•	•	•	•	•
PCK732	021814 (Assy.) ¹	HANDPIECE HANGER ASSEMBLY (BLACK)	•	•	•		•
PCK741	021813 (Assy.) ²	SYRINGE HANGER ASSEMBLY (BLACK)	•	•	•		•
PCK754	<i>(No OEM Part # Available)</i>	WATER SOLENOID REPAIR KIT	•	•	•		•
PCK759	15 32 535	HANDPIECE HANGER ASSEMBLY					7
PCK775	51 51 779	ANTI-SIPHON VALVE KIT				5	
PCK776	33 22 315 (Assy.) ³	HANDPIECE HANGER ASSEMBLY (WHITE)			•		6
PCK778	33 39 889 (Assy.) ⁴	SYRINGE HANGER ASSEMBLY (WHITE)			•		6
PCK780	<i>(No OEM Part # Available)</i>	FOOT CONTROL VALVE KIT	•	•	•	•	•
PCK781	<i>(No OEM Part # Available)</i>	FOOT CONTROLLER SPACER KIT	•	•	•	•	•
PCK783	33 39 970	HANDPIECE HANGER ASSEMBLY				•	
PCK784	<i>(No OEM Part # Available)</i>	FOOT CONTROLLER REPAIR KIT	•	•	•		•
PCL768	30 03 402	AIR INDICATOR (GREEN)				•	
PCS750	33 02 879	CUSHION SPRING	•	•	•		•
PCS779	33 03 067	STEM	•	•	•		•
PCS788	021833	AIR VALVE STEM	•	•	•		•
PCS789	020177	WATER VALVE STEM	•	•	•		•
PCS791	94 47 897	WATER VALVE STEM			•		•
PCT785	023725	FOOT CONTROLLER TUBING ASSEMBLY (8')	•	•	•		•
PCT786	<i>(No OEM Part # Available)</i>	FOOT CONTROLLER TUBING (25')	•	•	•		•
PCV755	33 03 000	CHECK VALVE ASSEMBLY	•	•	•		•
PCV765	<i>(No OEM Part # Available)</i>	VALVE STEM				•	
PCV766	15 33 764	VALVE TOGGLE SWITCH				•	
PCV767	18 98 902	CHECK VALVE				•	
PCV787	33 07 431	AUTOHOLDER VALVE				•	
RPK628	<i>(No OEM Part # Available)</i>	COMPRESSION NUT & SLEEVE-1/8" (10/pkg)	•	•	•		•
RPO302	<i>(No OEM Part # Available)</i>	O-RING (12/pkg)	•	•	•		•
RPO305	18 97 151	O-RING (12/pkg)	•	•	•		•
RPO310	33 03 075	O-RING (12/pkg)	•	•	•	•	•
RPO312	<i>(No OEM Part # Available)</i>	O-RING (12/pkg)	•	•	•		•
RPO314	51 74 052	O-RING (12/pkg)	•	•	•		•
RPO617	70 32 741	O-RING (12/pkg)	•	•	•	•	•

1 = 020292 (Top); 020185 (Bottom) 2 = 020547 (Top); 020185 (Bottom) 3 = 15 35 546 (Top); 15 35 827 (Bottom) 4 = 15 24 813 (Top); 15 35 827 (Bottom)
 5 = Fits Spirit AC Serial #4839 and below 6 = Spirit EC prior to October, 1993 7 = Spirit EC after October, 1993



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com