



NEW PARTS TO FIT TUTTNAUER

TUTTNAUER

3mm				
	RPI Part #TUV081 OEM Part #01810101 SOLENOID VALVE	RPI Part #TUK082 OEM Part #01810903 PLUNGER VALVE KIT	RPI Part #TUC083 OEM Part #01810102 MAGNETIC COIL (10W)	RPI Part #TUC084 OEM Part #01810909 & 01810910 WIRE CONNECTOR WITH GASKET

6mm				
	RPI Part #TUV085 OEM Part #01810103 SOLENOID VALVE	RPI Part #TUK086 OEM Part #01810904 PLUNGER VALVE KIT	RPI Part #TUC087 OEM Part #01810905 MAGNETIC COIL (12W)	RPI Part #TUC084 OEM Part #01810909 & 01810910 WIRE CONNECTOR WITH GASKET

RPI Part #TUP089 OEM Part #PUM055-0005 WATER PUMP (120VAC)	RPI Part #TUP090 OEM Part #PUM055-0007 WATER PUMP (230VAC)	RPI Part #TUK091 OEM Part #SKR203-0006* PUMP MOUNTING KIT <small>(*OEM # designates 1 Rubber Mount only)</small>	RPI Part #TUE080 OEM Part #CB930040 WATER FILL ELECTRODE <i>Also Sold Separately!</i> Ferrule (PTFE) RPI Part #TUF096 OEM Part #CB847140

RPI Part #TUS088 OEM Part #01630302 WATER LEVEL SWITCH	RPI Part #RPF644 OEM Part #03140036 HEPA FILTER	RPI Part #TUC067 OEM Part #02550019 WATER RESERVOIR/FILTER COVER <i>Also Available!</i> Water Reservoir Cover with Dipstick (RPI Part #TUK075) and Dipstick (RPI Part #TUS068)	RPI Part #TUA066 OEM Part # CT312036 DOOR SWITCH ACTIVATOR

PART NAME	RPI PART #	OEM #	FITS MODELS
HEPA FILTER	RPF644	03140036	EZ9, EZ10, EZ10k, 2340EA/EKA, 2540EA/EHS/EKA & 3870EA/EHS
DOOR SWITCH ACTIVATOR	TUA066	CT312036	EZ9, EZ10, EZ10k, ValueKlave™, 1730E/EK/M/MK, 2340E/EA/EK/EKA/M/MK, 2540E/EA/EH/EHS/EK/EKA/M/MK & 3870E/EA/EH/EHS/M
WATER RESERVOIR/ FILTER COVER	TUC067	02550019	EZ9, EZ10, EZ10k, ValueKlave™, 1730E/EK/M/MK, 2340E/EA/EK/EKA/M/MK, 2540E/EA/EK/EKA/M/MK & 3870E/EA/EHS/M
MAGNETIC COIL (10W)	TUC083	01810102	EZ9, EZ10, EZ10k, 1730E/EK, 2340E/EK, 2540E/EK & 3870E (See footnote ◆ below)
WIRE CONNECTOR WITH GASKET	TUC084	01810909 & 01810910	EZ9, EZ10, EZ10k, 1730E/EK, 2340E/EK, 2540E/EK & 3870E
MAGNETIC COIL (12W)	TUC087	01810905	EZ9, EZ10, EZ10k, 1730E/EK, 2340E/EK, 2540E/EK & 3870E (See footnote ▲ below)
WATER FILL ELECTRODE	TUE080	CB930040	EZ9, EZ10, EZ10k, 1730E/EK, 2340E/EA/EK/EKA, 2540E/EA/EH/EK/EKA & 3870E/EA/EH
FERRULE (PTFE)	TUF096	CB847140	EZ9, EZ10, EZ10k, 1730E/EK, 2340E/EA/EK/EKA, 2540E/EA/EH/EK/EKA & 3870E/EA/EH
PLUNGER VALVE KIT (3mm)	TUK082	01810903	EZ9, EZ10, EZ10k, 1730E/EK, 2340E/EK, 2540E/EK & 3870E (See footnote ☉ below)
PLUNGER VALVE KIT (6mm)	TUK086	01810904	EZ9, EZ10, EZ10k, 1730E/EK, 2340E/EK, 2540E/EK & 3870E (See footnote ▲ below)
PUMP MOUNTING KIT	TUK091	SKR203-0006 (1 Rubber Mount Only)	EZ9, EZ10, EZ10k, 1730E/EK, 2340E/EA/EK/EKA, 2540E/EA/EK/EKA & 3870E/EA
WATER PUMP (120VAC)	TUP089	PUM055-0005	EZ9, EZ10, EZ10k, 1730E/EK, 2340E/EA/EK/EKA & 2540E/EA/EK/EKA
WATER PUMP (230VAC)	TUP090	PUM055-0007	1730E/EK, 2340E/EA/EK/EKA, 2540E/EA/EK/EKA & 3870E/EA
WATER LEVEL SWITCH	TUS088	01630302	EZ9, EZ10, EZ10k, 1730E, 2340E/EA/EKA, 2540E/EA/EH/EHS/EKA & 3870E/EA/EH/EHS
SOLENOID VALVE (3mm)	TUV081	01810101	EZ9, EZ10, EZ10k, 1730E/EK, 2340E/EK, 2540E/EK & 3870E (See footnote * below)
SOLENOID VALVE (6mm)	TUV085	01810103	EZ9, EZ10 & EZ10k

◆ RPI PART #TUC083: Fits Air-Vent/Exhaust Valve, Compressor and Vent Valves on Models EZ9/EZ10/EZ10k.

Fits Air-Vent/Exhaust Valve, and Fill and Vent Valves on Manifold on Models 1730E/EK, 2340E/EK, 2540E/EK AND 3870E.

▲ RPI PART #TUC087/TUK086: Fits Water Pump Valve on Models EZ9/EZ10/EZ10k. Fits Manifold Fill Valve on Models 1730E/EK, 2340E/EK, 2540E/EK AND 3870E.

☉ RPI PART #TUK082: Fits Air-Vent/Exhaust Valve, Compressor and Vent Valves on Models EZ9/EZ10/EZ10k.

Fits Air-Vent/Exhaust Valve and Vent in Manifold Block on Models 1730E/EK, 2340E/EK, 2540E/EK AND 3870E.

* RPI PART #TUV081: Fits Air-Vent/Exhaust Valve, Fill and Compressor on Models EZ9/EZ10/EZ10k. Fits Air-Vent/Exhaust Valve on Models 1730E/EK, 2340E/EK, 2540E/EK AND 3870E.