



Dental Compressor PM Kits to fit Pelton & Crane PCAir3, PCAir5 & PCAir7!


RPI Part #CMK290
OEM Part #057328



COMPRESSOR PM KIT

- For oil-less dual headed compressors
- All parts also sold separately
- Includes: (2) Filter Element (RPI Part #CME288), (1) Coalescing Filter Element (RPI Part #CME043), (1) Coalescing Filter Element (RPI Part #CME044), PM Sticker, and PM Checklist
- Also available: Air Intake Filter (RPI Part #CMF282)

Models: PCAir3




RPI Part #CMK291
OEM Part #057329



COMPRESSOR PM KIT

- For oil-less triple headed compressors
- All parts also sold separately
- Includes: (3) Filter Element (RPI Part #CME288), (1) Coalescing Filter Element (RPI Part #CME043), (1) Coalescing Filter Element (RPI Part #CME044), PM Sticker, and PM Checklist
- Also available: Air Intake Filter (RPI Part #CMF282)

Models: PCAir5




RPI Part #CMK292
OEM Part #057330



COMPRESSOR PM KIT

- For oil-less quad headed compressors
- All parts also sold separately
- Includes: (4) Filter Element (RPI Part #CME288), (1) Coalescing Filter Element (RPI Part #CME043), (1) Coalescing Filter Element (RPI Part #CME044), PM Sticker, and PM Checklist
- Also available: Air Intake Filter (RPI Part #CMF282)

Models: PCAir7




RPI Part #CME043
OEM Part #004016

COALESCING FILTER ELEMENT

- Wilkerson® brand
- Size: 1.38" OD x 2.20" lg.
- Type B; .5 micron
- Includes: Filter Element, and (1 pc) Bowl O-Ring (RPI Part #RPO685)
- See **Note** below

Models: PCAir3, PCAir5 & PCAir7




RPI Part #CME044
OEM Part #003846

COALESCING FILTER ELEMENT

- Wilkerson® brand
- Size: 1.38" OD x 2.20" lg.
- Type C; 0.01 micron
- Includes: Filter Element, and (1 pc) Bowl O-Ring (RPI Part #RPO685)
- See **Note** below

Models: PCAir3, PCAir5, & PCAir7



RPI Part #CME288
OEM Part # (OEM Part # Not Available)

AIR FILTER ELEMENT

- 50 microns
- Size: 1.30" OD x .62" ID x 1.00" lg.
- See **Note** below

Models: PCAir3, PCAir5 & PCAir7




RPI Part #CMF282
OEM Part # (OEM Part # Not Available)

AIR INTAKE FILTER

- 50 microns
- Size: 2" dia. x 3.4" lg. x 1/4" MPT
- Includes: Air Filter Element (RPI Part #CME288)
- See **Note** below

Models: PCAir3, PCAir5 & PCAir7



Note: It's the RPI Advantage! All RPI PM Kits, Filters and Filter Elements shown above include a PM Sticker that can be attached to the compressor for recording the date that the filter was replaced.



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@riparts.com • WEBSITE www.riparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit the equipment noted. All product names used in this document are trademarks or registered trademarks of their respective holders. 2005CM Show 06/20