




# New Parts to fit Dental Compressors

FOR MORE INFORMATION ABOUT THE PARTS LISTED IN THIS FLYER, SEE RPI'S WEBSITE: [www.rpiparts.com](http://www.rpiparts.com)


## START BOXES



**RPI Part #CMA029**  
**OEM Part #**  
See Chart on Back

### START BOX (1HP, 208/230V)


- Includes Start Capacitor (RPI Part #CMC030), Run Capacitor (RPI Part #CMC031) & Relay (RPI Part #CMR032)



**RPI Part #CMA033**  
**OEM Part #**  
See Chart on Back

### START BOX (1HP, 230V)


- Includes Start Capacitor (RPI Part #CMC030), Run Capacitor (RPI Part #CMC034) & Relay (RPI Part #CMR032)



**RPI Part #CMA035**  
**OEM Part #**  
See Chart on Back

### START BOX (1HP, 115V)


- Includes Start Capacitor (RPI Part #CMC036), Run Capacitor (RPI Part #CMC037) & Relay (RPI Part #CMR032)



**RPI Part #CMA039**  
**OEM Part #**  
See Chart on Back

### START BOX (3/4HP, 208/230V)

- Includes Start Capacitor (RPI Part #CMC030) & Relay (RPI Part #CMR032)




**RPI Part #CMA041**  
**OEM Part #**  
See Chart on Back

### START BOX (3/4HP, 115V)

- Includes Start Capacitor (RPI Part #CMC042) & Relay (RPI Part #CMR032)


## START CAPACITORS



**RPI Part #CMC030**  
**OEM Part #**  
See Chart on Back

### CAPACITOR


- 145-174µf, 250VAC, 50/60Hz.
- Includes all parts as shown
- Also available, Start Box (RPI Part #'s CMA029, CMA033 & CMA039)



**RPI Part #CMC036**  
**OEM Part #**  
See Chart on Back

### CAPACITOR

- 430-516µf, 115VAC, 50/60Hz.
- Includes all parts as shown
- Also available, Start Box (RPI Part #CMA035)




**RPI Part #CMC042**  
**OEM Part #**  
See Chart on Back

### CAPACITOR

- 540-648µf, 115VAC, 50/60Hz.
- Includes all parts as shown
- Also available, Start Box (RPI Part #CMA041)


## RUN CAPACITORS



**RPI Part #CMC031**  
**OEM Part #**  
See Chart on Back

### CAPACITOR


- 20µf, 440VAC, 50/60Hz.
- Also available, Start Box (RPI Part #CMA029)



**RPI Part #CMC034**  
**OEM Part #**  
See Chart on Back

### CAPACITOR

- 10µf, 370VAC, 50/60Hz.
- Includes all parts as shown
- Also available, Start Box (RPI Part #CMA033)



**RPI Part #CMC037**  
**OEM Part #**  
See Chart on Back

### CAPACITOR

- 40µf, 370VAC, 50/60Hz.
- Includes all parts as shown
- Also available, Start Box (RPI Part #CMA035)

## THERMAL PROTECTOR ASSY



**RPI Part #CMP061**  
**OEM Part #**  
See Chart on Back

### THERMAL PROTECTOR ASSY

- For use with 1HP 208-230VAC or 1HP 230VAC Compressor Heads
- Nominal opens at 135°C ± 5°C and closes at 57°C ± 9°C
- Includes all parts as shown



**RPI Part #CMP062**  
**OEM Part #**  
See Chart on Back

### THERMAL PROTECTOR ASSY

- For use with 1HP 115VAC Compressor Heads
- Nominal opens at 120°C ± 5°C and closes at 78°C ± 9°C
- Includes all parts as shown



**RPI Part #CMP063**  
**OEM Part #**  
See Chart on Back

### THERMAL PROTECTOR ASSY

- For use with 3/4HP 208-230VAC or 3/4HP 230VAC Compressor Heads
- Nominal opens at 105°C ± 5°C and closes at 61°C ± 9°C
- Includes all parts as shown




**RPI Part #CMP064**  
**OEM Part #**  
See Chart on Back

### THERMAL PROTECTOR ASSY

- For use with 3/4HP 115VAC Compressor Heads
- Nominal opens at 120°C ± 5°C and closes at 69°C ± 9°C
- Includes all parts as shown

## RELAY (UNIVERSAL)



**RPI Part #CMR032**  
**OEM Part #**  
See Chart on Back

### RELAY (UNIVERSAL)

- Includes all parts as shown
- 110-270VAC single phase
- Max voltage contact rating: 502Vac
- Motor power rating: up to 10HP
- Also available, Start Box (RPI Part #'s CMA029, CMA033, CMA035, CMA039 & CMA041)

DENTAL COMPRESSORS –

Air Techniques, Apollo/Midmark®, DentalEZ®, Matrix/Midmark®, Midmark® & Tech West/Turbine, Inc.

SECTION

4

PAGE  
408

UPDATED  
09/09



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@rpiparts.com](mailto:order@rpiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)

The above parts are manufactured by Replacement Parts Industries, Inc. to fit the equipment noted. All product names used in this document are trademarks or registered trademarks of their respective holders. Web 09/09



# New Parts to fit Dental Compressors

FOR MORE INFORMATION ABOUT THE PARTS LISTED IN THIS FLYER, SEE RPI'S WEBSITE: [www.rpiparts.com](http://www.rpiparts.com)

RPI PART #	DESCRIPTION	AIR TECHNIQUES	APOLLO/MIDMARK®	DENTALEZ®	MATRIX/MIDMARK®	MIDMARK®	TECH WEST/TURBINE, IND.
<b>NEW RELAY (UNIVERSAL)</b>							
CMR032	RELAY (UNIVERSAL)	60501 62501 64501 66501	HCA30652 HCA30653 HCA30654 HCA30655	645011010* 645011011*	77001344* 77001346*	77001344* 77001346*	PR-115 PR-208
<b>NEW RUN CAPACITORS</b>							
CMC031	CAPACITOR (20µf, 440V, 50/60Hz.)	66516	HCA30622		77001346*	77001346*	RCC-208
CMC034	CAPACITOR (10µf, 370V, 50/60Hz.)			645011011*			
CMC037	CAPACITOR (40µf, 370V, 50/60Hz.)	62503	HCA30620	645011010*	77001346*	77001346*	RCC-115
<b>NEW START BOXES</b>							
CMA029	START BOX (1HP, 208/230V)	66515	SPD95291		77001346	77001346	CRA-208
CMA033	START BOX (1HP, 230V)	66500	SPD95290	645011011	77001344		CRA-230
CMA035	START BOX (1HP, 115V)	62500	SPD95280	645011010	77001345		CRA-115
CMA039	START BOX (3/4HP, 208/230V)	64515	SPD95285 SPD95286				
CMA041	START BOX (3/4HP, 115V)	60500	SPD95275				
<b>NEW START CAPACITORS</b>							
CMC030	CAPACITOR (145-174µf, 250V, 50/60Hz.)	64502	HCA30625		77001344* 77001345* 77001346*	77001344* 77001345* 77001346*	SCC-230
CMC036	CAPACITOR (430-516µf, 115V, 50/60Hz.)	62502	HCA30628		77001345*	77001345*	SCC-115
CMC042	CAPACITOR (540-648µf, 115V, 50/60Hz.)	60502	HCA30626				
<b>NEW THERMAL PROTECTOR ASSEMBLIES</b>							
CMP061	THERMAL PROTECTOR ASSEMBLY	66504	HCA30615	No OEM Part # available	61894802	61894802	KOPC-230
CMP062	THERMAL PROTECTOR ASSEMBLY	62504	HCA30610	No OEM Part # available	61894801	61894801	KOPC-230
CMP063	THERMAL PROTECTOR ASSEMBLY	64503	HCA30605				
CMP064	THERMAL PROTECTOR ASSEMBLY	64503	HCA30600				

\* The OEM does not list an individual part number, however the part is included in the OEM Part Number noted.

## Don't Forget ...

RPI offers even more of the dental compressor parts you need to fit Air Techniques, Apollo/Midmark®, DentaleZ®, Matrix/Midmark®, Midmark®, Tech West/Turbine, & more!

- Air Intake Filters
- Coalescing Filter Assemblies
- Cooling Fans
- Pressure Gauges
- Pressure Switches and Valves
- Purge Filters and Mufflers
- Safety Valves and more

