More New Parts to fit Dental Vacuum Units!


For more information about the parts listed in this flyer, see our website www.rpiparts.com

RPI Part #VPV119 OEM Part # See Chart Below

CHECK VALVE (1")

- Cracking Pressure: 0.5 PSI
- Port: 1" FPT
- Material: Bronze body; RPTFE seal
- Size: 1.93" dia. x 3.50" lg.

Fits: Motor Assembly




RPI Part #VPV120 OEM Part # See Chart Below

CHECK VALVE (1-1/2")

- Ports: 1-1/2" FPT
- · Material: Clear PVC
- Size: 3.38" dia. x 5.75" lq.

Fits: Drain Assembly



RPI Part #VPV122 OEM Part # See Chart Below

CHECK VALVE (1-1/2")

- Ports: 1-1/2" pipe compression fitting
- Material: PVC; Nitrile (Buna N) seal
- Size: 2.94" dia. x 5.25" lg.



RPI Part #VPC136 OEM Part # See Chart Below

SOLIDS COLLECTOR (1-1/2")

- Ports: 1-1/2" FPT
- Material: Black plastic body with clear plastic bowl
- Includes: Bowl, Bowl Gasket (RPI Part #VPG081) and Screen (20 Mesh) (RPI Part #VPS140)
- Also available: Screen (50 Mesh) (RPI Part #VPS139)

Fits: Separator Assembly



RPI Part #VPK137 OEM Part # See Chart Below

SOLIDS COLLECTOR REPLACEMENT KIT

• Includes all parts as shown Fits: Solids Collector (1-1/2")




RPI Part #VPK138 OEM Part # See Chart Below

SOLIDS COLLECTOR REPLACEMENT KIT (CASE of 8)

 Includes: (8) Solids Collector Replacement Kits (RPI Part #VPK137)





RPI Part #VPG081 0EM Part # See Chart Below

BOWL GASKET

- 3 per package
- Material: Buna N (Nitrile)
- Size: 3.13" OD x .13" Square
 Fits: Solids Collector Assembly

Fits: Solids Collector Assembly (1", 1-1/4" & 1-1/2")



RPI Part #VPS140 OEM Part # See Chart Below

SCREEN (20 MESH)


- Material: Stainless steel screen: Red plastic cage
- Size: 2.47" OD x 5.46" lg. Fits: Solids Collector (1-1/2")



RPI Part #VPS139 OEM Part # See Chart Below

SCREEN (50 MESH)

- Material: Stainless steel screen: Blue plastic cage
- Size: 2.47" OD x 5.46" lg. Fits: Solids Collector (1-1/2")



RPI Part #	DESCRIPTION	AIR Techniques	APOLLO/ Midmark®	DENTAL <i>ez</i> ®/Custom Air/Ramvac	MATRX/ Midmark®	MIDMARK®
VPC136	Solids Collector (1-1/2")	H5181 & H5216	PVC50700 & AVAVPF2		62937303	62937303
VPG081	Bowl Gasket - 3/PKG	No OEM Part # Available	MRP70936	No OEM Part # Available	77005158	77005158
VPK137	Solids Collector Replacement Kit	H5217	AVA60031		77005169	77005169
VPK138	Solids Collector Replacement Kit (Case of 8)	No OEM Part # Available	AVA60030		77005179	77005179
VPS139	Screen (50 Mesh)	No OEM Part # Available	No OEM Part # Available		No OEM Part # Available	No OEM Part # Available
VPS140	Screen (20 Mesh)	No OEM Part # Available	No OEM Part # Available		77005079	No OEM Part # Available
VPV119	Check Valve (1")					014-0552-01
VPV120	Check Valve (1-1/2")	54075				014-0551-00
VPV122	Check Valve (1-1/2")	54139				

rpi


replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

SECTION 4

PAGE

UPDATED 03/16



RPI Part #VPB100 **OEM Part # See Chart On Page 3** Industry Part #2/7M580JB

DRIVE BELT

• Effective Belt Length: 22.83"

• Size: .63" @ top x .31" thk Fits: Motor and Vacuum Pump




VPB100

CROSS SECTION SHOWN ACTUAL SIZE

RPI Part #VPB096 **OEM Part** # See Chart On Page 3 Industry Part #4L370-A35F

DRIVE BELT

- V-Belt
- Effective Belt Length: 37.0"
- Size: .50" @ top x .31" thk Fits: Motor and Vacuum Pump




CROSS SECTION SHOWN ACTUAL SIZE

RPI Part #VPB097 **OEM Part # See Chart On Page 3** Industry Part #4L380-A36F

DRIVE BELT


- V-Belt
- Effective Belt Length: 38.0"
- Size: .50" @ top x .31" thk Fits: Motor and Vacuum Pump



RPI Part #VPB098 **OEM Part # See Chart On Page 3** Industry Part #2/7M650JB

DRIVE BELT

- Effective Belt Length: 25.59"
- Size: .63" @ top x .31" thk Fits: Motor and Vacuum Pump



RPI Part #VPB099 **OEM Part # See Chart On Page 3** Industry Part #2/7MB670JB

DRIVE BELT

- Effective Belt Length: 26.38"
- Size: .63" @ top .31" thk Fits: Motor and Vacuum Pump



RPI Part #VPB101 **OEM Part** # See Chart On Page 3 Industry Part #2/7M600JB

DRIVE BELT

- Effective Belt Length: 23.62"
- Size: .63" @ top x .31" thk

Fits: Motor and Vacuum Pump



RPI Part #VPB102 **OEM Part # See Chart On Page 3** Industry Part #2/7M545JB

DRIVE BELT

- Effective Belt Length: 21.46"
- Size: .63" @ top x .31" thk Fits: Motor and Vacuum Pump



VPB102 **CROSS SECTION** SHOWN ACTUAL SIZE

SHOWN ACTUAL SIZE



Page 2



RPI Part #VPB103 **OEM Part # See Chart Below Industry Part #A41**



DRIVE BELT

- V-Belt; Effective Belt Length: 43.0"
- Size: .50" @ top x .31" thk
- All units require 2 Belts

Fits: Motor and Vacuum Pump






RPI Part #VPB104 **OEM Part # See Chart Below** Industry Part #A42



DRIVE BELT

- V-Belt; Effective Belt Length: 44.0"
- Size: .50" @ top x .31" thk
- All units require 2 Belts

Fits: Motor and Vacuum Pump



RPI Part #VPB105 **OEM Part # See Chart Below Industry Part #A40**




. ********

VPB105

DRIVE BELT

- V-Belt; Effective Belt Length: 42.0"
- Size: .50" @ top x .31" thk
- All units require 2 Belts

Fits: Motor and Vacuum Pump




RPI Part #VPK128 **OEM Part # See Chart Below**



- Includes all parts as shown
- Size: .25" x .25" x .69" lg.
- Radius: .38"; Spring Length: 1-1/16"
- OEM recommends replacement every 300 hours of use
- Also available: Pump Repair Kit (RPI Part #VPK127)


Fits: Liquid Evacuation Pump



RPI Part #VPK127 **OEM Part # See Chart Below**


PUMP REPAIR KIT

- Includes: (1) Impeller and (1) gasket
- Also available: Brush Kit (RPI Part #VPK128) Fits: Liquid Evacuation Pump



RPI Part #VPL131 **OEM Part # See Chart Below**





• 15W-50 (ISO grade 32)

Fits: Vacuum Unit

ALSO AVAILABLE BY THE CASE!

RPI Part #VPL132 **OEM Part # See Chart Below**

SYNTHETIC OIL (MOBIL 1®) (CASE of 6)

• Case includes: (6) 1 Quart (.946 Liter) bottles



RPI Part #VPL133 **OEM Part # See Chart Below**

OIL-SYN GEARBOX (1 Qt.)

- Size: 1 Quart (.946 Liter) bottle
- Mobil SHC 630 (ISO grade 220)

Fits: Vacuum Unit

ALSO AVAILABLE BY THE CASE!

RPI Part #VPL134 **OEM Part # See Chart Below**

OIL-SYN GEARBOX (CASE of 12)

· Case includes: (12) 1 Quart (.946 Liter) bottles



RPI Part #	DESCRIPTION	DENTAL <i>ez</i> ®/Custom Air/Ramvac	MATRX/ Midmark®	MIDMARK®	PELTON & Crane
RPT954	Belt Tensioner Tester	No OEM Part # Available		No OEM Part # Available	No OEM Part # Available
VPB096	Drive Belt			016-1055-00 & 2021110	
VPB097	Drive Belt			016-1056-00	
VPB098	Drive Belt	000273SP			
VPB099	Drive Belt	000274SP			
VPB100	Drive Belt	002035SP			
VPB101	Drive Belt	002027SP			
VPB102	Drive Belt	002034SP			
VPB103	Drive Belt	0002163SP			
VPB104	Drive Belt	0002164SP			
VPB105	Drive Belt	0002360SP			
VPK127	Pump Repair Kit		77001507	77001507	
VPK128	Brush Kit		77001506	77001506	
VPL131	Synthetic oil (Mobil 1®)	No OEM Part # Available			
VPL132	Synthetic oil (Mobil 1®) (Case of 6)	380706			
VPL133	Oil-Syn Gearbox (1 Qt.)			064-0028-01	
VPL134	Oil-Syn Gearbox			No OEM Part # Available	

RPI offers these parts and more to fit Dental Vacuum Units!

Don't Forget ...

- Adaptors and Connectors
- Anti-Syphon Valves
- Anti-Syphon Valve Repair Kits
- Check Valves
- Cuspidor Drain Valve (Vada Trap)
- Diaphragm
- Filter Elements
- Flow Regulators
- Manifold Gaskets

- Sink Drain Valve Assembly
- Solenoid Valve Assemblies
- Solenoid Valve Repair Kits
- Strainers
- Strainer Elements
- Vacuum Gauges
- Vacuum Relief Valves
- Water Control Assemblies
- and even more new parts to come!



replacement parts industries, inc.

20338 Corisco Street • Chatsworth, California • 91311

ADDRESS SERVICES REQUESTED

PRESORTED STANDARD U.S. POSTAGE PAID VAN NUYS, CA PERMIT NO. 2399

More New Parts to fit Dental Vacuum Units!



RPI Part #VPK123 OEM Part # See Chart Below

TANK WASHOUT KIT

Includes: Adaptor (3/8" Push-In x 1/2" MPT) and (10 ft.) Tubing (1/4" ID x 3/8" OD Clear Polyethylene) (RPI Part #RPT959)
 Fits: CAS Tank (Cyclonic Action Separator)





RPI Part #RPT959 OEM Part # See Chart Below

TUBING (1/4" ID x 3/8" OD CLEAR POLYETHYLENE)

• Sold by the foot (reel not included)

Fits: 3/8" Polyflo Fitting or 1/4" Barb Fitting

RPI PART#	DESCRIPTION	AIR Techniques	MATRX/ Midmark®	MIDMARK®
RPT959	Tubing (1/4" ID x 3/8" OD Clear Poly)	54601	77002004	77002004
VPK123	Tank Washout Kit		77008009 & 77002004	77008009 & 77002004