




# Foot Controls to fit Dentsply Cavitrion<sup>®</sup> Ultrasonic Scalers!



3 PIN CONFIGURATION




RPI Part #DSF039  
OEM Part #77885

**FOOT CONTROL**

- Single switch position
- 8 foot long gray #24 AWG shielded cord with 3 pin round plug
- Size: 2.63" W x 3.53" L x .85" H
- Non-skid base and pedal for positive traction

**Fits:** Cavitrion G115

**Models:** Cavitrion<sup>®</sup> BOBCAT<sup>®</sup> 25 KHz - G115




RPI Part #DSF012  
OEM Part #565162005

**FOOT CONTROL**

- 2 Switch boost feature
- 8 foot long black #22 AWG shielded cord with .097" diameter pin molded plug
- Size: 3.44" W x 4.00" L x 1.09" H
- Non-skid base and raised ribs on pedal for positive traction

**Fits:** Cavitrion G124

**Models:** Cavitrion<sup>®</sup> Select<sup>™</sup> SPS<sup>™</sup> 30K - G124




RPI Part #DSF040  
OEM Part #565162002

**FOOT CONTROL**

- Single switch position
- 8 foot long black #22 AWG shielded cord with .097" diameter pin molded plug
- Size: 2.63" W x 3.53" L x .85" H
- Non-skid base and pedal for positive traction

**Fits:** Cavitrion G121

**Models:** Cavitrion<sup>®</sup> Select<sup>™</sup> Scaler 25K - G121



RPI Part #DSF013  
OEM Part #565138002, 565138003


**FOOT CONTROL**

- 2 Switch boost feature
- 8 foot long gray #24 AWG shielded cord with 4 pin round plug
- Size: 3.44" W x 4.00" L x 1.09" H
- Non-skid base and raised ribs on pedal for positive traction

**Fits:** Cavitrions G119, G120, & G122

**Models:** Cavitrion<sup>®</sup> SPS<sup>™</sup> Ultrasonic Scaler - G119, Cavitrion<sup>®</sup> JET with SPS<sup>™</sup> - G120, and Cavitrion<sup>®</sup> Prophy-JET<sup>®</sup> Air Polishing System - G122

4 PIN CONFIGURATION



RPI Part #DSF042  
OEM Part #565162006

**FOOT CONTROL**

- Single switch position
- 8 foot long gray #24 AWG shielded cord, three wire with 3/16" quick connect terminals
- Size: 2.63" W x 3.53" L x .85" H
- Non-skid base and pedal for positive traction

**Fits:** Cavitrion G130

**Models:** Cavitrion<sup>®</sup> BOBCAT<sup>®</sup> Pro - G130



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@ripiparts.com](mailto:order@ripiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Dentsply equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Flyer Code 1705DS Rd #2 Show 09/18



**Complete list of RPI Parts to fit  
Dentsply Cavitrone<sup>®</sup>  
Ultrasonic Scalers!**

RPI Part #	OEM Part #	Description	BOBCAT <sup>®</sup> 25 KHz - G115	BOBCAT <sup>®</sup> Pro - G130	DualSelect <sup>™</sup> Dispensing System - G118	JET Plus <sup>™</sup> - G132	JET Plus <sup>™</sup> - G137	JET with SPS <sup>™</sup> - G120	Plus <sup>™</sup> Scaler - G131	Plus <sup>™</sup> Scaler - G136	Prophecy-JET <sup>®</sup> Air Polishing System - G122	Prophecy-JET <sup>®</sup> - G138	Select <sup>™</sup> Scaler 25K - G121	Select <sup>™</sup> SPS <sup>™</sup> 30K - G124	SPS <sup>™</sup> Surface Mount Ultrasonic Scaler - G98a	SPS <sup>™</sup> Ultrasonic Scaler - G119	Touch <sup>™</sup> Scaler - G1000
ADB012	No OEM Part # Available	Barb - 10/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
ADC034	No OEM Part # Available	Uni-Clamp - 20/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
ADC035	No OEM Part # Available	Sleeve Clamp - 20/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
ADC058	No OEM Part # Available	Cord Clip - 10/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
ADS032	No OEM Part # Available	Polyflo Sleeve/Ferrule (1/4") - 25/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
ADW038	No OEM Part # Available	Flat Washer - 50/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSA002	60367	Duckbill Filter Assembly	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSA008	63856	Hose Assembly & Filter (Air)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSA046	No OEM Part # Available	Water Supply Hose Assembly (Braided)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSA047	90125	Water Supply Line Assembly	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSC017	No OEM Part # Available	Coil (Solenoid Valve)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSF005	90158	Disc Filter (Water) - 10/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSF012	565162005	Foot Control	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSF013	565138002, 565138003	Foot Control	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSF039	77885	Foot Control	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSF040	565162002	Foot Control	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSF042	565162006	Foot Control	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSF048	623072001	Union Connector (3/16" Push-In)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSH003	60388	Water Hose (Braided Blue)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSH009	80786	Air Hose (Braided Gray)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSK001	90172	Retrofit Water Filter Kit	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSK016	No OEM Part # Available	Solenoid Valve Repair Kit	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSK018	60367 (complete Assy.)	Repair Kit (Duckbill Filter Assembly)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSK020	No OEM Part # Available	Solenoid Valve Repair Kit	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSK021	No OEM Part # Available	Water Filter Kit	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSK022	81167	Check Valve Kit	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSL024	No OEM Part # Available	Label (Air) - 3/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSL025	No OEM Part # Available	Label (Water) - 3/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSL026	No OEM Part # Available	Label (Caution) - 3/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DST014	574335001, 61350	Transistor	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSV006	571070003	Solenoid Valve (Water)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSV011	629194002	Solenoid Valve (Pinch)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSV015	571070002	Solenoid Valve (Water)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSV019	81258	Solenoid Valve (Water)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSV023	629211001	Check Valve	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DSV027	629202001	Needle Valve (Water)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPA915	90088	Air Filter Assembly	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPC557	554102003	Hospital Grade Power Cord (13A @ 125VAC, 8 ft.)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPC868	554102003	Hospital Grade Power Cord (13A @ 125VAC, 1.5 ft.)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPC869	554102003	Hospital Grade Power Cord (13A @ 125VAC, 3 ft.)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPE916	No OEM Part # Available	Filter Element	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF524	557034018	Fuse (1/8A, 250V) 5mm x 20mm Time Delay - 5/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF883	No OEM Part # Available	Polyflo Nut (1/4" Brass) - 5/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF907	621066001	Luer Lock (Male x 5/32" Barb) - 10/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF908	621065001	Luer Lock (Female x 5/32" Barb) - 10/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF913	No OEM Part # Available	Polyflo Fitting (1/4" Tube x 1/8" MPT) - 5/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF914	No OEM Part # Available	Adaptor (1/16" Barb x 1/8" MPT) - 5/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF917	No OEM Part # Available	Polyflo Fitting (1/4" Tube x #10-32 Male) - 5/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF918	79521	Quick Disconnect, Male (1/4" Poly Auto Shut-Off)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF920	No OEM Part # Available	Quick Disconnect, Female (Panel Mount 1/4" Poly)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF922	60537	Fuse (630MA, 250V) 5mm x 20mm Time Delay - 5/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF923	63827	Adaptor Elbow (1/16" Barb X 1/8" MPT) - 5/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF942	623064002	Cross Fitting (1/8" Barb x #10-32 Male) - 5/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF944	621068001	Adaptor Fitting (1/16" Barb x 1/8" Barb) - 10/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF948	No OEM Part # Available	Adaptor (1/16" Barb x 1/8" MPT) - 5/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF953	No OEM Part # Available	Polyflo Nut (1/4" Chrome) - 5/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF974	No OEM Part # Available	Luer Lock (Male x 1/8" Barb) - 10/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF975	621065001	Luer Lock (Female x 1/8" Barb) - 10/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF976	No OEM Part # Available	Male Connector (3/16" Push-In x #10-32 Male) - 2/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF977	6390253	Elbow Fitting (1/16" Barb x #10-32 Male) - 10/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPF978	No OEM Part # Available	Swivel Elbow (1/8" Barb x #10-32 Male) - 2/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPO304	No OEM Part # Available	O-Ring - 12/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPO342	62351	O-Ring (Black) - 12/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPO924	62605	O-Ring (Green) - 12/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPT296	No OEM Part # Available	Universal Sleeve Tool	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPT909	No OEM Part # Available	Braided Tubing (1/4" OD Blue) - per ft.	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPT910	No OEM Part # Available	Braided Tubing (1/4" OD Blue) (Box of 100 ft.)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPT911	No OEM Part # Available	Braided Tubing (1/4" OD Gray) - per ft.	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPT912	No OEM Part # Available	Braided Tubing (1/4" OD Gray) (Box of 100 ft.)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPT927	No OEM Part # Available	Tubing (.160" ID x 1/4" OD Blue) - per ft.	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPT940	61963	Tubing (1/16" ID x 1/8" OD Clear Polyurethane) - per ft.	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPT941	625075001	Tubing (.107" ID x 3/16" OD Blue) - per ft.	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPT973	625072001	Tubing (1/16" ID x 3/16" OD Black) - per ft.	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RPT999	No OEM Part # Available	Tubing (.160" ID x 1/4" OD Blue) (Box of 100 ft.)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RXT001	625075001	Tubing (.107" ID x 3/16" OD Blue) (Box of 100 ft.)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RXT008	712-00004	Cable Tie (.270" W x 8.75" lg. White) - 25/pkg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
RXT016	61631	Pinch Tubing (1/16" ID x 1/8" OD Red) - per ft.	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

**Call Us Toll-Free Monday-Friday, 8:00 AM-4:30 PM (Pacific Time), (800) 221-9723**