



New Parts To Fit Air Techniques, Pelton & Crane, and SciCan

SciCan: Statim 5000

RPI Part #SCS029
OEM Part #01-101649S




CASSETTE SEAL, 5000


- Black Silicone
- Includes: RPS287 (4) Pouches of Liquid Soap

Fits: Cassette Lid

Models: STATIM 5000



Air Techniques: Peri-Pro I, II, & III




RPI Part #ATF649
OEM Part #94278

THERMO FUSE (98°C)

Fits: Peri-Pro I & II

The Peri-Pro I & II uses (OEM Part #94278)
RPI Part#ATF649, Rated @ 98° C.



RPI Part #ATF651
OEM Part #90278



THERMO FUSE (109°C)

Fits: Peri-Pro III

The Peri-Pro III uses (OEM Part #90278)
RPI Part #ATF651, Rated @ 109°C.

Pelton & Crane: Chairman & Coachman

RPI Part #PCB730
OEM Part #040988




FOOT SWITCH ASSY. PCB

- Includes Actuator Spacer

Fits: Foot Switch Assy.

Models: Chairman 5000



Actuator Spacer
RPI Part# PCA726


RPI Part #PCB723
New OEM Part #018326



SWITCH BACK PCB ASSEMBLY

Fits: Chair back assembly

Models: Coachman




RPI Part #PCB736
New OEM Part #94 34 424



UP INTERCONNECT (PCB)

Fits: Incline/recline controls

Models: Chairman 5000 & 5010



RPI Part #PCA726
OEM Part #5428292



ACTUATOR SPACER

- Closed Cell Neoprene
- 4 per package

Fits: Footswitch PCB Assy.

OEM#3040988, RPI# PCB730,

Models: Chairman 5000



RPI Part #PCB733
New OEM Part #30 06 629



LIMIT SWITCH (PCB)

Fits: Incline/recline and upper traverse

Models: Chairman 5090

