




New Parts to fit Hill-Rom Stabilet Infant Warmer

RPI Part #AIB116
OEM Part #BDZ107/BDZ108

RPI Added Value!

FLUORESCENT BALLAST

• Generic part; can be used to replace right hand or left hand ballast
Fits: Overhead Light Assy
Models: Stabilet



RPI Part #AIH118
OEM Part #BCR730

FLUORESCENT TUBE HOLDER


Fits: Overhead Light Assy.
Models: Stabilet



RPI Part #AIE115
OEM Part #BDJ002 & BDJ003

ENCLOSURE, FLUORESCENT TUBE


Fits: Overhead Light Assy.
Models: Stabilet



RPI Part #LMP021
OEM Part #BDJ001

FLUORESCENT LAMP

Fits: Overhead Light Assy
Models: Stabilet




RPI Part #AIR112
OEM Part #BDZ205

RPI Added Value!

LIGHT RELAY ASSY


Fits: Overhead Light Assy.
• Includes installed timing resistor
Models: 2000A/B, 2000C, 1250



RPI Part #AIR124
OEM Part #BCZ952

RECEPTACLE ASSY, PROBE


Fits: Chassis Assy
Models: Stabilet



RPI Part #AIB109
OEM Part #BKK010

BASSINET BRACKET/HINGE

• Beige
Fits: Bassinet
Models: Stabilet 2000A/B, 2000C, 1250



RPI Part #AIB110
OEM Part #KB010

BASSINET BRACKET/HINGE

• Black
Fits: Bassinet
Models: Stabilet 2000A/B, 2000C, 1250



RPI Part #AIB111
OEM Part #BKK015

CORNER BUMPER


Fits: Base Assy.
Models: 2000A/B, 2000C



RPI Part #AIC119
OEM Part #BPM002

CASTER - 5" SWIVEL WITH BRAKE


Fits: Base Assembly
Models: Stabilet



RPI Part #AIC120
OEM Part #BPM002

CASTER - 5" SWIVEL (NO BRAKE)


Fits: Base Assembly
Models: Stabilet



RPI Part #AIS113
OEM Part #BC0003

THERMAL SWITCH


Fits: Overhead Light Assy.
Models: 2000A/B, 2000C, 1250



RPI Part #AIT114
OEM Part #BCP425 & BEZ504

FAST PAK TRIAC ASSY WITH CABLE


Fits: Chassis Assy.
Models: 2000A/B, 2000C, 1250



RPI Part #RPC464
OEM Part # NA

HEAT SINK COMPOUND

• MSDS and Technical Data sheet available on RPI website
• Non-Silicone based thermal grease
• 5 per package
• Size: .07 fl. oz (2ml)




RPI Part #AIS122
OEM Part #BDB001

RPI Added Value!

SWITCH, MAIN POWER


• Includes AIL123 and LMP022
Fits: Front Panel
Models: Stabilet 2000C



RPI Part #AIL123
OEM Part #BDB002

LENS, MAIN POWER SWITCH


Fits: Main Power Switch - RPI Part #AIS122
Models: Stabilet 2000C



RPI Part #LMP022
OEM Part #BDB003

LAMP

• 14 VDC, 1.1W, Base: Wedge, T-1 3/4
• Incandescent
Fits: Main Power Switch - RPI Part #AIS122
Models: Stabilet 2000C



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Hill-Rom Air Shields equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. WEBSITE 02/05