





## New Line of Parts to fit CASTLE/GETINGE BULK STERILIZERS




**SOL-1  
SOLENOID VALVE  
ASSEMBLY**

RPI Part #CSV004  
OEM Part #517339  
• Steam to Chamber  
• 3522, 3525 & 3533




**SOLENOID VALVE REPAIR KIT**  
RPI Part #AMK177  
OEM Part #61301601739




**SOL-1  
SOLENOID VALVE  
ASSEMBLY**

RPI Part #CSV005  
OEM Part #517345  
• Steam to Chamber  
• 3523




**SOLENOID VALVE REPAIR KIT**  
RPI Part #AMK183  
OEM Part #54190




**SOL-3  
SOLENOID VALVE  
ASSEMBLY**

RPI Part #CSV007  
OEM Part #517347  
• Air-In  
• 3523 & 3533




**SOLENOID VALVE REPAIR KIT**  
RPI Part #CSK018  
OEM Part #517295




**SOL-4  
SOLENOID VALVE  
ASSEMBLY**

RPI Part #CSV004  
OEM Part #517339  
• Multitherm Bypass  
• 3523, 3623 & 3633




**SOLENOID VALVE REPAIR KIT**  
RPI Part #AMK177  
OEM Part #61301601739




**SOL-5  
SOLENOID VALVE  
ASSEMBLY**

RPI Part #CSV009  
OEM Part #517354  
• Slow Exhaust  
• 3523




**SOLENOID VALVE REPAIR KIT**  
RPI Part #CSK020  
OEM Part #517298

**SOL-5 & 6  
SOLENOID VALVE ASSEMBLY**  
• Slow Exhaust & Vac-Water  
*Sol. Valve Assy. not yet available;  
Repair Kit is available, see below.*




**DUAL SOLENOID VALVE REPAIR KIT**  
RPI Part #CSK019  
OEM Part #61301601738  
• 3522, 3533, 3622, 3623 & 3633




**SOL-6  
SOLENOID VALVE  
ASSEMBLY**

RPI Part #CSV010  
OEM Part #517349  
• Vac-Water  
• 3523 & 3525




**SOLENOID VALVE REPAIR KIT**  
RPI Part #AMK174  
OEM Part #56161

**SOL-7  
SOLENOID VALVE ASSEMBLY**  
• Multitherm Steam  
*Sol. Valve Assy. not yet available;  
Repair Kit is available, see below.*




**SOLENOID VALVE REPAIR KIT**  
RPI Part #CSK022  
OEM Part #517297  
• 3523




**SOL-7  
SOLENOID VALVE  
ASSEMBLY**

RPI Part #CSV004  
OEM Part #517339  
• Multitherm Steam  
• 3622, 3623 & 3633




**SOLENOID VALVE REPAIR KIT**  
RPI Part #AMK177  
OEM Part #61301601739




**SOL-8  
SOLENOID VALVE  
ASSEMBLY**

RPI Part #CSV004  
OEM Part #517339  
• Jacket  
• 3522, 3523 & 3533




**SOLENOID VALVE REPAIR KIT**  
RPI Part #AMK177  
OEM Part #61301601739




**SOL-8  
SOLENOID VALVE  
ASSEMBLY**

RPI Part #CSV014  
OEM Part #517342  
• Jacket  
• 3622, 3623 & 3633




**SOLENOID VALVE REPAIR KIT**  
RPI Part #CSK025  
OEM Part #54186



**SOL-10 & 11  
SOLENOID VALVE  
ASSEMBLY**

RPI Part #CSV012  
OEM Part #541582  
• Door  
• 3522, 3523, 3525  
& 3533



**SOLENOID VALVE REPAIR KIT**  
RPI Part #CSK023  
OEM Part #61301600442



**replacement parts industries, inc.**


CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@riparts.com](mailto:order@riparts.com) • WEBSITE [www.riparts.com](http://www.riparts.com)

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Castle/Getinge equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Web 02/08

**QUICK  
REFERENCE!**

**SOLENOID VALVE ASSEMBLIES & REPAIR KITS  
TO FIT CASTLE/GETINGE BULK STERILIZERS!**

SOLENOID VALVE ASSY. FITS	MODELS	SOLENOID VALVE ASSY. RPI PART #	REPAIR KIT RPI PART #
Sol-1 Steam to Chamber	3522, 3525 & 3533	CSV004	AMK177
Sol-1 Steam to Chamber	3523	CSV005	AMK183
Sol-3 Air-In	3523 & 3533	CSV007	CSK018
Sol-4 Multitherm Bypass	3523, 3623 & 3633	CSV004	AMK177
Sol-5 Slow Exhaust	3523	CSV009	CSK020
Sol-5 & 6 Slow Exhaust & Vac Water	3522, 3533, 3622, 3623 & 3633	<i>Sol. Valve Assy. not yet available</i>	CSK019
Sol-6 Vac-Water	3523 & 3525	CSV010	AMK174
Sol-7 Multitherm Steam	3523	<i>Sol. Valve Assy. not yet available</i>	CSK022
Sol-7 Multitherm Steam	3622, 3623 & 3633	CSV004	AMK177
Sol-8 Jacket	3522, 3523 & 3533	CSV004	AMK177
Sol-8 Jacket	3622, 3623 & 3633	CSV014	CSK025
Sol-10 & 11 Door	3522, 3523, 3525 & 3533	CSV012	CSK023



**See us on the web for more information...[www.rpiparts.com](http://www.rpiparts.com)**