


Models: 100 Series, 400 Series, 500 Series, 700 Series, 3522, 3523, 3525, 3533, 6300 Series, 6400 Series, 6900 Series, & MC3600

RPI is Your Best Source for Door Gaskets to fit Getinge/Castle Bulk Sterilizers!

GETINGE/CASTLE




RPI Part #CSG102
OEM Part #61301607691

DOOR GASKET (17.50" x 17.50")

- Material: Blue silicone
 - Custom formulated for greater lubricity!
 - .65" OD cross section
- Fits:** Chamber
Models: 400 Series



CSG102
CROSS SECTION
NOT SHOWN ACTUAL SIZE



RPI Part #CSG002
OEM Part #61301605017 & 58397


DOOR GASKET (16" X 16")

- Material: Blue silicone
- New material with greater lubricity!
- .43" OD cross section

Fits: 16" x 16" Door
Models: 3522, 3525 & 3533



CSG002
CROSS SECTION
NOT SHOWN ACTUAL SIZE



RPI Part #CSG103
OEM Part #61301607692


DOOR GASKET (21" x 21")

- Material: Blue silicone
- Custom formulated for greater lubricity!
- .65" OD cross section

Fits: 21.50" x 21.50" Door
Models: 500 Series



CSG103
CROSS SECTION
NOT SHOWN ACTUAL SIZE



RPI Part #CSG030
OEM Part #513966


DOOR GASKET (16" x 16")

- Material: Blue silicone
- New material with greater lubricity!
- .40" x .40" "U" shaped cross section

Fits: 16" x 16" Door
Models: 3522, 3525 & 3533



CSG030
CROSS SECTION
NOT SHOWN ACTUAL SIZE



RPI Part #CSG104
OEM Part #61301609906


DOOR GASKET (26.50" x 36")

- Material: Blue silicone
- Custom formulated for greater lubricity!
- .65" OD cross section

Fits: 26.5" x 36.00" Door
Models: 700 Series



CSG104
CROSS SECTION
NOT SHOWN ACTUAL SIZE



RPI Part #CSG001
OEM Part #61301605016 & 58398


DOOR GASKET (20" X 20")

- Material: Blue silicone
- New material with greater lubricity!
- .43" OD cross section

Fits: 20" x 20" Door
Models: 100 Series & 3523



CSG001
CROSS SECTION
NOT SHOWN ACTUAL SIZE



RPI Part #CSG031
OEM Part #513967

DOOR GASKET (20" x 20")

- Material: Blue silicone
- New material with greater lubricity!
- .40" x .40" "U" shaped cross section

Fits: 20" x 20" Door
Models: 100 Series & 3523



CSG031
CROSS SECTION
NOT SHOWN ACTUAL SIZE

See reverse side for more!


Warning: These gaskets are not for use in Ethylene Oxide (E.O. Gas) sterilizers

rpi replacement parts industries, inc.
CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@riparts.com • WEBSITE www.riparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Getinge/Castle equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. 1109CS Show 05/17

SECTION 10
PAGE 310
UPDATED 05/17

RPI is Your Best Source for Door Gaskets to fit Getinge/Castle Bulk Sterilizers!



RPI is pleased to introduce 3 new Door Gaskets to fit the following Getinge/Castle Bulk Sterilizers:

- 400 Series
- 500 Series
- 700 Series

Plus, RPI also offers Door Gaskets that are now custom formulated for greater lubricity to fit the following Models:

- 100 Series
- 6300 Series
- 6900 Series
- 3522, 3523, 3525 & 3533
- 6400 Series
- MC3600

RPI is Your Best Source for Door Gaskets to fit Getinge/Castle Bulk Sterilizers!


ADDRESS SERVICES REQUESTED

PREPARED STANDARD
U.S. POSTAGE PAID
VAN NUYS, CA
PERMIT NO. 2399

replacement parts industries, inc. 
20338 Corisco Street • Chatsworth, California • 91311




RPI Part #CSG003
OEM Part #40441
DOOR GASKET (24" x 36")
• Material: Blue silicone
• New material with greater lubricity!
• .43" OD cross section
Fits: 24" x 36" Door
Models: 6300 Series, 6400 Series & 6900 Series




CSG003
CROSS SECTION
NOT SHOWN ACTUAL SIZE



RPI Part #CSG032
OEM Part #536651
DOOR GASKET (24" x 36")
• Material: Blue silicone
• New material with greater lubricity!
• 1" OD cross section
Fits: 24" x 36" Door
Models: MC3600




CSG032
CROSS SECTION
NOT SHOWN ACTUAL SIZE

Warning: These gaskets are not for use in Ethylene Oxide (E.O. Gas) sterilizers