
# OHMEDA • GI



# **New parts to fit**

# Giraffe Infant Incubators and OmniBed!





### RPI Part #0MK028


OEM Part #6600-1406-500 (Inside Latch Cover only)/6600-1403-500 (Side Wall Latch only)/6600-1405-500 (Latch Spring only)/6600-1404-500 (Latch Cover only)

### SIDE WALL LATCH KIT

• Includes all parts as shown

Fits: Side Wall

Models: Giraffe OmniBed & Giraffe Incubator





### RPI Part #0MK029

OEM Part #6600-1239-500 (Hinge Bottom only)/6600-1242-500 (Hinge Cover only)/6600-1240-500 (Middle Hinge only)/6600-1041-400 (Hinge Pin only)/6600-1040-400 (Left Spring only)/6600-1026-400 (Right Spring only)

### **PORTHOLE HINGE KIT**

• Includes all parts as shown

Fits: Side Panels

Models: Giraffe OmniBed & Giraffe Incubator

### RPI Part #0MK027 0EM Part #6600-0738-700

### **PORTHOLE LATCH ASSY KIT**

• Includes all parts as shown **Fits:** Side and Front Panels

Models: Giraffe OmniBed & Giraffe Incubator





# RPI Part #0MG017 0EM Part #6600-1248-500

### **CORNER GROMMET**

Material: Translucent thermoplastic rubber
 Fits: End Panel

Models: Giraffe OmniBed & Giraffe Incubator






RPI Part #0MP021 0EM Part #6600-1007-700

### **TOUCH PANEL ASSEMBLY**

Size: ~7.12" x 13"For use with <u>OmniBed</u>

Fits: Display Module Models: Giraffe OmniBed For Use With OmniBed!



RPI Part #0MP022 0EM Part #6600-1254-600

### **TOUCH PANEL ASSEMBLY**

• Size: ~7.12" x 13"

• For use with <u>Incubator</u>

Fits: Display Module
Models: Giraffe Incubator



Incubator!



RPI Part #0MS025 0EM Part #6600-0736-701

### **CANOPY LIFT SWITCH (LEFT)**

Fits: Upper Chassis Models: Giraffe OmniBed





RPI Part #0MS026 0EM Part #6600-0736-702

**CANOPY LIFT SWITCH (RIGHT)** 

Fits: Upper Chassis Models: Giraffe OmniBed



# replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

SECTION 7


PAGE **212** 

UPDATED 01/14

### **Infant Incubators**

### Model: Giraffe® Incubator & Giraffe® OmniBed®








### **BREAK PLATE**

- · Includes all parts as shown
- Included in: Compartment Air Probe (RPI Part #OMP018)

Fits: Compartment Air Probe

Models: Giraffe OmniBed & Giraffe Incubator





### RPI Part #0MP018 0EM Part #6600-1512-700

### **COMPARTMENT AIR PROBE**

· Includes all parts as shown

· Also available: Humidity Sensor (RPI Part #0MS038)

Fits: Rear Wall Panel

Models: Giraffe OmniBed & Giraffe Incubator

### RPI Part #0MP037 **0EM Part #6600-1474-500**

### **BACK PLATE**

- Includes all parts as shown
- Included in: Air Compartment Probe (RPI Part #0MP018)

Fits: Compartment Air Probe

Models: Giraffe OmniBed & Giraffe Incubator

RPI Part #0MS038 **0EM Part #6600-0784-700** (Sensor Only)/6600-0540-400 (0-ring Only)/6600-1475-500 (Retainer Only)

### **HUMIDITY SENSOR**

- Includes all parts as shown
- Also available: Calibration Kit (RPI Part #0MK040) and Compartment Air Probe (RPI Part #0MP018)

Fits: Compartment Air Probe

Models: Giraffe OmniBed & Giraffe Incubator





### RPI Part #0MK040 OEM Part #6600-0048-850

### **CALIBRATION KIT**

- Use to calibrate Humidity Sensor (RPI Part #0MS038)
- · Solution usable for up to 1 year after mixing
- Includes: Vial with sodium chloride and rubber elbow

Fits: Humidity Sensor

Models: Giraffe OmniBed & Giraffe Incubator



# RPI Part #LMP027 OEM Part #6600-0912-200

### **LAMP**


- 12V, 50W
- Base: 2-Pin
- Type: Halogen; EXT

Models: Giraffe Exam Light, Giraffe OmniBed



RETAINER

**RP0304** 



### RPI Part #0MS023 **OEM Part #6600-1154-600** (Switch Only)/ 6600-1165-400 (Locking Nut Only)

### **DOOR SWITCH**

- SPDT
- 1A@125 VAC
- 1/8" Tabs
- · Includes all parts as shown

Fits: Heater Doors Models: Giraffe OmniBed





### RPI Part #0MS034 OEM Part #6600-1014-602

### **POWER SW / CIRCUIT BREAKER**

- Rocker switch (on-off)
- 12A @ 125/250 VAC; DPST / Trip amps 16.2
- Black switch body with white legend: I/O
- 1/4" Tabs

Fits: Lower Electrical Enclosure


Models: Giraffe OmniBed & Giraffe Incubator



- Fits: Exam Light

& Giraffe Incubator





### RPI Part #0MB030

**OEM Part #6600-1372-500** (Lift Belt Only)/ 6600-1071-400 (Belt Tension Spring Only)/ 6600-1187-400 (Grommet Only)



### **LIFT BELT**

- Length: 75"; Pitch: 0.1475"; Width: 0.23"
- Material: Molded polyurethane with twin stainless steel cores
- Includes: Belt Tension Spring (RPI Part #OMS031) and heater door opening rod grommet


Fits: Lift Rail Belt Channel Models: Giraffe OmniBed



### **BELT TENSIONER SPRING**

- Size: ~.54"0D x .75" lg. x .06" wire
- Material: Stainless steel
- Included in: Lift Belt (RPI Part #OMB030)

Fits: Lift Belt/Lift Rail Channel Models: Giraffe OmniBed



OEM Part #6600-0718-700/ 6600-1187-400 (Grommet Only)



### **HOOD HARNESS**

• Includes all parts as shown Fits: Lift Rail Cable Channel Models: Giraffe OmniBed

### RPI Part #0MS039

**OEM Part #6600-0718-700 (Complete Hood** harness)/6600-1187-400 (Grommet Only)

### **MICROSWITCH**

- SPDT; 0.1 A@48VAC
- Includes: Mounting hardware, shrink tubing. and heater door opening rod grommet
- Included in: Hood Harness (RPI Part #0MH032)
- Replace the complete Hood Harness (RPI Part #OMH032) or just replace the Microswitches, you decide - it's the RPI advantage!

Fits: Hood Harness/Lift Rail Channel

Models: Giraffe OmniBed




### RPI Part #0MS019 **OEM Part #6600-0820-700 (Sensor** Only)/6600-0709-405 (Screw Only)

### **AIR FLOW SENSOR**

- IR Emitter / Detector pair
- Includes all parts as shown

Fits: Lower Chassis

Models: Giraffe OmniBed & Giraffe Incubator




### RPI Part #0MR033 **OEM Part #6600-1003-600**

### **SOLID STATE RELAY**

- Size: ~1.2" x 1" x .9"; 1/4" Tabs
- SPST-N.O.; 12 A@24-280 VAC
- Control voltage 3-30 VDC;
- · Includes all parts as shown

Fits: Electrical Enclosure

Models: Giraffe OmniBed & Giraffe Incubator






### RPI Part #0MF024 0EM Part #6600-1523-700

### **COOLING FAN**

- Size: 1.575" x 1.575" x .41" fan housing
- 6.0-13.8 VDC; 5800 RPM
- · Includes all parts as shown

Fits: Heat Sink Vent

Models: Giraffe OmniBed & Giraffe Incubator



## RPI Part #0MB035 OEM Part #6600-0562-603

### **CIRCUIT BREAKER (3.5A)**

- 3.5 Amps
- 1/4" Tabs


RPC464

Includes all parts as shown

Fits: Connection Panel

Models: Giraffe OmniBed & Giraffe Incubator















PRESORTED STANDARD U.S. POSTAGE PAID VAN NUYS, CA **PERMIT NO. 2399** 

INSTALLATION

ADDRESS SERVICES REQUESTED

# Don't Forget ... RPI also carries these great Parts to fit Giraffe Incubators & OmniBed!

### RPI Part #0MH005 0EM Part #6600-0220-850

### **HEATING ASSEMBLY**

- 460W, 99VAC, ~21Ω
- Includes Connector and **Ground Wire Assembly**

Fits: Hood Assembly

Models: Giraffe OmniBed & Giraffe Incubator

### RPI Part #0MF015 OEM Part #6600-0207-850

### MICRO FILTER

- 10 per package
- ~1-1/2" x 5-1/2"
- Synthetic media contains no glass fibers
- Latex free

Fits: Filter Housing

Models: Giraffe OmniBed & Giraffe Incubator

### RPI Part #0MP009 OEM Part #6600-0873-700

### SKIN TEMPERATURE PROBE (DISPOSABLE)

- · Single use only
- 3 pin connector
- Includes: (10) probes and (10 pieces) Reflective Probe Cover (RPI Part #RPC619)

Fits: Probe Panel

Models: Giraffe OmniBed & Giraffe Incubator



RPI Part #0MP010 0EM Part #6600-0874-700

SKIN TEMPERATURE PROBE (DISPOSABLE) - Case

 Includes: (50) probes and (50 pieces) Reflective Probe Cover (RPI Part #RPC619)



• Also available: Reflective Probe Covers (RPI Part #RPC619)

Fits: Probe Panel

Models: Giraffe OmniBed & Giraffe Incubator



RPI Part #RPC619 OEM Part #0203-1980-300

### **REFLECTIVE PROBE COVER**

- 100 per package; Single use only
- Metalized Mylar with PSA backing; Latex Free

Fits: Skin Temperature Probe

**Models:** All Infant Warmer and Incubator Skin Temperature Probes