



Reflectors and More to fit A-DEC Cascade 6300 Dental Lights!

RPI Part #ADR185
OEM Part # (No OEM Part # available)
REFLECTOR with Free Lens Splash Shield
 • New glass - not re-coated
 • Includes Reflector, and Lens Splash Shield
 (RPI Part #ADL126)
Fits: Light head
Model: All 6300 models

**FREE
LENS SPLASH SHIELD**




RPI Part #ADS129
OEM Part #041-336-00
POWER SWITCH
 • Includes all parts as shown
 • SPST
 • 20 A @ 125 VAC;
 10 A @ 250 VAC
Fits: Head Assembly and Swivel Block
Model: All 6300 models



RPI Part #ADL126
OEM Part #28-0503-01
LENS SPLASH SHIELD
 • Clear plastic
 • Mounting hardware not included
Fits: Light head
Model: All 6300 models



RPI Part #ADS125
OEM Part #28-0633-00
LAMP SOCKET ASSEMBLY
 • Includes all parts as shown
 • Uses replacement Lamp
 (RPI Part #LMP600)
Fits: Inside of Focus Tube
Model: All 6300 models



RPI Part #LMP600
New OEM Part #041-179-01
Old OEM Part #041-179-00
LAMP
 • 24 V; 150 W
 • Base: 2-pin
Fits: Light socket
Model: All 6300 models



	RPI PART#	OEM PART #	DESCRIPTION
NEW!	ADL126	28-0503-01	LENS SPLASH SHIELD
NEW!	ADR185	(No OEM Part #)	REFLECTOR WITH FREE LENS SPLASH SHIELD
	ADS125	28-0633-00	LAMP SOCKET ASSEMBLY
	ADS129	041-336-00	POWER SWITCH
	ADS130	041-177-00	INTENSITY SWITCH
	ADT127	28-0566-00	LOCKING TAB (PART NOT SHOWN ABOVE)
	ADT132	28-0724-00	TRANSFORMER (PART NOT SHOWN ABOVE)
	ADT133	28-0822-00	TRANSFORMER (PART NOT SHOWN ABOVE)
	ADT134	28-0720-00	TRANSFORMER (PART NOT SHOWN ABOVE)
	LMP600	041-179-01	LAMP
	RPH243	046-123-00	CIRCUIT BREAKER (PART NOT SHOWN ABOVE)