

Models: Belmont X-Calibur and Healthco X-Calibur

NEW PARTS TO FIT BELMONT & HEALTHCO DENTAL LIGHTS


RPI Part #BEH030
OEM Part #HLU 1130

WIRING HARNESS IN BALANCE ARM

- Includes connector
- 6 ft. long

Fits: Balance arm assembly

Models: Belmont X-Calibur 046/BLW & 047/BLC;
Healthco X-Calibur 046/HLW & 047/HLC



RPI Part #BEH031
OEM Part #HLU 0350

WIRE HARNESS IN BALANCE ARM

- Includes connector
- 10 ft. long

Fits: Balance arm assembly

Models: Belmont X-Calibur HLPW;
Healthco X-Calibur BLPW



RPI Part #BEH032
OEM Part #HLU 1140

WIRE HARNESS IN TRACK MOUNT

- Includes connector
- 8 ft. long

Fits: Track mount BLT

Models: Belmont X-Calibur 048/BLT;
Healthco X-Calibur 048/HLT




RPI Part #BES023
OEM Part #8000-1001

POWER SWITCH

- Includes all parts as shown
- SPDT (On-Off-On); Nylon toggle
- 10 A @ 250 VAC; 15 A @ 125 VAC
- 3/4 hp, 125-250 VAC

Fits: Head assembly

Models: All Belmont X-Calibur models;
All Healthco X-Calibur models




RPI Part #BES034
New OEM Part #8000-02005
Old OEM Part #800-02005

FRONT SHIELD

- Clear Lexan
- 6.250" Dia. x 4.870" Flats x .06" thk.

Fits: All oval handles

Models: All Belmont X-Calibur with oval handles



RPI Part #BES033
OEM Part #HLU 2005

FRONT SHIELD

- Clear Lexan
- 6.36" x 5.38" x .06" thk

Fits: All square handles

Models: All Belmont X-Calibur with square handles;
All Healthco X-Calibur with square handles




RPI Part #LMP601
Industry Part #H-3

LAMP

- 12 VDC; 55 W
- Base:Wire; Halogen

Fits: Light socket

Models: All Belmont X-Calibur models;
All Healthco X-Calibur models



RPI Part #BES002
OEM Part #HLU 2054

RUBBER SHEEL

- 2 per package; Buna-N
- 5/32" dia. x 4-5/8" lg. with angled end

Fits: Yolk assembly

Models: All Belmont X-Calibur models;
All Healthco X-Calibur models




RPI Part #BES006
OEM Part #HLU 2004

SCREW (FRAME BACK COVER)

- Includes all parts as shown
- Stainless
- Special thread: #8-36 UNF

Fits: Back frame cover

Models: All Belmont X-Calibur models;
All Healthco X-Calibur models




NEW PARTS TO FIT BELMONT & HEALTHCO DENTAL LIGHTS

RPI Part #BEK028
OEM Part # (No OEM Part # Available)

**WASHER KIT
(JOINT BOLT ASSEMBLY)**


• Includes all parts as shown
Fits: Joint bolt assembly
Models: All Belmont X-Calibur models;
All Healthco X-Calibur models



RPI Part #BEA001
OEM Part #8000-2023

SPRING CLIP HOLDER ASSEMBLY


Fits: Light Socket
Models: All Belmont X-Calibur models;
All Healthco X-Calibur models



RPI Part #BEH013
OEM Part #8000-2027

SOCKET HOLDER


• Includes all parts as shown
• Brass
Fits: Light head
Models: All Belmont X-Calibur models;
All Healthco X-Calibur models



RPI Part #BES012
OEM Part #HLU 2025

LIGHT SOCKET

• Includes all parts as shown
Fits: Light Head Assembly
Models: All Belmont X-Calibur models;
All Healthco X-Calibur models



RPI Part #BEW010
New OEM Part #HLU 2026

**GROUNDING WIRE
ASSEMBLY**

Fits: Light socket assembly
Models: All Belmont X-Calibur models;
All Healthco X-Calibur models



RPI Part #BEW016
New OEM Part #HLU 2011

RUBBER WASHER


• 2 per package
• 1.210" ID x 1.780" OD x .087" thk
• Clear silicone
Fits: Both sides of reflector assembly
Models: All Belmont X-Calibur models;
All Healthco X-Calibur models



RPI Part #BEN014
OEM Part # (No OEM Part # Available)

NUT (SOCKET HOLDER)

• Anodized black
Fits: Light socket holder
Models: All Belmont X-Calibur models;
All Healthco X-Calibur models



**RPI PARTS ARE
GUARANTEED
TO YOUR 100%
COMPLETE SATISFACTION!**