

Models: Belmont X-Calibur and Healthco X-Calibur

**NEW PARTS TO FIT BELMONT & HEALTHCO DENTAL LIGHTS**


RPI Part #BEH030  
OEM Part #HLU 1130

**WIRING HARNESS IN BALANCE ARM**

- Includes connector
- 6 ft. long

Fits: Balance arm assembly

Models: Belmont X-Calibur 046/BLW & 047/BLC;  
Healthco X-Calibur 046/HLW & 047/HLC



RPI Part #BEH031  
OEM Part #HLU 0350

**WIRE HARNESS IN BALANCE ARM**

- Includes connector
- 10 ft. long

Fits: Balance arm assembly

Models: Belmont X-Calibur HLPW;  
Healthco X-Calibur BLPW



RPI Part #BEH032  
OEM Part #HLU 1140

**WIRE HARNESS IN TRACK MOUNT**

- Includes connector
- 8 ft. long

Fits: Track mount BLT

Models: Belmont X-Calibur 048/BLT;  
Healthco X-Calibur 048/HLT




RPI Part #BES023  
OEM Part #8000-1001

**POWER SWITCH**

- Includes all parts as shown
- SPDT (On-Off-On); Nylon toggle
- 10 A @ 250 VAC; 15 A @ 125 VAC
- 3/4 hp, 125-250 VAC

Fits: Head assembly

Models: All Belmont X-Calibur models;  
All Healthco X-Calibur models




RPI Part #BES034  
New OEM Part #8000-02005  
Old OEM Part #800-02005

**FRONT SHIELD**

- Clear Lexan
- 6.250" Dia. x 4.870" Flats x .06" thk.

Fits: All oval handles

Models: All Belmont X-Calibur with oval handles




RPI Part #BES033  
OEM Part #HLU 2005

**FRONT SHIELD**

- Clear Lexan
- 6.36" x 5.38" x .06" thk

Fits: All square handles

Models: All Belmont X-Calibur with square handles;  
All Healthco X-Calibur with square handles




RPI Part #LMP601  
Industry Part #H-3

**LAMP**

- 12 VDC; 55 W
- Base:Wire; Halogen

Fits: Light socket

Models: All Belmont X-Calibur models;  
All Healthco X-Calibur models



RPI Part #BES002  
OEM Part #HLU 2054

**RUBBER SHEEL**

- 2 per package; Buna-N
- 5/32" dia. x 4-5/8" lg. with angled end

Fits: Yolk assembly

Models: All Belmont X-Calibur models;  
All Healthco X-Calibur models




RPI Part #BES006  
OEM Part #HLU 2004

**SCREW (FRAME BACK COVER)**

- Includes all parts as shown
- Stainless
- Special thread: #8-36 UNF

Fits: Back frame cover

Models: All Belmont X-Calibur models;  
All Healthco X-Calibur models



**DENTAL LIGHTS - BELMONT & HEALTHCO**

SECTION  
4

PAGE  
1202


UPDATED  
06/04

**NEW PARTS TO FIT BELMONT & HEALTHCO DENTAL LIGHTS**

RPI Part #BEK028  
OEM Part # (No OEM Part # Available)

**WASHER KIT  
(JOINT BOLT ASSEMBLY)**


• Includes all parts as shown  
Fits: Joint bolt assembly  
Models: All Belmont X-Calibur models;  
All Healthco X-Calibur models



RPI Part #BEA001  
OEM Part #8000-2023

**SPRING CLIP HOLDER ASSEMBLY**


Fits: Light Socket  
Models: All Belmont X-Calibur models;  
All Healthco X-Calibur models



RPI Part #BEH013  
OEM Part #8000-2027

**SOCKET HOLDER**


• Includes all parts as shown  
• Brass  
Fits: Light head  
Models: All Belmont X-Calibur models;  
All Healthco X-Calibur models



RPI Part #BES012  
OEM Part #HLU 2025

**LIGHT SOCKET**

• Includes all parts as shown  
Fits: Light Head Assembly  
Models: All Belmont X-Calibur models;  
All Healthco X-Calibur models



RPI Part #BEW010  
New OEM Part #HLU 2026

**GROUNDING WIRE  
ASSEMBLY**

Fits: Light socket assembly  
Models: All Belmont X-Calibur models;  
All Healthco X-Calibur models



RPI Part #BEW016  
New OEM Part #HLU 2011

**RUBBER WASHER**

• 2 per package  
• 1.210" ID x 1.780" OD x .087" thk  
• Clear silicone  
Fits: Both sides of reflector assembly  
Models: All Belmont X-Calibur models;  
All Healthco X-Calibur models



RPI Part #BEN014  
OEM Part # (No OEM Part # Available)

**NUT (SOCKET HOLDER)**

• Anodized black  
Fits: Light socket holder  
Models: All Belmont X-Calibur models;  
All Healthco X-Calibur models



**RPI PARTS ARE  
GUARANTEED  
TO YOUR 100%  
COMPLETE SATISFACTION!**