

Models: Chairman, Chairman 5000, Chairman 5010 and Coachman

Capacitors, Switches & Fuse To Fit Pelton & Crane Chairman & Coachman Chairs!

RPI Part #PCC728
OEM Part #94-416-43

LIFT CAPACITOR (56-60 MFD)

- Plastic; Electrolytic medium

Fits: Base motor/lift motor

Models: Chairman 5000 & Coachman




RPI Part #PCS716
OEM Part #007429

TRAVERSE SWITCH

- Momentary contact (On/Off/On)
- 6A @ 125 V; 3A @ 250 V
- Snap in mounting with low profile paddle (dark grey paddle)

Fits: Back assembly and foot switch assembly

Model: Chairman




RPI Part #PCS717
OEM Part #007428

TILT SWITCH

- Momentary contact (On/Off/On)
- 6A @ 125 V; 3A @ 250 V
- Snap in mounting with low profile paddle (white paddle)

Fits: Back assembly

Model: Chairman



RPI Part #PCC734
OEM Part #77-31-904/019518

CAPACITOR (21-25 MFD)

- Plastic, Electrolytic medium

Fits: Incline/recline and upper traverse

Models: Chairman 5000 and Coachman




RPI Part #PCS718
OEM Part #007427

LIFT SWITCH

- Momentary contact (On/Off/On)
- 6A @ 125 V; 3A @ 250 V
- Snap in mounting with low profile paddle (black paddle)

Fits: Back assembly

Model: Chairman




RPI Part #PCS713
OEM Part #007431

ROLLER MICRO SWITCH

- .250" (6.38mm) roller lever
- 5A @ 250 VAC

Fits: Seat wiring ass

Model: Chairman




RPI Part #PCS724
OEM Part #018420

MICRO SWITCH

- 5A @ 250 VAC

Fits: Chair back asser

Model: Coachman




RPI Part #PCS731
OEM Part #15-25-737

CORRECTION LIMIT SWITCH

- 10A @ 125 VAC

Fits: Lift motor

Model: Chairman 5000



RPI Part #RPF367
OEM Part #18-81-163

FUSE (2A, 250V)

- 5 per package
- Size: 5 mm x 20 mm (.197" x .796")
- MDL Glass, Fast-Acting

Fits: Inner base assembly 240 V only

Models: Chairman 5000 & 5010

