



INTRODUCING THE RPI FIELD SERVICE CALIBRATION KIT TO FIT MIDMARK•RITTER M9 & M11

RPI Part #MIK074


OEM Part # (No OEM Part # Available)

FIELD SERVICE CALIBRATION KIT

- To check and adjust pressure and temperature
- Includes all parts as shown
- All parts in kit are sold separately - see below

Models:

- M9 (S/N Series: CZ, DA, DB, DX, DY, FD & OM)
- M9D (S/N Series: FF & FG)
- M11 (S/N Series: ES, ET, FP, FR & GB)
- M11D (S/N Series: GC & GD)



RPT480 (x25)

RPT113

RPT460

Carrying Case

*See RPI Part #RPC476 below

All Parts in Kit Also Sold Separately – See Below!


RPI Part #RPG461

OEM Part #002-0372-00

TEST PRESSURE GAUGE

- Includes Gauge, Tubing (RPI Part #RPT462) and Kwik Clamps (RPI Part #RPC463)
- Dual gauge display: 0-60 psi/0-400 KPA
- Gauge has rubber cover protector
- Can be used on most sterilizers for in-line pressure check

Models: M9 & M11




RPI Part #RPT460

OEM Part # (No OEM Part # Available)

TRIM POT ADJUSTMENT TOOL

- 5" Long
- Plastic with stainless steel blades
- Includes shirt pocket clip

Fits: Most trim pots




RPI Part #RPT480

OEM Part #015-0013-05 (1/pkg)

CABLE TIE (HIGH TEMP)

- 25 per package
- Blue
- 1/8" x 4"
- Temperature rating: Up to 338°F




RPI Part #RPC463

OEM Part # (No OEM Part # Available)

KWIK™ CLAMPS

- 6 per package
- Temperature range: -40°F to 302°F
- Nylon®
- Fits:** 3/8" OD tubing




RPI Part #RPT113

OEM Part # (No OEM Part # Available)

MAX-REGISTER THERMOMETER

- Mercury
- Length: 5-1/2"
- Accuracy Rating: ±1 Division (±2°F)
- Range: 150°F to 350°F
- Includes Thermometer and metal protective case




RPI Part #RPC476

OEM Part # (No OEM Part # Available)

CARRYING CASE

- Plastic case with foam inserts*
- Inside dimensions: 9.5"(L) x 6"(W) x 2.69 (D); Outside dimensions: 10.5"(L) x 7"(W) x 3"(D)
- *Note: When carrying case is purchased separately, the "pick-n-pluck" foam inserts can be customized by you to create your own storage design and a blank case label is included for custom description to fit your needs.




RPI Part #RPT462

OEM Part # (No OEM Part # Available)

TUBING (SILBRADE)

- Sold per foot
- 1/8" ID x 3/8" OD
- Silicone tubing with reinforced braid
- Max temp 350°F
- Working Pressure: 233 psi @ 70°F
- Fits:** Barbed Fittings (1/8" diameter)



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpi.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Midmark•Ritter equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Website 06/05