





INTRODUCING NEW PARTS TO FIT GENDEX GXP™ FILM PROCESSORS




RPI Part #GXF012
OEM Part #652-0012P1
HEATER FAN (115V)
 • 4.69" square x 1.5" wide
Fits: Heater/Fan Assembly




RPI Part #GXF011
OEM Part #122-0096G1
COOLING FAN (115V)
 • 3.15" square x 1" wide
Fits: Chassis



RPI Part #GXB009
OEM Part #420-0007P1
DRIVE BELT
 • 190XL025 Power Grip Timing Belt
Fits: Developer, fixer, wash and dryer racks




RPI Part #GXS014
OEM Part #632-0045P1
PUSH BUTTON SWITCH (RED)
 • Endo switch
 • Includes Switch Cover (RPI Part #GXC015)
Fits: Switch panel




RPI Part #GXS013
OEM Part #632-0046P1
PUSH BUTTON SWITCH (WHITE)
 • Run/Stand-by and On/Off switches
 • Includes Switch Cover (RPI Part #GXC015)
Fits: Switch panel




RPI Part #GXC015
OEM Part #433-0060P1
SWITCH COVER
 • Transparent
Fits: All push button switches and lamp



RPI Part #GXB004
OEM Part #652-0015P1/Brush
MOTOR BRUSH (2/pkg)
 • Cap sold separately (RPI Part #RPC038)
Fits: Brush holder assembly on main drive motor



RPI Part #RPH638
OEM Part # 646-0047P1
FUSE HOLDER (MAIN)
 • 3/16" terminals
 • Fits standard 1/4" x 1-1/4" fuse
 • Includes mounting hardware
Fits: Rear panel



RPI Part #GXF024
OEM Part #112-0629G1
LEVELING FOOT
 • 4 per package
 • Non-skid PVC base
Fits: Bottom frame



RPI Part #RPC038
OEM Part #652-0015P1/Cap
BRUSH CAP (6/pkg)
 • Brush sold separately (RPI Part #GXB004)
Fits: Main drive motor




RPI Part #RPF070
OEM Part #646-0048P1
FUSE (MAIN)
 • 12A, 250V; 1/4" x 1-1/4"
 • Ceramic
Fits: Rear panel

PC BOARD FUSES: F1 TO F4

RPI Part #	OEM Part #	Position & Description
RPF045	640-0054P16	F1; 1A, 250V
RPF046	46-1700-21P34	F2; 1-1/2A, 250V
RPF053	46-1700-21P33	F3; 8A, 250V
RPF512	No OEM Part # Available	F4; 1A, 250V

TO HELP FIX THOSE ANNOYING LEAKS!

RPI Part #RPS639
OEM Part # (No OEM # Available)
SILICONE GASKET MAKER 598
 • 5 per package
 • .05 fl. oz. per unit
 • Provides a sealant for pump tubing connections



RPI Part #RPC637
OEM Part # 407-0159P1
KWIK™ CLAMP (6/pkg)
 • Range: 23/64"-25/64"
 • Use with pump tubing



RPI Part #RPO360
OEM Part #404-0036P2
O-RING (6/pkg)
 • .301" x .441" x .070"
 • Material: EPDM
Fits: Standpipe



LOOK FOR MORE PARTS TO FIT GENDEX FILM PROCESSORS AND X-RAY EQUIPMENT LATER THIS SUMMER!



replacement parts industries, inc.


CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@ripiparts.com • WEBSITE www.ripiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Gendex equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Flyer No Prices 06/07



RPI introduces a new line ...
parts to fit
GENDEX
GXP™ FILM PROCESSORS


...and many more to come later this summer!



**CUSTOMIZE YOUR WORK SPACE
 WITH RPI POWER CORDS!**

*Depending on where your outlet is located ...
 choose the power cord that best fits your needs!*

Right Angle Plug




RPI Part #RPC582
OEM Part #644-0009P1

INDUSTRIAL GRADE POWER CORD

- 14/3 SJT, 105°C; 15A@125V
- NEMA 6-15P; EN60320-C13R Right Angle

Straight Plug



RPI Part #RPC262


OEM Part # (No OEM # Available)

INDUSTRIAL GRADE POWER CORD

- 14/3 SJT, 60°C; 15A@125V
- NEMA 5-15P; EN60320-C13

**NEED THE
 POWER MODULE?**

*Check the power module.
 Is it time to replace it?*



RPI Part #HBM011
OEM Part #651-0068P1

POWER MODULE

- 15A@250V