

Table Top Sterilizers

Models: UltraClave® M9, M9D, M11 & M11D

MIDMARK® • RITTER®



**JUST IN TIME
FOR THE HOLIDAYS!**

- New Dam Gasket and New PM Kit ... to fit M9 & M9D door with metal flap!
- New and Improved Dam Gasket ... to fit M9 & M9D door without metal flap!
- New and Improved Dam Gasket ... to fit all M11 & M11D!

All parts in Kit also sold separately!

RPI Part #MIK158
OEM Part #002-0361-00
(Door & Dam Gasket only)

STERILIZER PM KIT

- Includes all parts as shown plus a PM Check List & PM Sticker
- Also available: Sterilizer PM Kit (RPI Part #MIK072) for door without metal flap

Fits: Door **with** metal flap
(OEM Part #030-0698-00)
Models: M9 & M9D

Kit fits M9 & M9D
(door with metal flap)

MIG157
Dam Gasket
• Only fits Door with metal flap

MIF062
Mesh Chamber Filter
• For newer style machines

MIF038
Sintered Bronze Filter
• For older style machines

MIG027
Door Gasket

MIS076
Door Springs

All parts in Kit also sold separately!

RPI Part #MIK072
OEM Part #002-0361-01
(Door & Dam Gasket only)

STERILIZER PM KIT

- Includes all parts as shown plus a PM Check List & PM Sticker
- Also available: Sterilizer PM Kit (RPI Part #MIK158) for door with metal flap

Fits: Door **without** metal flap
(OEM Part #030-1016-00)
Models: M9 & M9D

Kit fits M9 & M9D
(door without metal flap)

MIG035
Dam Gasket
• Only fits Door without metal flap

MIF062
Mesh Chamber Filter
• For newer style machines

MIF038
Sintered Bronze Filter
• For older style machines

MIG027
Door Gasket

MIS076
Door Springs

NEW New Design For Easy Install & Better Performance

All parts in Kit also sold separately!

RPI Part #MIK080
OEM Part #002-0504-00
(Door & Dam Gasket only)

STERILIZER PM KIT

- Includes all parts as shown plus a PM Check List & PM Sticker

Models: M11 & M11D

Kit fits all M11 & M11D

MIG036
Dam Gasket

MIF062
Mesh Chamber Filter
• For newer style machines

MIF038
Sintered Bronze Filter
• For older style machines

MIG028
Door Gasket

MIS076
Door Springs

NEW New Design For Easy Install & Better Performance



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Midmark®/Ritter® equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Show 12/11 updated 02/12


SECTION
10
PAGE
1018
UPDATED
12/11

JUST IN TIME FOR THE HOLIDAYS!



To fit M9 & M9D door
with metal flap ...

NEW DAM GASKET & NEW PM KIT!



To fit M9 & M9D door without metal flap
and all M11 & M11D ...

NEW & IMPROVED DAM GASKETS!

- Easier to install!
- Designed to channel water back into the chamber!



MIG035



MIG036

Just look at the difference ...
See for yourself!

Side view of
OLD gasket

Side view of
NEW gasket



replacement parts industries, inc.

P.O. Box 5019 • Chatsworth, California • 91313-5019

PRESORTED STANDARD
U.S. POSTAGE PAID
VAN NUYS, CA
PERMIT #2399

• Address Service Requested •



An ISO 9001 Certified Company



LIST OF PM KITS AND PARTS TO FIT MIDMARK® • RITTER® M9, M9D, M11 & M11D ULTRACLAVE® STERILIZERS

RPI PART #	OEM PART #	DESCRIPTION	MODELS
MIF038	002-0360-00 (Sintered Style)	FILTER (SINTERED STYLE)	Fits M9, M9D, M11 & M11D
MIF062	002-0360-00 (Mesh Style)	FILL/VENT MESH CHAMBER FILTER - 2/PKG	Fits M9, M9D, M11 & M11D
MIG027	053-0366-00 (Door Gasket only)	DOOR GASKET	Fits M9 & M9D
MIG028	053-0527-00 (Door Gasket only)	DOOR GASKET	Fits M11 & M11D
MIG035	053-0903-00 (Dam Gasket only)	DAM GASKET	Fits M9 & M9D with Door without metal flap
MIG036	053-0784-00 (Dam Gasket only)	DAM GASKET	Fits M11 & M11D
MIG157	053-0508-00 (Dam Gasket only)	DAM GASKET	Fits M9 & M9D with Door with metal flap
MIK072	002-0361-01 (Door & Dam Gaskets only)	STERILIZER PM KIT	Fits M9 & M9D with Door without metal flap
MIK080	002-0504-00 (Door & Dam Gaskets only)	STERILIZER PM KIT	Fits M11 & M11D
MIK158	002-0361-00 (Door & Dam Gaskets only)	STERILIZER PM KIT	Fits M9 & M9D with Door with metal flap
MIS076	050-2347-00	DOOR SPRINGS - 2/PKG	Fits M9, M9D, M11 & M11D