

FILTER, CHAMBER SEALS & PM KITS


TO FIT STERIS SYSTEM



RPI PART #SSF012 OEM Part #A1525

A1525 STERILE WATER FILTER


- 0.2 micron
- Silicone O-ring Model: System 1



RPI PART #SSS001 **OEM Part #200213**

INFLATABLE CHAMBER SEAL


- Includes chamber seal, stem, coupler, plug, and adapter Model: System 1



RPI PART #SSS037 OEM Part #200253

INFLATABLE CHAMBER SEAL

- EPDM
- · Includes chamber seal and stem Model: System 1






RPI PART #SSK038 OEM Part #200687

MAJOR PM KIT #1

- · Includes:
- (1) Check Valve (CK-8 Drain Check) Retro Fit Kit (RPI Part #SSK008),
- (1) Check Valve CK-1 (RPI Part #SSV002),
- (2) Check Valves CK-2 & CK-3 (RPI Part #SSV003), and
- (2) Pinch Valve Sleeves (RPI Part #SSS007)

Model: System 1



RPI PART #SSK039 OEM Part #200975

MAJOR PM KIT #2

- · Includes:

 - (1) Check Valve CK-1 (RPI Part #SSV002), (2) Check Valves CK-2 & CK-3 (RPI Part #SSV003), and
 - (2) Pinch Valve Sleeves (RPI Part #SSS007)

Model: System 1

