




# FILTER, CHAMBER SEALS & PM KITS TO FIT STERIS SYSTEM 1

**LOOK**  
They're Back!

**RPI PART #SSF012**  
**OEM Part #A1525**

### A1525 STERILE WATER FILTER


- 0.2 micron
  - Silicone O-ring
- Model:** System 1



**RPI PART #SSS001**  
**OEM Part #200213**

### INFLATABLE CHAMBER SEAL


- EPDM
  - Includes chamber seal, stem, coupler, plug, and adapter
- Model:** System 1



**RPI PART #SSS037**  
**OEM Part #200253**

### INFLATABLE CHAMBER SEAL

- EPDM
  - Includes chamber seal and stem
- Model:** System 1




**LOOK**  
What's New!

**RPI PART #SSK038**  
**OEM Part #200687**

### MAJOR PM KIT #1

- Includes:
    - (1) Check Valve (CK-8 Drain Check) Retro Fit Kit (RPI Part #SSK008),
    - (1) Check Valve CK-1 (RPI Part #SSV002),
    - (2) Check Valves CK-2 & CK-3 (RPI Part #SSV003), and
    - (2) Pinch Valve Sleeves (RPI Part #SSS007)
- Model:** System 1



**RPI PART #SSK039**  
**OEM Part #200975**

### MAJOR PM KIT #2

- Includes:
    - (1) Check Valve CK-1 (RPI Part #SSV002),
    - (2) Check Valves CK-2 & CK-3 (RPI Part #SSV003), and
    - (2) Pinch Valve Sleeves (RPI Part #SSS007)
- Model:** System 1

