



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com




# JUST IN TIME FOR THE HOLIDAYS ...

PC Board to fit New Style Midmark® UltraClave® M9 & M11 now available!



**Fits  
New Style  
M9 & M11**



RPI Part #MIB129  
OEM Part #022-0762-00


**CONTROL PC BOARD (REFURBISHED)**

- Includes: Mounting Bracket
- When returning the old PC Board, **YOU MUST** include the Mounting Bracket.

**Fits:** Chassis

**Models:** M9 -020 & -022 Serial # prefixes RN, RR & V; and M11 -020 & -022 Serial # prefixes RS, RV & V only

**Fits  
Old Style  
M9**



RPI Part #MIB130  
OEM Part #002-0434-00/002-0434-03


**CONTROL PC BOARD (REFURBISHED)**

- Mounting Bracket not included
- When returning the old PC Board, **DO NOT** include the Mounting Bracket. Keep the old Mounting Bracket for later use.

**Fits:** Chassis

**Models:** M9 -001 thru -019 with Serial # prefixes CZ, DA, DB, DX, DY, FD & OM

**Fits  
Old Style  
M11**



RPI Part #MIB131  
OEM Part #002-0501-00

**CONTROL PC BOARD (REFURBISHED)**

- Mounting Bracket not included
- When returning the old PC Board, **DO NOT** include the Mounting Bracket. Keep the old Mounting Bracket for later use.

**Fits:** Chassis

**Models:** M11 -001 thru -019 with Serial # prefixes ES, ET, FP, FR & GB


**Here's How the \$100 Credit Exchange Program Works.**

- 1) Place an order for a Control PC Board (PCB), and RPI will send a shipping box and prepaid UPS return label to you. You will not pay for shipping when using this box.
- 2) Using the shipping box and UPS return label, ship as follows:
  - Include the PCB and Mounting Bracket if shipping the *New Style* M9/M11 PCB (OEM Part #022-0762-00).
  - Do not include the Mounting Bracket if shipping the *Old Style* M9/M11 PCB (OEM Part #'s 002-0434-00/002-0434-03/002-0501-00).
- 3) Once the non-operational board has been received by RPI and verified that it is in reasonable condition without burn holes, a replacement PCB will be shipped to you. Your invoice will reflect the cost of the board plus the \$100 credit for the Core.

**Important Information.** PC Board performance verified and tested and Quality Control test printout is included with each PC Board. 90 Day Warranty.

**Warning.** To avoid damaging the replacement PC Board and possibly voiding the warranty, check the condition of the *Heater and/or Pulse Solenoid* prior to installing the replacement PC Board. DO NOT use PC Board in 220V or version "D" sterilizers.

**Fits  
All Styles  
M9 & M11**



RPI Part #MIK159  
OEM Part # (No OEM Part # Available)

**PCB MOUNTING KIT**


- Includes all parts as shown
- Fits:** PC Board to Mounting Bracket and Chassis
- Models:** M9 & M11



RPI Part #MIK117  
OEM Part # (No OEM Part # Available)

**PRESSURE TUBE KIT**

- Includes all parts as shown
- Translucent silicone tubing
- Models:** M9 & M11



**For more great parts to fit Midmark® UltraClave® M9 & M11 Autoclaves, check www.rpiparts.com**