



... Proud to be a Healthcare Components Group brand ...




HEALTHCARE  
Components Group

MIDMARK • RITTER

We've Got You Covered ... New Parts to fit

# MIDMARK • RITTER ULTRACLAVE® M9 & M11 STERILIZERS!




**RPI Part #MIK267**  
**OEM Part #002-0782-00/002-0782-01/  
002-0782-02/002-0782-03**

### TOP COVER KIT

- Includes: Touch Pad (RPI Part #MIP099), Pressure Relief Handle (RPI Part #MIH269), Cover Plate (Printer) (RPI Part #MIC270) and (6 pcs) Screw (#10 x 1-1/4" Phillips Pan Head) (RPI Part #RPH851)
- Also available: Steam Block Kit (RPI Part #MIK236), Display Assembly (RPI Part #MIA147), Hinge Cover (RPI Part #MIC271) and Speednut (RPI Part #PCN136)

**Fits:** Top of Sterilizer

**Models:** M9-020 thru -022 and M9D-020 thru -022




**RPI Part #MIK268**  
**OEM Part # 002-0780-00/002-0780-01**

### TOP COVER KIT

- Includes: Touch Pad (RPI Part #MIP099), Pressure Relief Handle (RPI Part #MIH269), Cover Plate (Printer) (RPI Part #MIC270) and (6 pcs) Screw (#10 x 1-1/4" Phillips Pan Head) (RPI Part #RPH851)
- Also available: Steam Block Kit (RPI Part #MIK237), Display Assembly (RPI Part #MIA147), Hinge Cover (RPI Part #MIC271) and Speednut (RPI Part #PCN136)

**Fits:** Top of Sterilizer

**Models:** M11-020 thru -022 & M11D-020 thru -022




**RPI Part #MIK265**  
**OEM Part #002-0783-00/  
002-0783-01/  
053-1251-00**

### DOOR PANEL KIT

- Includes: Door Handle (RPI Part #MIH163), Hinge Cover (RPI Part #MIC271) and (4 pcs) Screw (#8 x 7/8" Phillips Pan Head) (RPI Part #RPH855) - *It's the RPI Advantage!*
- Also available: Steam Block Kit (RPI Part #MIK236)

**Fits:** Door Frame

**Models:** M9-020 thru -022 & M9D-020 thru -022



**RPI Part #MIK266**  
**OEM Part #002-0781-00/  
002-0781-01/  
053-1255-00**

### DOOR PANEL KIT

- Includes: Door Handle (RPI Part #MIH163), Hinge Cover (RPI Part #MIC271) and (4 pcs) Screw (#8 x 7/8" Phillips Pan Head) (RPI Part #RPH855) - *It's the RPI Advantage!*
- Also available: Steam Block Kit (RPI Part #MIK237)

**Fits:** Door Frame

**Models:** M11-020 thru -022 & M11D-020 thru -022



replacement parts industries, inc. Add this page to your RPI catalog – Section 10, Page 1026


CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@ripiparts.com](mailto:order@ripiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)

The above parts are manufactured by Replacement Parts Industries, Inc. to fit the equipment noted. All product names used in this document are trademarks or registered trademarks of their respective holders. Flyer Code 1703MI RD2 Web 08/22

SECTION  
10

PAGE  
1026

UPDATED  
08/22



**RPI Part #MIH269**  
**OEM Part #053-1252-00**

**PRESSURE RELIEF HANDLE**

- Includes: (1 pc) Screw (#10 x 1-1/4" Phillips Pan Head) (RPI Part #RPH851)
- Included in: Top Cover Kit (RPI Part #'s MIK267 & MIK268)
- Also available: Safety Valve (40 PSI) (RPI Part #'s MIV122 & MIV090)

**Fits:** Top Cover to Safety Valve  
**Models:** M9-020 thru -022, M9D-020 thru -022, M11-020 thru -022 & M11D-020 thru -022




**RPI Part #MIC271**  
**OEM Part #053-1277-00**

**HINGE COVER**

- Included in: Door Panel Kit (RPI Part #'s MIK265 & MIK266)

**Fits:** Hinge Assembly  
**Models:** M9-020 thru -022, M9D-020 thru -022, M11-020 thru -022 & M11D-020 thru -022



**RPI Part #MIC270**  
**OEM Part #053-0933-00**

**COVER PLATE (PRINTER)**

- Included in: Top Cover Kit (RPI Part #'s MIK267 & MIK268)

**Fits:** Top Cover  
**Models:** M9-020 thru -022, M9D-020 thru -022, M11-020 thru -022 & M11D-020 thru -022


# Make Your Midmark •Ritter® UltraClave® Look as Good as New with RPI's Decals!

**RPI Part #MIK329**  
**OEM Part #002-10231-00/061-0311-00**

**DECAL KIT**

- Ideal for bringing an old sterilizer back to life.
- Includes: (2 pcs) Decal (Hot Surface) (RPI Part #MID330), Decal (Distilled Water) (RPI Part #MID331), Decal (Water Level) (RPI Part #MID332), and Decal (Caution) (RPI Part #MID335)

**Fits:** Front and Rear Panels  
**Models:** M9, M9D, M11, M11D



**RPI Part #MID330**  
**OEM Part #061-0726-00**

**DECAL (HOT SURFACE)**

- 2 per package
- Size: 1.50" x 1.56"
- Included in: Decal Kit (RPI Part #MIK329)

**Fits:** Front Panel and Inside of Door Panel  
**Models:** M9, M9D, M11 & M11D




**RPI Part #MID331**  
**OEM Part #061-0827-00**

**DECAL (DISTILLED WATER)**

- Size: 2.00" x 2.75"
- Included in: Decal Kit (RPI Part #MIK329)

**Fits:** Front Panel  
**Models:** M9, M9D, M11 & M11D



**RPI Part #MID332**  
**OEM Part #061-0311-00**

**DECAL (WATER LEVEL)**

- Size: 1" x 4.25"
- Included in: Decal Kit (RPI Part #MIK329)

**Fits:** Front Panel Behind the Drain Hose  
**Models:** M9, M9D, M11 & M11D





**RPI Part #MID335**  
**OEM Part #MPM60455**

**DECAL (CAUTION)**

- Size: .75" x 1.50"
- Included in: Decal Kit (RPI Part #MIK329)

**Fits:** Rear Panel  
**Models:** M9, M9D, M11 & M11D






**RPI Part #MIH327**  
**OEM Part #015-1639-00**

**HEATER WIRE HARNESS**

- Size: #14 AWG, 14.00" lg. wire, 1/4" QC terminals, with wire marker #11 and #12
- Also available: Jumper Wire (RPI Part #MIJ328) and Main Harness (RPI Part #MIH238)

**Fits:** PC Board to Heater (Newer Style Units Only)

**Models:** M9-020 thru -022, M9D-020 thru -022, M11-020 thru -022 & M11D-020 thru -022



**RPI Part #MIJ328**  
**OEM Part #015-1640-00**

**JUMPER WIRE**

- Size: #14 AWG, 2.85" lg. wire, 1/4" QC flag terminals, with wire marker #14
  - Included in: Main Harness (RPI Part #MIH238)
  - Also available: Heater Wire Harness (RPI Part #MIH327)
- Fits:** Hi Limit Thermostat to Thermostat (Newer Style Units Only)
- Models:** M9-020 thru -022, M9D-020 thru -022, M11-020 thru -022 & M11D-020 thru -022

**The Grommet is now included in the Steam Trap Kit!**



**RPI Part #RXG099**  
**OEM Part #053-0394-02**

**GROMMET**

- 6 per package
- Size: .250" ID x .500" OD x .250" Thk.; Groove Size: .375" dia. x .125" wide
- Material: Black SBR Rubber, 60 Duro

**Fits:** Water Reservoir to 1/4" Copper Tubing, Groove Size: .375" dia. x .125" wide (Older Style Units Only)

**Models:** M9-001 thru -019, M9D-001 thru -019, M11-001 thru -019 & M11D-001 thru -019



**RPI Part #MIK078**  
**OEM Part #002-0654-00**

**STEAM TRAP KIT (New Style/Upgrade Kit)**

- Includes: Steam Trap Valve (RPI Part #MIV081), (1 pc) Grommet (RPI Part #RXG099), (3 pcs) Compression Nut (RPI Part #RCN059), and copper tubing
- Fits:** Chamber to Condensating Tank (Older Style Units Only)
- Models:** M9-001 thru -019, M9D-001 thru -019, M11-001 thru -019 & M11D-001 thru -019

**Drain Hose Kits to fit M9 & M11 Sterilizers with two different Hose Clip options!**

*Check Serial Numbers for which Kit you need!*

**Fits Newer Models**



**RPI Part #MIC333**  
**OEM Part #015-1695-00**


**HOSE CLIP**

- 2 per package
- Size: 11/16" diameter capacity, 1/2" wide, black plastic
- Fits mounting hole size: .252" thru .267" dia.
- Included in: Drain Hose Kit (RPI Part #MIK334)

**Fits:** Reservoir Drain Tubing to Front Panel (Newer Style Units Only)

**Models:** M9/M9D (-020 thru -022 Serial #'s V1464210-Present) & M11/M11D (-020 thru -022 Serial #'s V1479341-Present)

**Fits Newer Models**



*Includes Push-In Mount Hose Clips!*

**RPI Part #MIK334**  
**OEM Part # (No OEM Part # Available)**

**DRAIN HOSE KIT**

- 13" lg. drain hose can be cut to fit as needed.
- Includes: (2 pcs) Hose Clip (RPI Part #MIC333), (1 pc) Kwik Clamp (RPI Part #RPC286), (13" Long) Tubing (3/8" ID x 5/8" OD Clear) (RPI Part #RPH285)
- Also available: Drain Hose Kit (RPI Part #MIK061) for older style units

**Fits:** Reservoir to Front Panel with Hole Diameter Size: .252" thru .267" for the "Revised Hose Clip Replacement"

**Models:** M9/M9D (-020 thru -022 Serial #'s V1464210-Present) & M11/M11D (-020 thru -022 Serial #'s V1479341-Present)

**Fits Older Models**



*Includes Rivot Mount Hose Clips!*

**RPI Part #MIK061**  
**OEM Part # (No OEM Part # Available)**

**DRAIN HOSE KIT**

- 13" lg. drain hose can be cut to fit as needed.
- Includes: (2 pcs) Hose Clip (RPI Part #RPC283), (2 pcs) Pop Rivet (RPI Part #RPR284), (1 pc) Kwik Clamp (RPI Part #RPC286), (13" Long) Tubing (3/8" ID x 5/8" OD Clear) (RPI Part #RPH285)
- Also available: Drain Hose Kit (RPI Part #MIK334) for newer style units

**Fits:** Reservoir to Front Panel with Hole Diameter Size: .122" thru .126" (Older Style Units Only)

**Models:** M9/M9D (-001 thru -019), M11/M11D (-001 thru -019), M9/M9D (-020 thru -022 Serial #'s V2200-V1464209) & M11/M11D (-020 thru -022 Serial #'s V2200-V1479340)