



More New Parts to fit Midmark • Ritter 7 & M7 SpeedClave® Sterilizers




INDEPENDENT BUZZER-STYLE TIMER


- Use this Kit to replace or upgrade a Timer with an independent buzzer.
(Use *Time Switch (RPI Part #RCT048)* if you need a timer with a built-in bell.)
- Timer Knob sold separately – see RPI Part #'s RCK111 and RCK092.
- The Wire Harness in this Kit may be needed to extend the length of the original Wire Harness when upgrading to this style timer.

RPI Part #RCK116
OEM Part #002-0373-00

TIMER KIT (30 MIN)
Fits: Front Panel
Model: M7



Includes all parts as shown





RPI Part #RCK111
OEM Part #053-1008-01

TIMER KNOB
Fits: Timer
Model: M7

RPI Part #RCK092
OEM Part #H98062

TIMER KNOB
Fits: Timer
Models: 7 & M7




RPI Part #RCT107
OEM Part # (No OEM Part # Available)

TIMER (30 MIN)
• See "Independent Buzzer-Style" note above with Timer Kit (RPI Part #RCK116)
Fits: Front Panel
Model: M7



RPI Part #RCB096
New OEM Part #015-1813-00
Old OEM Part #015-0745-02

BUZZER
Fits: Machine Base
Model: M7



RPI Part #RCH117
OEM Part # (No OEM Part # Available)


WIRE HARNESS
• See "Independent Buzzer-Style" note above with Timer Kit (RPI Part #RCK116)
Fits: Original Wire Harness
Model: M7

BUILT-IN BELL-STYLE TIMER

- Use this Timer to replace a Timer with a built-in bell. This timer will not work in conjunction with machines that utilize an independent buzzer. (Use *Timer Kit (RPI Part #RCK116)* if you need a timer with an independent buzzer.)

RPI Part #RCT048
OEM Part #H284667-1

TIME SWITCH
Fits: Front Panel
Models: 7 & M7-without buzzer




Includes all parts as shown



RPI Part #RPK213
OEM Part # (No OEM Part # Available)


KNOB
Fits: Time Switch (RPI Part #RCT048) only
Models: 7 & M7



RPI Part #RCA099
OEM Part #H98064


FLEX CABLE ASSEMBLY
Fits: Front Panel and Relay Assembly
Models: 7 & M7

Includes all parts as shown



RPI Part #RCK110
OEM Part #053-1008-00

TEMPERATURE KNOB
Fits: Flex Cable Assembly
Model: M7



RPI Part #RCK093
OEM Part #H98063

CONTROL KNOB
Fits: Flex Cable Assembly
Models: 7 & M7



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Midmark-Ritter equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. WEB flyer with newsletter 04/09




More New Parts to fit Midmark•Ritter 7 & M7 SpeedClave® Sterilizers

RPI Part #RCH118
Old RPI Part #RCH024
OEM Part #002-0839-00, H284687-1
& 002-0245-00

HEATER ELEMENT ASSEMBLY


- Includes all parts as shown
- See **Service Tip**
- Fits:** Chamber
- Models:** 7 & M7

Attaching Hardware (RPI Part #RCH119)
 for the Heater Element Assembly (RPI Part #RCH024) is also sold separately. Includes (2) each Fiber Washer, Flat Washer, Internal Tooth Washer and Hex Nut



Service Tip

To help prevent clogging of the 7 or M7, install RPI Filter (RPI Part #MIF062) into the chamber end of the Water Fill/Exhaust line.




RPI Part #MIF062
New OEM Part #002-0360-00 (Mesh Style)
Old OEM Part #016-0790-00

FILL/VENT MESH CHAMBER FILTER

- See **Service Tip**
- Fits:** Chamber
- Models:** 7 & M7


RPI Part #RCK106
NEW OEM Part #002-0248-00
OLD OEM Part #053-0513-00

DOOR HANDLE KIT
Fits: Door Assembly
Models: 7 & M7



RPI Part #RCP120
OEM Part #042-0005-16

CAM PIN
Fits: Door Handle Assembly
Models: 7 & M7




RPI Part #RCN103
OEM Part #K906034

PUSH CAP NUT (5/PKG)
Fits: Cam Pin (RPI Part #RCP120)
Models: 7 & M7




RPI Part #MIT047
New OEM Part #002-0370-00
Old OEM Part #015-0647-00

OVERTEMP THERMOSTAT
Fits: Bracket that mounts to under side of chamber
Model: M7




RPI Part #RCT050
OEM Part #H225639

MANUAL THERMOSTAT
Fits: Bracket that mounts to under side of chamber
Models: 7 & M7



RPI Part #RCL121
OEM Part #015-1304-00

PILOT LIGHT (GREEN)
Fits: Front Panel
Model: M7



Don't forget, RPI also carries these parts to fit Midmark•Ritter 7 & M7 sterilizers and more!

RPI PART #	OEM PART #	DESCRIPTION	RPI PART #	OEM PART #	DESCRIPTION
PCC003	H98159	INDUSTRIAL GRADE POWER CORD (7 & M7)	RCS012	040-0010-103, H96673 & 3046	HINGE SCREW (7 & M7)
RCB002	H284709	BELLOWS ASSEMBLY (7 Only)	RCS042	002-0263-00 & 33053	RELAY SWITCH ASSEMBLY (7 & M7)
RCB089	002-0375-00	AIR VENT BELLOWS ASSEMBLY (7 & M7) (Old OEM Part #014-0193-00)	RCS086	H225730	SLEEVE (7 & M7)
RCB100	002-0375-00	AIR VENT BELLOWS ASSEMBLY (7 & M7) (Old OEM Part #014-0193-00)	RCS098	002-0247-00	VALVE STEM KIT (7 & M7)
RCC010	34646	CAM (7 Only)	RCS108	H106330	DOOR SPRING (7 & M7)
RCC095	H98009	FILL COVER (7 & M7)	RCS115	015-1791-00	FILL/VENT SWITCH (M7 Only)
RCC114	014-0236-10	COIL (FILL/VENT) (M7 Only)	RCV051	H98089 & 11683	SAFETY VALVE (31 PSI) (7 & M7)
RCF016	H98021	SUCTION CUP FOOT (7 & M7)	RCV112	014-0236-00	SOLENOID VALVE (FILL/VENT) (M7 Only)
RCG019	H98134 & 33166	DOOR GASKET (7 & M7)	RCW009	H96771 & 3059	DOOR WASHER (7 & M7)
RCG085	002-0242-00	TEMPERATURE GAUGE (DIAL THERMOMETER) (7 & M7)	RPB167	H98160	STRAIN RELIEF BUSHING (7 & M7)
RCG087	H98135	RESERVOIR GASKET (7 & M7)	RPF426	053-0990-00	FEMALE QUICK CONNECT FITTING (GREY) (M7 Only)
RCG104	053-1579-00	FIBER GASKET (7 & M7)	RPF427	053-0990-01	MALE QUICK CONNECT FITTING (GREY) (M7 Only)
RCK113	014-0236-11	REPAIR KIT (FILL/VENT) (M7 Only)	RPF679	015-0346-35	FUSE (12A, 250VAC) - FAST ACTING (M7)
RCL090	H98036	PILOT LIGHT (RED) (7 & M7)	RPH674	015-1259-04	FUSE HOLDER (M7 Only)
RCN013	041-0010-08 & 3047	HINGE SCREW NUT (7 & M7)	RPK431	002-0653-00	DRAIN TUBE ASSEMBLY KIT (M7)
RCP011	042-0007-07 & H97989	CAM PIN (7 & M7)	RPO306	H98136 & 36906	O-RING (7 & M7)

See the RPI Catalog or website for more parts to fit the 7 & M7 sterilizers!

Note: Parts fit model(s) noted to the right of the descriptions above.