


TELLE


New Parts to fit Midmark • Ritter Tables and Chairs


WARNING: When attaching hydraulic hoses to the cylinders, torque the hose fittings between 55 - 60 Inch Pounds only. Additional torque could damage the threads in the cylinder.
















POWER HOSE Fits: Pump to Base Models: 106 & 306







Fits: Pump to Base Models: 106 & 306





* See Chart on Back for Models and Serial Number Breakdowns







replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

SECTION

PAGE **304**

RPI PART #	OEM PART #	DECRIPTION	MODELS with Prefix and Serial # Breakdown Note: Part fits all prefixes and all serial #'s for models listed if not otherwise noted.
MIC089	002-0094-00	BASE CYLINDER	105-001 thru -006; 107; 305-001 thru -004; 307-001 thru -007
MIC111	002-0100-00	BACK CYLINDER	105-001 thru -006; 305-001 thru -004
MIC112	002-0137-00	TILT CYLINDER	112-001 ; 113-001 ; 117-001 thru -003 serial #'s AJ1000 thru AJ3143
MIC113	002-0271-00	TILT CYLINDER	117-001 thru -003 serial #'s AJ3144 & FE1000 thru present; 317-001
MIC114	002-0346-00	BACK CYLINDER	117-001 thru -003 serial #'s FE1000 thru present; 119-001 thru -016 serial #'s BY2006, CA1191, DF2695, EM1000, EN1000, EY1000, EZ1000 & HK1000 thru present; 311-005; 317-001; 318-001; 319-003 thru -006 with prefixes HN, KK, KL & V
MIH100	002-0144-00 002-0054-00	POWER HOSE	106-001 thru -004; 306-001 thru -003
MIH101	002-0143-00 002-0053-00	RETURN HOSE	106-001 thru -004; 306-001 thru -003
MIH102	002-0095-00	RETURN HOSE	107; 307-001 thru -007
MIH103	002-0096-00	POWER HOSE	107; 307-001 thru -007
MIH104	002-0004-00	POWER HOSE	110-001
MIH105	002-0109-00	POWER HOSE	105-001 thru -002 serial #'s G1000 thru G1324 and H1000 thru H1019
MIH106	002-0005-00	RETURN HOSE	110-001
MIH107	002-0098-00	RETURN HOSE	105-001 thru -006 serial #'s G1325, H1020, V1000, CF1000, CG1000 & CA1000 thru present; 305-001 thru -004
MIH108	002-0099-00	POWER HOSE	105-001 thru -006 serial #'s G1325, H1020, V1000, CF1000, CG1000 & CA1000 thru present; 305-001 thru -004
MIH109	002-0110-00	RETURN HOSE	105-001 thru -002 serial #'s G1000 thru G1324 and H1000 thru H1019



replacement parts industries, inc.

P.O. Box 5019 • Chatsworth, California • 91313-5019

PRESORTED STANDARD U.S. POSTAGE PAID VAN NUYS, CA PERMIT #2399

• Address Service Requested •

An ISO 9001 Certified Company

