





**INTRODUCING A NEW LINE OF PARTS TO FIT
GE/CRITIKON DINAMAP™ 8100 & 8100T NIBP MONITORS!**




PUMP ASSEMBLY
RPI Part #GEP004
OEM Part #320413

Fits: Chassis
Models: 8100, 8100T,
9300 & 9340



CHECK VALVE
RPI Part #RPV642
OEM Part #754116


Fits: Pump Assembly
Models: 8100, 8100T,
9300 & 9340



TUBING
RPI Part #RPT645
OEM Part #740107


- 1/8" ID x 1/4" OD
- Opaque; Silicone blend
- Sold by the foot (Reel not included)

Fits: Chassis
Models: 8100, 8100T,
9300 & 9340



FRONT PANEL OVERLAY
RPI Part #GEK002
OEM Part #701252


Fits: Front Panel
Model: 8100



Includes
Temperature
Probe Overlay

FRONT PANEL OVERLAY
RPI Part #GEK001
OEM Part #300918


Fits: Front Panel
Model: 8100T



BATTERY
RPI Part #GEB006
OEM Part #320378

- 12VDC, 2.8AMP Hr.

Fits: Chassis
Model: 8100, 8100T, 8110, 8120,
9300, 9340, 9350, 9400 & 9500XL




BATTERY CABLE
RPI Part #GEC013
OEM Part #316341

Fits: Power Supply PCB
Model: 8100, 8100T, 9300 & 9340



CALIBRATION PLUG
RPI Part #GEP014
OEM Part #320438

Fits: Temp Probe Connector
Model: 8100T




QUICK CONNECT
RPI Part #GEC012
OEM Part #754505

Fits: Dual Tubing Set
Models: 8100 & 8100T



O-RING (12/pkg)
RPI Part #RPO313
OEM Part #300593


Fits: Quick Connect
Models: All models using dual
tubing set



Board Mounted

TEMP PROBE CONNECTOR
RPI Part #GEC008
OEM Part #300681


Fits: Display PCB
Model: 8100T



Wired to Board

TEMP PROBE CONNECTOR
RPI Part #GEC009
OEM Part #316339


Fits: Case Bracket
Model: 8100T



POWER INPUT MODULE
RPI Part #GEM015
OEM Part #320431

- Includes fuse holder and (2 pieces) 1/2A@250VAC Fuse (RPI Part #RPF043)

Fits: Chassis
Models: 8100 & 8100T



HOSPITAL GRADE POWER CORDS
RPI Part #RPC557 - 8' long

RPI Part #RPC558 - 18' long
RPI Part #RPC559 - 25' long