



The Motor Plus More New Parts to fit PELTON & CRANE CHAIRMAN DENTAL CHAIRS!




RPI Part #PCG739
OEM Part #022701

TRAVERSE GEAR (SINGLE)

- 17 Teeth
- Includes Grease (RPI Part #RPL717)

Fits: Traverse Motor
Model: Chairman




RPI Part #PCG762
OEM Part #024770

LIFT GEAR

- 39 Teeth, Helical
- Includes Grease (RPI Part #RPL717)

Fits: Lift Motor
Model: Chairman




RPI Part #PCG740
OEM Part #024766

TRAVERSE GEAR (Double)

- 16 Teeth/35 Teeth
- Includes Grease (RPI Part #RPL717)

Fits: Traverse Motor
Model: Chairman




RPI Part #PCM745
OEM Part #017628 Serial #1001-4189
/007881 Serial #4190 and above

LIFT MOTOR

- Includes Lift Motor Pulley Kit (RPI Part #PCK792) and (3 pieces) Pig Tail Connector (RPI Part #RPT013)

Fits: Base and Lift Assembly
Model: Chairman



RPI Part #RPL717
OEM Part # (No OEM Part # Available)

GREASE

- 30cc per package
- Renolit G 193 Grease
- MSDS and Technical Data Sheet available on RPI Website

Fits: Lift Motor
Model: Chairman

List of RPI Parts to fit Pelton & Crane Chairman Dental Chairs

| RPI Part # | OEM Part # | Description | Notes |
|---------------|-------------------------|-------------------------|-------|
| PCA747 | 007660 | AUTO RECLINE ACTUATOR | |
| PCB729 | 017850 | PCB ASSEMBLY | |
| PCC746 | 007903 | COILING WIRE ASSEMBLY | |
| PCG739 | 022701 | TRAVERSE GEAR (SINGLE) | |
| PCG740 | 024766 | TRAVERSE GEAR (DOUBLE) | |
| PCG762 | 024770 | LIFT GEAR | |
| PCK738 | 012476 | GEAR AND RACK | |
| PCK742 | 011195 | BRAKE REPAIR KIT | |
| PCK792 | 007575 | LIFT MOTOR PULLEY KIT | |
| PCM745 | 017628* & 007881** | LIFT MOTOR | |
| PCP721 | 007595 | FOOT SWITCH PEDAL | |
| PCS713 | 007431 | ROLLER MICRO SWITCH | |
| PCS714 | 007506 | FORMED LEVER SWITCH | |
| PCS715 | 007596 | TRAVERSE LIMIT SWITCH | |
| PCS716 | 007429 | TRAVERSE SWITCH | |
| PCS717 | 007428 | TILT SWITCH | |
| PCS718 | 007427 | LIFT SWITCH | |
| PCS719 | 007898 | OPERATE SWITCH ASSEMBLY | |
| PCS720 | 007899 | EXIT SWITCH ASSEMBLY | |
| RPL717 | No OEM Part # Available | GREASE | |

*017628 Serial #1001-4189 **007881 Serial #4190 and above



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Pelton & Crane equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Web 05/10