

# 13 WAYS TO BE PREPARED WHEN SERVICING PELTON & CRANE STERILIZERS

Introducing thirteen sterilizer PM kits exclusively from RPI. Don't simply put in a new gasket when you can replace the bellows, filter and other critical PM parts all at the same time with the RPI Sterilizer PM Kits. What's more, the kits are less expensive than purchasing the parts individually!

## 1

### RPI Sterilizer PM Kit (RPI Part #PCK216)

#### Fits Magna-Clave®

- Gasket (RPI Part #PCG019)
- Bellows (RPI Part #PCB001)
- Air Valve O-ring (RPI Part #RP0354)
- Filter Fill Line (RPI Part #PCF009)
- Filter with Cotter Pin (RPI Part #PCK226)

## 2

### FILTER with COTTER PIN (RPI Part #PCK226)

#### Fits Magna-Clave®

- Includes filter and cotter pin

## 3

### RPI Sterilizer PM Kit (RPI Part #PCK217)

#### Fits Sentry

- Gasket (RPI Part #PCG020)
- Bellows (RPI Part #PCB001)
- Air Valve O-ring (RPI Part #RP0354)
- Filter Fill Line (RPI Part #PCF009)

## 4

### RPI Sterilizer PM Kit (RPI Part #PCK218)

#### Fits OCM

- Gasket (RPI Part #PCG014)
- Bellows (RPI Part #PCB001)
- Air Valve O-ring (RPI Part #RP0354)
- Filter Fill Line (RPI Part #PCF009)

## 5

### RPI Sterilizer PM Kit (RPI Part #PCK219)

#### Fits OCR & OCR+

- Gasket (RPI Part #PCG015)
- Bellows (RPI Part #PCB001)
- Air Valve O-ring (RPI Part #RP0354)
- Filter Fill Line (RPI Part #PCF009)

## 6

### RPI Sterilizer PM Kit (RPI Part #PCK220)

#### Fits Validator® 8

- Gasket (RPI Part #PCG065)
- Bellows (RPI Part #PCB001)
- Air Valve O-ring (RPI Part #RP0354)
- Filter Fill Line (RPI Part #PCF009)
- Door Spacer (RPI Part #PCS067)

## 7

### RPI Sterilizer PM Kit (RPI Part #PCK221)

#### Fits Validator® 10

- Gasket (RPI Part #PCG066)
- Bellows (RPI Part #PCB001)
- Air Valve O-ring (RPI Part #RP0354)
- Filter Fill Line (RPI Part #PCF009)
- Door Spacer (RPI Part #PCS067)

## 8

### RPI Sterilizer PM Kit (RPI Part #PCK222)

#### Fits Validator® 8 Plus

- Gasket (RPI Part #PCG065)
- Filter Fill Line (RPI Part #PCF009)
- Reservoir Filter Kit (RPI Part #PCK167)
- Door Spacer (RPI Part #PCS067)

## 9

### RPI Sterilizer PM Kit (RPI Part #PCK223)

#### Fits Validator® 10 Plus

- Gasket (RPI Part #PCG066)
- Filter Fill Line (RPI Part #PCF009)
- Reservoir Filter Kit (RPI Part #PCK167)
- Door Spacer (RPI Part #PCS067)

## 10

### RPI Bellows Kit (RPI Part #PCK224)

#### Fits OCM, OCR, OCR+, Sentry, Validator® 8 & Validator® 10

- Bellows Cap (RPI Part #PCC005)
- Bellows (RPI Part #PCB001)
- Air Valve O-ring (RPI Part #RP0354)

## 11

### RPI Filter & Tube Kit (RPI Part #PCK225)

#### Fits Delta 8, Delta 10, OCM, OCR, OCR+, Sentry & all Validator® series

- Filter Fill Line (RPI Part #PCF009)
- Fill Chamber Tube (RPI Part #PCT010)

## 12

### RPI Sterilizer PM Kit (RPI Part #PCK231)

#### Fits Delta 8


- Gasket (RPI Part #PCG065)
- Hepa Filter (RPI Part #RPF644)
- Compressor Filter (RPI Part #SCF002)
- Reservoir Filter Kit (RPI Part #PCK167)
- Filter Fill Line (RPI Part #PCF009)
- Door Spacer (RPI Part #PCS067)
- Filter In-Line Assembly (RPI Part #PCF233)

## 13

### RPI Sterilizer PM Kit (RPI Part #PCK232)

#### Fits Delta 10

- Gasket (RPI Part #PCG066)
- Hepa Filter (RPI Part #RPF644)
- Compressor Filter (RPI Part #SCF002)
- Reservoir Filter Kit (RPI Part #PCK167)
- Filter Fill Line (RPI Part #PCF009)
- Door Spacer (RPI Part #PCS067)
- Filter In-Line Assembly (RPI Part #PCF233)



**replacement parts industries, inc.**

CALL (800) 221-9723 • FAX (818) 882-7028 • E-MAIL [order@rpiparts.com](mailto:order@rpiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)