
# RPI IS YOUR BEST SOURCE FOR

# RACKS AND TRAYS TO FIT TABLE TOP STERILIZERS!



# MIDMARK® • RITTER® UltraClave® M9 (M9-001 thru -019 & M9-020 thru -022)



#### RPI Part #MIK204

OEM Part #030-0710-01 (Tray Rack), 050-3691-00 (Tray Rack), 002-0253-00 (Instrument Tray (Small)), 002-0374-00 (Instrument Tray (Large)) & 050-1773-00 (Tray Plate)

#### **RACK AND TRAY KIT**

- Includes all parts as shown; All parts also sold separately!
- Material: Stainless steel

Fits: Chamber

RPI Part #MIR208

0EM Part #030-0710-01 & 050-3691-00

TRAY RACK

Accommodates old and new style Instrument Trays!

RPI Part #MIT206

**OEM Part #002-0253-00** 

**INSTRUMENT TRAY (SMALL)** 

Use with wire or sheet metal style racks!

• Size: 6" x 12"

Also fits: SpeedClave® M7 (M7-001 thru M7-022)

RPI Part #MIT207 OEM Part #002-0374-00

#### **INSTRUMENT TRAY (LARGE)**

• Size: 7" x 12"

Use with wire or sheet metal style racks!

RPI Part #MIP205

OEM Part #050-1773-00 (Old style) &


050-3750-00 (New style)

TRAY PLATE

Use with wire or sheet metal style racks!

# Fits

# MIDMARK® • RITTER® UltraClave® M11 (M11-001 thru -019 & M11-020 thru -022)



#### RPI Part #MIK209

OEM Part #050-3920-00 (Tray Rack), 050-2644-00 (Instrument Tray (Small)), 050-4260-00 (Instrument Tray (Small)), 050-4259-00 (Instrument Tray (Large)), 050-2643-00 (Instrument Tray (Large)) & 050-3919-00 (Tray Plate)

#### **RACK AND TRAY KIT**

- Includes all parts as shown; All parts also sold separately!
- Material: Stainless steel

Fits: Chamber

RPI Part #MIR213 0EM Part #030-0887-00 & 050-3920-00 (Tray Rack)

**TRAY RACK** 

Accommodates old and new style Instrument Trays!

RPI Part #MIT211 0EM Part #050-2644-00 & 050-4260-00

**INSTRUMENT TRAY (SMALL)** 

• Size: 7" x 15"

Use with wire or sheet metal style racks!

RPI Part #MIT212 0EM Part #050-2643-00 & 050-4259-00

INSTRUMENT TRAY (LARGE)

• Size: 9" x 15"

Use with wire or sheet metal style racks!

RPI Part #MIP210
0EM Part #050-2600-00 (Old Style)
& 050-3919-00 (New Style)

TRAY PLATE

Use with wire or sheet metal style racks!



# replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

## PELTON & CRANE Delta AF10, Q10 & XL10


## RPI Part #PCT142 0EM Part #1539357/004429

#### **INSTRUMENT TRAY (SMALL)**

- Size: 6-1/2" x 17" x 1"
- · Material: Stainless steel
- Tray can be used with larger tray (RPI Part #PCT143)

Fits: Chamber

Models: AF10, Q10 & XL10




### RPI Part #PCT143 0EM Part #1539340/044159/004430

#### **INSTRUMENT TRAY (LARGE)**

- Size: 8" x 17" x 2-1/2"
- · Material: Stainless steel
- Tray can be used with smaller tray (RPI Part #PCT142)

Fits: Chamber

Models: AF10, Q10 & XL10



## Fits

## PELTON & CRANE OCM & Sentry


#### RPI Part #PCT140 OEM Part #004040

#### **INSTRUMENT TRAY (SMALL)**

- Size: 5" x 13" x 1"
- Material: Stainless steel
- Tray can be used with larger tray (RPI Part #PCT141)

Fits: Chamber

Models: OCM & Sentry




#### RPI Part #PCT141 OEM Part #004141

#### **INSTRUMENT TRAY (LARGE)**

- Size: 6-1/2" x 13" x 1"
- · Material: Stainless steel
- Tray can be used with smaller tray (RPI Part #PCT140)

Fits: Chamber Model: OCM




#### RPI Part #PCR188 OEM Part #004286

# TRAY REST & SUPPORT ASSEMBLY

- Size: 4.5" x 13"
- · Material: Stainless steel

Fits: Chamber Model: OCM



## PELTON & CRANE OCR & OCR+


#### RPI Part #PCT142 0EM Part #1539357/004429

#### **INSTRUMENT TRAY (SMALL)**

- Size: 6-1/2" x 17" x 1"
- · Material: Stainless steel
- Tray can be used with larger tray (RPI Part #PCT143)

Fits: Chamber

Models: OCR & OCR+




#### RPI Part #PCT143 OEM Part #1539340/044159/004430

#### **INSTRUMENT TRAY (LARGE)**

- Size: 8" x 17" x 2-1/2"
- Material: Stainless steel
- Tray can be used with smaller tray (RPI Part #PCT142)

Fits: Chamber Models: OCR & OCR+




#### RPI Part #PCR189 OEM Part #004426

# TRAY REST & SUPPORT ASSEMBLY

- Size: 6" x 17"
- Material: Stainless steel

Fits: Chamber

Models: OCR & OCR+



## PELTON & CRANE Validator 10 & Plus 10


#### RPI Part #PCT142 0EM Part #1539357/004429

#### **INSTRUMENT TRAY (SMALL)**

- Size: 6-1/2" x 17" x 1"
- Material: Stainless steel
- Tray can be used with larger tray (RPI Part #PCT143)

Fits: Chamber

Model: Validator 10 (AA)



#### RPI Part #PCT143 OEM Part #1539340/044159/004430


#### **INSTRUMENT TRAY (LARGE)**

- Size: 8" x 17" x 2-1/2"
- Material: Stainless steel
- Tray can be used with smaller tray (RPI Part #PCT142)

Fits: Chamber

Model: Validator 10 (AA) & Validator Plus 10 (AD)





RPI Part #TUK163 OEM Part #TRH173-0005 (Tray Holder)/ TRY173-0008 Tray (Wire)

#### **HOLDER AND TRAY KIT**

- Includes all parts as shown
- All parts also sold separately!
- · Material: Stainless steel

Fits: Chamber

Fits


Models: Valueklave (1730MKV) & 1730E/EK/M/MK

RPI Part #TUH164 OEM Part #TRH173-0005 (Replaces CU510010) TRAY HOLDER

RPI Part #TUT165 OEM Part #TRY173-0008 (Replaces CU520010)

TRAY (WIRE)

• Size: 5.46" x 12.25"



RPI Part #TUK166 OEM Part #TRH234-0012 (Tray Holder)/ TRY254-0003 Tray (Wire)

#### **HOLDER AND TRAY KIT**

- Includes all parts as shown
- All parts also sold separately!
- Material: Stainless steel

Fits: Chamber Models: EZ9, EZ9k & 2340E/EA/EK/EKA/M/MK

Fits

RPI Part #TUH167 OEM Part #TRH234-0012 (Replaces CT510010)


**TRAY HOLDER** 

RPI Part #TUT168 OEM Part #TRY254-0003 (Replaces CT520010)

TRAY (WIRE)

• Size:17.39" x 7.60"

## **TUTTNAUER** EZ10 & 2540 Series



RPI Part #TUK169 OEM Part #TRH254-0036 (Tray Holder)/ TRY254-0003 Tray (Wire)

#### **HOLDER AND TRAY KIT**

- Includes all parts as shown
- All parts also sold separately!
- · Material: Stainless steel

Fits: Chamber

Models: EZ10, EZ10k & 2540E/EK/EKA/M/MK


RPI Part #TUH170 OEM Part #TRH254-0036 (Replaces CV510010) TRAY HOLDER

RPI Part #TUT168 OEM Part #TRY254-0003 (Replaces CT520010)

TRAY (WIRE)

• Size:17.39" x 7.60"

## **TUTTNAUER** 3870 Series



RPI Part #TUK171

OEM Part #CC510010 (Tray Holder)/
CC520010 (Large Tray)/
CC520020 (Small Tray)

#### **HOLDER AND TRAY KIT**

- Includes all parts as shown
- All parts also sold separately!
- · Material: Stainless steel

Fits: Chamber

Models: 3870E/EA/EK/EKA/M/MK

RPI Part #TUH172 OEM Part #CC510010 TRAY HOLDER

RPI Part #TUT173 OEM Part #CC520020

TRAY (SMALL)

• Size: 27.21" x 11.87"

RPI Part #TUT174 OEM Part #CC520010

TRAY (LARGE)

• Size: 27.21" x 14.43"

## Keep Your Sterilizer Chambers Clean And Performing At Their Best With RPI Sterilizer Cleaners!



· Sterilizer cleaning solution

• (1) 16 oz. bottle Models: Midmark® • Ritter® sterilizers



 Also available by the case of 12 bottles! Just ask for MID-CLEAN™ (CASE) (RPI Part #MIC144)

#### Fits **PELTON & CRANE**

RPI Part #PCC249 OEM Part #047508 (Omni Cleaner Plus™)

#### PEL-CLEAN™

- · Sterilizer cleaning solution
- (1) 16 oz. bottle

Models: Pelton & Crane sterilizers



 Also available by the case of 12 bottles! Just ask for PEL-CLEAN™ (CASE) (RPI Part #PCC250)

# Fits

### **TUTTNAUER**

RPI Part #TUC094 OEM Part #CB0010 (Chamber Brite™)

#### TUTT-CLEAN™

- Sterilizer cleaner
- (1) Box of 10 packets Models:

Most Tuttnauer® autoclaves including 1730, 2340, 2540, 3850 & 3870 Series



 Also available by the case of 12 boxes! Just ask for TUTT-CLEAN™ (CASE) (RPI Part #TUC095)

# RPI IS Your Best Source For Table Top Sterilizer Racks And Trays!

ADDRESS SERVICES REQUESTED

PERMIT NO. 2399 АЭ (SYUN NAV U.S. POSTAGE PAID ORAGNATE GETROSER

625 Cochran Street • Simi Valley, California • 93065

# replacement parts industries, inc.



# Keep Your Sterilizer Racks and Trays Clean and Shiny With RPI's Sci-Dry™ And Cleaning Kit!

Repels moisture and condensation on Sterilizer Racks and Trays leaving them clean and bright!

RPI Part #SCA054 OEM Part #20ZPLUS (Stat-Dri Plus™) / 80ZPLUST (Stat-Dri PlusT™)

#### SCI-DRY™ (2 oz.)

- Drying and rinse agent
- (1) 2 oz. spray bottle

Models: Most Midmark® • Ritter®, Pelton & Crane, SciCan®, and Tuttnauer table top sterilizers



· Also available by the case of 12 bottles! Just ask for SCI-DRY™ (2 oz. - CASE) (RPI Part #SCA055)

Everything you need to clean the sterilizer chamber, racks and trays, the groove where the door gasket installs, and more in one convenient Cleaning Kit!

## RPI Part #RPK791

**OEM Part #** (No OEM Part # Available)

