



... Proud to be a **Healthcare Components Group** company ...



**HEALTHCARE
Components Group**

STERIS

New PM Packs and more to fit ...

Steris[®]

V-Pro[®]

Low Temperature Sterilization Systems



**RPI Part #SSP069
OEM Part #P764335-557**

PM PACK (VACUUM)

- Includes: Vacuum Pump Oil (Leybold Pump) (RPI Part #RXLO34), Cartridge Filter (RPI Part #SSF067), Filter Element (RPI Part #SSE070), (1 pc) Disposal Bag (RPI Part #RPB838), Diaphragm (Valve Umbrella) (RPI Part #SSD071), (1 pkg) Sealing Washer (RPI Part #SSW072), and (3 ft.) Tubing (5/16" ID x 9/16" OD Natural) (RPI Part #RXT077)

Models: V-PRO 60 & V-PRO s2



**RPI Part #SSP063
OEM Part #P764335-558**

PM PACK (GENERAL)

- Includes: HEPA Filter (RPI Part #SSF059), (1 pc) Male Nut (1/8" OD Tubing) (RPI Part #RXF073), (1 pc) Ferrule (1/8" OD Tubing) (RPI Part #RXF074), and (1 pkg) Filter (RPI Part #SSF064)

Models: V-PRO 60 & V-PRO s2



**RPI Part #SSP075
OEM Part #P764337-176**

PM PACK (VACUUM)

- Includes: (3) Vacuum Pump Oil (Leybold Pump) (RPI Part #RXLO34), (2) Cartridge Filter (RPI Part #SSF067), (2) Filter Element (RPI Part #SSE070), (1 pc) Disposal Bag (RPI Part #RPB838), (2) Diaphragm (Valve Umbrella) (RPI Part #SSD071), (2 pkgs) Sealing Washer (RPI Part #SSW072), Drain Plug Gasket (RPI Part #SSG068), Exhaust Filter (RPI Part #SSF060), (5 pcs) Absorbent Pad (RPI Part #RPP958) and (3 ft.) Tubing (5/16" ID x 9/16" OD Natural) (RPI Part #RXT077)

Models: V-PRO Max 2



**RPI Part #SSP058
OEM Part #P764333-354**

PM PACK (GENERAL)

- Includes: HEPA Filter (RPI Part #SSF059), (4) Exhaust Filter (RPI Part #SSF060), (1 pc) Male Nut (1/8" OD Tubing) (RPI Part #RXF073), (1 pc) Ferrule (1/8" OD Tubing) (RPI Part #RXF074), Tee (1/8" OD Tubing) (RPI Part #RXF075), Battery (Lithium) (RPI Part #RXB076), Air Intake Filter (RPI Part #CMF282), Filter Element (RPI Part #SSE061), Solenoid Valve Repair Kit For 1/4" Port Valves (RPI Part #AMK175), and (2) Isolation Valve (3-Way) (RPI Part #SSV062)

Models: V-PRO 1, V-PRO 1+ & V-PRO Max



**RPI Part #SSP065
OEM Part #P764334-481**

PM PACK (VACUUM)

- Includes: Air Intake Filter (Metal Style) (RPI Part #CMF002), Exhaust Filter (RPI Part #SSF060), Filter Element (RPI Part #SSE066), (2) Cartridge Filter (RPI Part #SSF067), (1 pc) Disposal Bag (RPI Part #RPB838), (5 pcs) Absorbent Pad (RPI Part #RPP958), (3) Vacuum Pump Oil (Leybold Pump) (RPI Part #RXLO34), and Drain Plug Gasket (RPI Part #SSG068)

Models: V-PRO 1, V-PRO 1+ & V-PRO Max



replacement parts industries, inc. A Healthcare Components Group Company

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@ripiparts.com • WEBSITE www.ripiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit the equipment noted. All product names used in this document are trademarks or registered trademarks of their respective holders. 1908SS Rd 2 SHOW 11/22

SECTION

10

PAGE

507

UPDATED

11/22

RPI Part #SSF059
OEM Part #P101006-172

HEPA FILTER

- 1.2 micron
- Size: 2.50" OD x 7.60" lg.
- Ports: 3/8" FPT with 1/2" Barbs
- Included in: PM Pack (General) (RPI Part #'s SSP058 and SSP063)

Fits: Air Intake

Models: V-PRO 1, V-PRO 1+, V-PRO 60, V-PRO Max, V-PRO Max 2 & V-PRO s2




RPI Part #SSV062
OEM Part #P093940-125

ISOLATION VALVE (3-WAY)

- 24VDC, 1.5W; 19.70" lg.
- #24 AWG wires
- Ports: 1/4-28; 3-Way
- Included in: PM Pack (General) (RPI Part #SSP058)

Fits: Injection Assembly

Models: V-PRO 1, V-PRO 1+ & V-PRO Max



RPI Part #AMK175
OEM Part #P910009-024

SOLENOID VALVE REPAIR KIT FOR 1/4" PORT VALVES

- Includes all parts as shown
- Included in: PM Pack (General) (RPI Part #SSP058)

Fits: Solenoid Valve

Models: V-PRO 1, V-PRO 1+ & V-PRO Max



RPI Part #RXB076
OEM Part #P764332-350

BATTERY (LITHIUM)

- 3VDC Lithium; Size: .62" dia x 1.00" lg
- Included in: PM Pack (General) (RPI Part #SSP058)

Fits: Control Box Assembly

Models: V-PRO 1, V-PRO 1+ & V-PRO Max



RPI Part #RXF075
OEM Part #P387348-974

TEE (1/8" OD TUBING)

- 1/4-28 for 1/8" OD tubing fitting; Material: PEEK®
- Included in: PM Pack (General) (RPI Part #SSP058)

Fits: Hydrogen Peroxide Delivery System

Models: V-PRO 1, V-PRO 1+, V-PRO 60, V-PRO Max, V-PRO Max 2 & V-PRO s2



RPI Part #SSE061
OEM Part #P764333-215

FILTER ELEMENT

- 20 micron; Size: .81" OD x .55" ID x 1.22" lg.
- Included in: PM Pack (General) (RPI Part #SSP058)

Fits: Air Intake

Models: V-PRO 1, V-PRO 1+ & V-PRO Max



RPI Part #RXF073
OEM Part #P387349-030

MALE NUT (1/8" OD TUBING)

- 6 per package
- 1/4-28 x .61" lg for 1/8" OD tubing; Material: PEEK®
- Included in: PM Pack (General) (RPI Part #'s SSP058 and SSP063)

Fits: Hydrogen Peroxide Delivery System

Models: V-PRO 1, V-PRO 1+, V-PRO 60, V-PRO Max, V-PRO Max 2 & V-PRO s2



RPI Part #CMF282
OEM Part #P764333-087

AIR INTAKE FILTER

- 50 Microns
- Size: 2" dia. x 3.40" lg. x 1/4" MPT
- Includes: Air Filter Element (RPI Part #CME288)

• Included in: PM Pack (General) (RPI Part #SSP058)

Fits: Compressor

Models: V-PRO 1, V-PRO 1+ & V-PRO Max




RPI Part #RXF074
OEM Part #P387349-031

FERRULE (1/8" OD TUBING)

- 6 per package
- For 1/8" OD tubing; Material: Tefzel®
- Included in: PM Pack (General) (RPI Part #'s SSP058 and SSP063)

Fits: Hydrogen Peroxide Delivery System

Models: V-PRO 1, V-PRO 1+, V-PRO 60, V-PRO Max, V-PRO Max 2 & V-PRO s2



RPI Part #CME288
OEM Part # (No OEM Part # Available)

AIR FILTER ELEMENT

- 50 Microns; Size: 1.30" OD x .62" ID x 1.00" lg.
- Included in: Air Intake Filter (RPI Part #CMF282)

Fits: Air Intake Filter Housing

Models: V-PRO 1, V-PRO 1+ & V-PRO Max



RPI Part #SSF064
OEM Part #P129385-263

FILTER

- 2 per package
- Size: .30" dia. x .90" lg; Port: M5 male
- Included in: PM Pack (General) (RPI Part #SSP063)

Fits: Injection Manifold (60), 3-way valve (1, 1+, and Max)

Models: V-PRO 1, V-PRO 1+, V-PRO 60, V-PRO Max, V-PRO Max 2 & V-PRO s2



RPI Part #SSF060
OEM Part #P129385-276

EXHAUST FILTER

- Size: .67" dia x 1.50" lg; 1/4" MPT
- Included in: PM Pack (General) (RPI Part #SSP058) and PM Pack (Vacuum) (RPI Part #'s SSP065 & SSP075)

Fits: Auto Drain Assembly

Models: V-PRO 1, V-PRO 1+, V-PRO Max & V-PRO Max 2



RPI Part #SSE070
OEM Part #P093940-310



FILTER ELEMENT

- Size: 6.83" OD x 7.12" lg.
- Oil vapor and odor removing activated carbon
- Included in: PM Pack (Vacuum) (RPI Part #'s SSP069 & SSP075)

Fits: Vacuum Exhaust Assembly

Models: V-PRO 60, V-PRO s2 & V-PRO Max 2



RPI Part #RXT077
OEM Part #P764335-465

TUBING (5/16" ID X 9/16" OD CLEAR)

- Sold by the foot
- Material: Industrial Grade PVC (Non-Food Grade)
- Included in: PM Pack (Vacuum) (RPI Part #'s SSP069 & SSP075)

Fits: Vacuum Pump Drain Use Only

Models: V-PRO 60, V-PRO s2 & V-PRO Max 2



RPI Part #SSD071
OEM Part #P129603-002

DIAPHRAGM (VALVE UMBRELLA)

- Material: Blue Silicone
- Size: 18mm (.709") OD x 14.3mm (.563") lg.
- Included in: PM Pack (Vacuum) (RPI Part #'s SSP069 & SSP075)

Fits: Oil Mist Eliminator

Models: V-PRO 60, V-PRO s2 & V-PRO Max 2




RPI Part #SSW072
OEM Part #P129603-014

SEALING WASHER

- 2 per package
- Material: Red Silicone
- Size: .290" ID x .688" OD x .093" thk
- Included in: PM Pack (Vacuum) (RPI Part #'s SSP069 & SSP075)

Fits: Oil Mist Eliminator

Models: V-PRO 60, V-PRO s2 & V-PRO Max 2



x2

RPI Part #SSE066
OEM Part #P764334-368

FILTER ELEMENT

- 10 micron; Size: 7.87" OD x 9.62" lg.
- Included in: PM Pack (Vacuum) (RPI Part #SSP065)

Fits: Air Intake

Models: V-PRO 1, V-PRO 1+ & V-PRO Max



RPI Part #CMF002
OEM Part #P129385-485

AIR INTAKE FILTER (METAL STYLE)

- Size: 1.84" OD x 3.63" lg.
- Auto Muffler MO2
- Included in: PM Pack (Vacuum) (RPI Part #SSP065)

Models: V-PRO 1, V-PRO 1+ & V-PRO Max



RPI Part #SSF067
OEM Part #P764334-369

CARTRIDGE FILTER

- 0.3 micron
- Size: 2.78" OD x 7.87" lg.
- Includes: Filter and (2) O-rings
- Included in: PM Pack (Vacuum) (RPI Part #'s SSP065, SSP069 & SSP075)

Fits: Oil Mist Eliminator

Models: V-PRO 1, V-PRO 1+, V-PRO 60, V-PRO Max, V-PRO Max 2 & V-PRO s2



RPI Part #SSG068
OEM Part #P764332-377

DRAIN PLUG GASKET

- Material: Viton®
- Size: .615" ID x .875" OD x .080" thk
- Included in: PM Pack (Vacuum) (RPI Part #'s SSP065 & SSP075)

Fits: Oil Mist Eliminator

Models: V-PRO 1, V-PRO 1+, V-PRO Max & V-PRO Max 2



RPI Part #RXLO34
OEM Part #P764334-644/P764335-529

VACUUM PUMP OIL (LEYBOLD PUMP)

- 1 liter (32 fl. oz.) bottle
- Molykote L-1668 FM oil
- Includes: Yorker spout
- Included in: PM Pack (Vacuum) (RPI Part #'s SSP065, SSP069 & SSP075)
- Also available: Bottle (1 liter / 32 fl. oz.) (RPI Part #RPB939), Disposal Bag (RPI Part #RPB838) and Absorbent Pad (RPI Part #RPP958)

Fits: Vacuum Pump

Models: V-PRO 1, V-PRO 1+, V-PRO 60, V-PRO Max, V-PRO Max 2 & V-PRO s2



RPI Part #RPB838
OEM Part # (No OEM Part # Available)

DISPOSAL BAG

- 3 per package
- Used oil collection and disposal bag
- Included in: PM Pack (Vacuum) (RPI Part #'s SSP065, SSP069 & SSP075)

Models: V-PRO 1, V-PRO 1+, V-PRO 60, V-PRO Max, V-PRO Max 2 & V-PRO s2




RPI Part #RPP958
OEM Part # (No OEM Part # Available)

ABSORBENT PAD

- 3 per package; Size: 9-1/4" x 14-1/2"
- Can be used for oil drip collection and spills
- Included in: PM Pack (Vacuum) (RPI Part #'s SSP065, SSP069 & SSP075)

Models: V-PRO 1, V-PRO 1+, V-PRO 60, V-PRO Max, V-PRO Max 2 & V-PRO s2



RPI Parts to fit Steris® V-Pro® Low Temperature Sterilization Systems!

RPI Part #	OEM Part #	Description	V-Pro 1	V-Pro 1+	V-Pro 60	V-Pro Max	V-Pro Max 2	V-Pro s2
AMK175	P910009-024	SOLENOID VALVE REPAIR KIT FOR 1/4" PORT VALVES	•	•		•		
CME288	<i>(OEM Part # Not Available)</i>	AIR FILTER ELEMENT	•	•		•		
CMF002	P129385-485	AIR INTAKE FILTER (METAL STYLE)	•	•		•		
CMF282	P764333-087	AIR INTAKE FILTER	•	•		•		
RPB838	<i>(OEM Part # Not Available)</i>	DISPOSAL BAG	•	•	•	•	•	•
RPP958	<i>(OEM Part # Not Available)</i>	ABSORBENT PAD	•	•	•	•	•	•
RXB076	P764332-350	BATTERY (LITHIUM)	•	•		•		
RXF073	P387349-030	MALE NUT (1/8" OD TUBING)	•	•	•	•	•	•
RXF074	P387349-031	FERRULE (1/8" OD TUBING)	•	•	•	•	•	•
RXF075	P387348-974	TEE (1/8" OD TUBING)	•	•	•	•	•	•
RXL034	P764334-644 & P764335-529	VACUUM PUMP OIL (LEYBOLD PUMP)	•	•	•	•	•	•
RXT077	P764335-465	TUBING (5/16" ID X 9/16" OD CLEAR)			•		•	•
SSD071	P129603-002	DIAPHRAGM (VALVE UMBRELLA)			•		•	•
SSE061	P764333-215	FILTER ELEMENT	•	•		•		
SSE066	P764334-368	FILTER ELEMENT	•	•		•		
SSE070	P093940-310	FILTER ELEMENT			•		•	•
SSF059	P101006-172	HEPA FILTER	•	•	•	•	•	•
SSF060	P129385-276	EXHAUST FILTER	•	•		•	•	
SSF064	P129385-263	FILTER	•	•	•	•	•	•
SSF067	P764334-369	CARTRIDGE FILTER	•	•	•	•	•	•
SSG068	P764332-377	DRAIN PLUG GASKET	•	•		•	•	
SSP058	P764333-354	PM PACK (GENERAL)	•	•		•		
SSP063	P764335-558	PM PACK (GENERAL)			•			•
SSP065	P764334-481	PM PACK (VACUUM)	•	•		•		
SSP069	P764335-557	PM PACK (VACUUM)			•			•
SSP075	P764337-176	PM PACK (VACUUM)					•	
SSV062	P093940-125	ISOLATION VALVE (3-WAY)	•	•		•		
SSW072	P129603-014	SEALING WASHER			•		•	•