



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

Just in time for the Holidays ...

New Pump Kits to fit SciCan® Stat/IM® 2000 & 5000 Cassette Sterilizers now available!




RPI Part #SCP050
OEM Part #01-100568S (Invensys Pump - STAT/IM 2000)/
01-101761S (Invensys Pump - STAT/IM 5000)

PUMP KIT

- 120Volts, 60Hz.
- This Invensys brand pump will accommodate STAT/IM 2000 and 5000 sterilizers with PCB Software Rev. 6.4 or older only
- Includes all parts as shown

Fits: Chassis

Models: 2000 & 5000




RPI Part #SCP051
OEM Part #01-110445S (Ulka Pump -STAT/IM 2000)/
01-110446S (Ulka Pump -STAT/IM 5000)

PUMP KIT

- 120Volts, 60Hz.
- This Ulka brand pump will accommodate STAT/IM 2000 and 5000 sterilizers with PCB Software Rev. 7 and above only
- Includes all parts as shown

Fits: Chassis

Models: 2000 & 5000




RPI Part #SCT053
OEM Part #01-108816S

METERING TUBE (290)

- Use as replacement pump tube or as pump recovery tube
- Size: 1/8" OD tubing x 5-1/2" lg.
- Also included in: Pump Kits (RPI Part #'s SCP050 & SCP051)

Fits: Pump to Boiler

Models: 2000 & 5000



For a complete listing of parts to fit
SciCan® Stat/IM® 2000 & 5000 Cassette Sterilizers,
 see www.rpiparts.com

**Plus even more new parts
 coming soon!**