

## NEW PARTS TO FIT STERIS SYSTEM 1

### STERILANT ASPIRATOR ASSEMBLY


RPI Part #SSA029

OEM Part #A1600

- Includes sterilant aspirator and tubing (RPI Part #SST030)

Fits: Instrument trays

Model: System 1



Tubing also sold separately... see RPI Part #SST030 below.

### CHECK VALVE with O-RINGS

RPI Part #SSV004

OEM Part #400390

- Includes check valve and (1) O-ring (RPI Part #RPO343) & (2) O-ring (RPI Part #RPO388)

Fits: Air Manifold Assembly

Model: System 1



O-rings also sold separately... see RPI Part #RPK036 below.

### SILICONE TUBING


RPI Part #SST030

OEM Part #N/A

- .250" OD x .125" ID x .063" Wall
- Sold by the foot
- Medical grade Silicone

Fits: Sterilant Aspirator Assembly (RPI Part #SSA029)

Model: System 1



### CHECK VALVE O-RING KIT


RPI Part #SSK036

OEM Part #N/A

- Includes o-rings: (1) RPI Part #RPO343 & (2) RPI Part #RPO388
- Material: Red Silicone
- Hardness: 70 Duro

Fits: Check Valve (RPI Part #SSV004)

Model: System 1



RPI Part #SSK036  
Includes o-rings: (1) RPO343 & (2) RPO388  
O-rings also sold separately... see RPI Part #RPO343 & RPO388

RPI Part #RPO343: ID: 5/16", OD: 7/16", C/S: 1/16"

RPI Part #RPO388: ID: 7/16", OD: 9/16", C/S: 1/16"

### FILTER SCREEN & O-RING KIT

RPI Part #SSK035

OEM Part #N/A

- Includes filter screen and O-ring (RPI Part #RPO389)

Fits: Drain section of instrument tray

Model: System 1



O-ring also sold separately... see RPI Part #RPO389

### S.S. REPAIR KIT (CK-8 CHECK VALVE)


RPI Part #SSK034

OEM Part #200696

- Includes stainless steel spring and (1) each of O-ring RPI Part #RPO303 & RPO318

Fits: CK-8 Check Valve

Model: System 1



O-rings also sold separately... see RPI Part #RPO303 & RPO318

OEM PART #	RPI PART #	DESCRIPTION
A1501	SSF013	PREFILTER CARTRIDGE A1501
A1504	SSF014	PREFILTER CARTRIDGE A1504
A1600	SSA029	STERILANT ASPIRATOR ASSEMBLY
200226	SSK006	CHECK VALVE KIT
200056	SSK010	GAS SPRING ASSEMBLY
200696	SSK034	S.S. REPAIR KIT (CK-8 CHECK VALVE)
400345	SSV002	CHECK VALVE (CK-1)
400368	SSV005	CHECK VALVE (S.S.)
400390	SSV004	CHECK VALVE with O-RINGS
400420	SSS007	PINCH VALVE SLEEVE
400510	SSV003	CHECK VALVE (S.S. CK-2 & CK-3)
400696	SSK008	CHECK VALVE (CK-8 DRAIN CHECK) RETRO FIT KIT
N/A	MDS017	FLOAT ASSEMBLY
N/A	SSK035	FILTER SCREEN AND O-RING
N/A	SSK036	CHECK VALVE O-RING KIT (6 of each o-ring per package)
N/A	SST030	SILICONE TUBING
2-011 (Industry part #)	RPO343	O-RING (12 per package)
2-012 (Industry part #)	RPO303	O-RING (12 per package)
2-013 (Industry part #)	RPO388	O-RING (6 per package)
2-017 (Industry part #)	RPO318	O-RING (12 per package)
2-121 (Industry part #)	RPO389	O-RING (6 per package)



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@rpiparts.com](mailto:order@rpiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Steris equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. w/corrections 09/03