



Service Bulletin

Keeping You Informed from the Leader in Replacement Parts for Healthcare Equipment since 1972



It's Back and Ready to Ship – the Heater Bar Assembly that fits Air Techniques A/T2000 & A/T2000 Plus !

RPI Part #ATH639
OEM Part #43290

HEATER BAR ASSEMBLY
Fits Air Techniques models
A/T2000 & A/T2000 Plus



Take a Moment to Update Your RPI Catalog – these Hospital Grade Power Cords are 13A !
Section 6, Page 300

RPI Part #RPC557 (8 ft. long)
RPI Part #RPC558 (18 ft. long)
RPI Part #RPC559 (25 ft. long)
HOSPITAL GRADE POWER CORD
• 16/3 SJT, 60°C
13A @ 125VAC



Redesigned for A Better Fit and Now In Stock – the Transfer Arm Gears that fit the Air Techniques Peri- Pro, II & III !

RPI Part #ATG613
OEM Part #90525
“L” TRANSFER ARM GEAR
Fits Air Techniques
Peri-Pro, II & III

RPI Part #ATG615
OEM Part #90518
TRANSFER ARM GEAR
Fits Air Techniques
Peri-Pro, II & III



Corrections to Note ... Update Your RPI Catalog Today!
Midmark • Ritter M9 & M11
Section 10, Pages 807 & 809

RPI Part #MIA077
OEM Part #050-3970-00 (Old #050-2623-00)
SPRING ARM
Models with correct Serial # prefixes:
M11 Serial #'s -001 thru -019 with prefixes ES, ET, FP, FR, GB & NP;
M11D Serial #'s -001 thru -019 with prefixes GC & GD



The Autoholder Valve in the Handpiece Hanger Assembly that fits the P&C Spirit is New and It Fits the Handpiece Hanger Perfectly !

RPI Part #PCK783
OEM Part #33 39 970
HANDPIECE HANGER ASSEMBLY
Includes Handpiece Hanger
and the **NEW** Autoholder Valve (RPI Part #PCV793)
Fits Spirit AC

RPI Part #MIA086
OEM Part #050-2362-00
SPRING ARM
Models with correct Serial # prefixes:
M9 Serial #'s -001 thru -019 with prefixes CZ, DA, DB, DX, DY, FK, FL, LA & OM;
M9D Serial #'s -001 thru -019 with prefixes FF & FG

RPI Part #PCV793
OEM Part #33 07 431
AUTOHOLDER VALVE
Fits Spirit AC Handpiece Hanger

RPI Part #PCV052
OEM Part #002-0359-00
PRESSURE RELIEF VALVE
Models with correct Serial # prefixes:
M9 Serial #CZ1000 thru CZ4334, DA1000 thru DA1004, DB1000 thru DB1237, DX1000 thru DX2555, DY1000 thru DY1176, FD1000 thru FD1399, FK1000 thru FK1699, FL1000 thru FL1474 & OM1000 thru OM145248;
M9D Serial #FF1000 thru FF1034 & FG1000 thru FG1112;
M11 Serial #ES1000 thru ES5751, ET1000 thru ET7341, FP1000 thru FP1075, FR1000 thru FR1246 & GB1000 thru GB1129;
M11D Serial #GC1000 thru GC1093 & GD1000 thru GD1102

RPI Part #PCV787
OEM Part # (No OEM Part # Available)
AUTOHOLDER VALVE
Fits Foot Controls and Handpiece Hangers
This part does not fit the Spirit AC Handpiece Hanger

RPI Part #MIS039
OEM Part #002-0362-00
DOOR SWITCH ASSEMBLY

RPI Part #MIC037
OEM Part #052-0173-00
CONDENSATION COIL

RPI Part #MIP050
OEM Part #002-0357-00
TEMPERATURE PROBE ASSEMBLY

RPI Part #MIT047
OEM Part #002-0370-00
OVERTEMP THERMOSTAT

Models with correct Serial # prefixes:
M9 Serial #'s -001 thru -019 with prefixes CZ, DA, DB, DX, DY, FK, FL, LA & OM;
M9D Serial #'s -001 thru -019 with prefixes FF & FG;
M11 Serial #'s -001 thru -019 with prefixes ES, ET, FP, FR, GB & NP;
M11D Serial #'s -001 thru -019 with prefixes GC & GD



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com