Models: NX, 100NX & 100S



# MORE NEW PARTS TO FIT STERRAD NX, 100NX & 100S

#### **HYDROGEN PEROXIDE STERILIZATION SYSTEMS!**




#### RPI Part #SDK050 0EM Part #02-52410-X-XXX

#### VAPORIZER/CONDENSER SERVICE KIT

- Includes: Insulation Block Kit (RPI Part #SDK093), Valve Manifold Assembly (RPI Part #SDA048), Wire Brush (Scrubbing) (RPI Part #RXB059), Wire Brush (1/2" Dia.) (RPI Part #RXB060), Cleanout Tool (.030") (RPI Part #RXT062), (1 pc) Heat Sink Compound (RPI Part #RPC464), (2) 0-Rings,
   (2) cleaning pads and mounting hardware
- Also available: Vacuum Test Gauge Kit (RPI Part #SDK112), Pliers (Soft Jaw) (RPI Part #RXT005) and Insulation Blanket (RPI Part #SDB094)

Fits: Vaporizer/Condenser Assembly

Models: NX





#### RPI Part #SDK051


OEM Part # (No OEM Part # Available)

#### VAPORIZER/CONDENSER SERVICE KIT

- Includes: Insulation Block Kit (RPI Part #SDK093), Valve Manifold Assembly (RPI Part #SDA082), Wire Brush (Scrubbing) (RPI Part #RXB059), Wire Brush (1/2" Dia.) (RPI Part #RXB060), Cleanout Tool (.030") (RPI Part #RXT062), (1 pc) Heat Sink Compound (RPI Part #RPC464), (6) 0-Rings,
   (2) Cleaning pads and mounting hardware
- Also available: Vacuum Test Gauge Kit (RPI Part #SDK112), Pliers (Soft Jaw) (RPI Part #RXT005) and Insulation Blanket (RPI Part #SDB097)

Fits: Vaporizer/Condenser Assembly

Models: 100NX



#### RPI Part #SDK093

OEM Part # (No OEM Part # Available)




#### INSULATION BLOCK KIT

- Includes: Insulation block, (1 pc) High Temperature lubricant (RPI Part #RPL090), 0-Ring Removal Tool (RPI Part #RXT002), (2 pcs) velcro pads (hook only), (4) 0-Rings and mounting hardware
- Also available: Insulation Blanket (RPI Part #SDB094), Insulation Blanket (RPI Part #SDB097),
   Vacuum Test Gage Kit (RPI Part #SDK112) and Pliers (Soft Jaw) (RPI Part #RXT005)

Fits: Vaporizer/Condenser Assembly

Models: NX & 100NX





## replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

10
PAGE

**536** 

UPDATED 07/20


#### RPI Part #SDK090

OEM Part # (No OEM Part # Available)

## INJECTOR ASSEMBLY REBUILD KIT

 Includes: Wire Brush (Scrubbing) (RPI Part #RXB059), (1 pkg) Needle H202 Transfer (RPI Part #SDN073), (1 pc) High Temp Lubricant (RPI Part #RPL090), (2) needle guides and bushings, (3) compression springs, diaphragm, sealing washer, (3) O-Rings, (2) cleaning pads and retaining ring

Fits: Injector Assembly Models: 100S



#### RPI Part #SDK091


**OEM Part** # (No OEM Part # Available)

#### INJECTOR ASSEMBLY TOOL KIT

- . Includes: Needle Aligner (RPI Part #SDT110), Bushing Retainer (RPI Part #SDT111), O-Ring Removal/Gage Tool (RPI Part #RXT082), Cleanout Tool (.067") (RPI Part #RXT061) and Spanner Wrench (RPI Part #RXT044)
- · Custom tools used to install the needles into the Injector Valve Assembly and disassemble/clean the Injector Base

Fits: Injector Valve Assembly

Models: 100S




#### RPI Part #RXT044

OEM Part # (No OEM Part # Available)

#### SPANNER WRENCH

- Size: 6.25" lg. x 2.63" span
- . Material: .120" dia. stainless steel spring wire

Fits: Injector Housing Models: 100S



#### RPI Part #RXT082


OEM Part # (No OEM Part # Available)

#### **0-RING REMOVAL/GAGE TOOL**

- Size: 3.75" lg. x 1/4" dia., plastic body
- Includes: (2) Protective vinyl caps
- . Used to check for existing 0-Rings and remove them from the Injector Assembly.

Fits: Injector Assembly

Models: 100S



#### RPI Part #RXT061

OEM Part # (No OEM Part # Available)

#### CLEANOUT TOOL (.067")

- Size: .067" dia. drill bit x 1/4" dia., plastic body
- · Includes: Protective vinyl cap Fits: Injector Valve Assembly

Models: 100S



#### RPI Part #SDT110

OEM Part # (No OEM Part # Available)

#### **NEEDLE ALIGNER**

- Size: 1.50" L x .75" W
- · Material: Stainless steel

Fits: Injector Assembly Needle Bushings

Models: 100S

#### RPI Part #SDT111 OEM Part # (No OEM Part # Available)

#### BUSHING RETAINER


- Size: 1.50" L x 1.12" W
- · Material: Stainless steel

Fits: Injector Assembly Models: 100S









#### RPI Part #SDK112 OEM Part # (No OEM Part # Available)

#### VACUUM TEST GAUGE KIT

- · Includes all parts as shown
- · All parts also sold separately Fits: Pressure Control Manifold Models: NX, 100NX & 100S

#### RPI Part #RXG039 0EM Part #35-04264-001

#### VACUUM TEST GAUGE

- . Includes: Gauge, 9VDC battery, and storage case
- · Brand: Inficon
- · Range: Atmosphere to 1 Micron
- . Display units: Microns, Pascal, Millibar and Torr
- Also available: Centering Ring (16KF) (RPI Part #RXR042) and Clamp (10KF & 16KF) (RPI Part #RXC063)

Fits: 16KF flange Models: NX, 100NX & 100S




#### RPI Part #RXR042 0EM Part #30-09399-0-001

### CENTERING RING (16KF)

- Size: 16KF (16mm, "KF" style)
- . Material: 303 Stainless steel hub with (1 pc) Viton® (FKM) 0-Ring (RPI Part #RX0043)
- . Also available: Clamp (10KF & 16KF) (RPI Part #RXC063)

Fits: K16 Flange

Models: NX, 100NX & 100S




### RPI Part #RX0043

OEM Part # (No OEM Part # Available)

#### 0-RING

- · 3 per package
- Size: 18mm (.709") ID x 28mm (1.102") OD x 5mm (.197") C/S
- . Material: Viton® (FKM); 75 Duro
- . Included in: Centering Ring (16KF) (RPI Part #RXR042)

Fits: 16KF Centering Ring Models: NX, 100NX, 100S



SHOWN ACTUAL SIZE

#### RPI Part #RXC063

OEM Part # (No OEM Part # Available)

#### **CLAMP (10KF & 16KF)**

- Size: fits 10KF and 16KF fittings
- Material: Aluminum body
- · Also available: Centering Ring (16KF) (RPI Part #RXR042)

Fits: Vacuum Pump Assembly Models: 100S, 100NX & NX



RP0972

(1 pc)

#### RPI Part #SDF092 0EM Part #30-09256-X-001

#### ADAPTOR FITTING

- Size: 16KF flange x 7/16-20 thread
- Material: Stainless steel with (1pc) Viton® (FKM) 0-Ring (RPI Part #RP0972)
- . Also available: Centering Ring (16KF) (RPI

Part #RXR042) and Clamp (10KF & 16KF) (RPI Part #RXC063)

Fits: Vacuum Manifold Models: NX & 100NX

#### RPI Part #RP0972

OEM Part # (No OEM Part # Available)

#### Industry Part #AS568-012

#### 0-RING

- 12 per package
- Size: .364" ID x .504" OD x .070" C/S
- . Material: Viton® (FKM); 75 Duro

Fits: Adapter Fitting and Hex Port Plug

Models: NX & 100NX

SHOWN ACTUAL SIZE

#### RPI Part #RXF078

OEM Part # (No OEM Part # Available)

#### VACUUM REDUCER TEE (1" x 1" x 3/4")

- . Connection: fits 25KF x 25KF x 16KF flange
- · Material: Stainless steel Fits: Vacuum Pump Assembly

Models: 100S



#### RPI Part #RXF079

OEM Part # (No OEM Part # Available)

#### VACUUM RADIUS ELBOW (3/4")

- · Connection: 16KF flange
- · Material: Stainless steel Fits: Vacuum Pump Assembly

Models: 100S



#### RPI Part #RPL615

OEM Part # (No OEM Part # Available)

#### **0-RING LUBRICANT**

- 2 oz. tube
- . Non-Toxic; Super-O-Lube® made by Parker
- Temperature Range: -65° to 400° F
- · Can be used on silicone and non-silicone rubber materials

Models: NX, 100NX & 100S







#### RPI Part #SDN073


**OEM Part #** (No OEM Part # Available)

#### **NEEDLE H202 TRANSFER**

- · 2 per package
- Includes: (4) 0-rings

Fits: Injector Valve Assembly

Model: 100S




#### RPI Part #SDF098 0EM Part #21-01104-040

#### **SCREEN FILTER**

• Size: 40KF (40mm, "KF" style)

 Material: Stainless steel centering ring with a 30 mesh stainless steel screen and (1 pc) Viton® (FKM) 0-Ring (RPI Part #RX0081)

Fits: Vacuum Assembly Models: 100NX & 100S



#### RPI Part #RXB059

**OEM Part** # (No OEM Part # Available)

#### WIRE BRUSH (SCRUBBING)

• Size: 6" body with 1" dia. bristle

· Material: Galvanized stainless steel bristles

Fits: Vaporizer Service Kit Models: NX, 100NX & 100S



#### RPI Part #RXB060

OEM Part # (No OEM Part # Available)

#### WIRE BRUSH (1/2" DIA)

• Size: 8" body with 1/2" dia. medium stiffness bristles

Material: Brass bristles
 Fits: Vaporizer Service Kit
 Models: NX & 100NX



#### RPI Part #RXT062


OEM Part # (No OEM Part # Available)

#### CLEANOUT TOOL (.030")

 Size: .030" dia. drill bit x 1/4" dia., plastic body

Includes: Protective vinyl cap
 Fits: Vaporizer/Condenser Assembly

Models: NX & 100NX



#### RPI Part #RPC464

**OEM Part #** (No OEM Part # Available)

#### **HEAT SINK COMPOUND**

• 5 per package; Size: .07 fl. oz (2ml)

Non-Silicone based thermal grease

Models: NX & 100NX



#### RPI Part #RPL090

OEM Part # (No OEM Part # Available)

#### HIGH TEMP LUBRICANT

• 5 per package; 1/2 gram per unit

• Temperature range: -40°F to +400°F (-40°C to +204°C)

Models: NX, 100NX & 100S





#### RPI Part #RX0081

**OEM Part** # (No OEM Part # Available)

#### 0-RING


• 3 per package

• Size: 42mm (1.653") ID x 52mm (2.047") OD x 5mm (.197") C/S

. Material: Viton® (FKM); 75 Duro

Fits: 40KF Centering Ring

Models: 100S



SHOWN ACTUAL SIZE



#### RPI Part #RX0080

**OEM Part** # (No OEM Part # Available)

#### 0-RING


· 3 per package

 Size: 28mm (1.020") ID x 38mm (1.496") OD x 5mm (.197") C/S

• Material: Viton® (FKM); 75 Duro

Fits: 25KF Centering Ring

Models: 100S



SHOWN ACTUAL SIZE

#### RPI Part #RXT002

OEM Part # (No OEM Part # Available)

#### **O-RING REMOVAL TOOL**

 Made from Nylon® to avoid scratching delicate O-Ring sealing surfaces

Size: 5.5"

Models: NX, 100NX & 100S



#### RPI Part #RXT005

OEM Part # (No OEM Part # Available)

#### PLIERS (SOFT JAW)

 These adjustable pliers are ideal for griping and protecting delicate surfaces without scarring.

• 1.75" max opening; 10" lg.

 Also available: Replacement Pads (RPI Part #RXT006)

Models: NX, 100NX & 100S



#### RPI Part #RXT020

OEM Part # (No OEM Part # Available)

#### SPANNER WRENCH

- Size: 1.25" Hex x 3.00" Ig.
- · Material: Blue anodized aluminum
- Note: Holes through hex body allow use of a screwdriver as a leverage handle.
- Used to install/rebuild Solenoid Valves (RPI Part #'s SDV031, SDV049, SDV096 and SDK052)

Fits: Vent Valve

Models: 100S, NX & 100NX





#### RPI Part #SDV031 0EM Part # (No OEM Part # Available)

#### **SOLENOID VALVE**

- 12VDC, 16W (10Ω)
- · Material: Viton® (FKM) seals
- Also available: Solenoid Valve Rebuild Kit (RPI Part #SDK052) and Spanner Wrench (RPI Part #RXT020)

Fits: Vent Valve Assembly Models: 100NX





## RPI Part #SDV049 0EM Part # (No OEM Part # Available)

#### SOLENOID VALVE

- 24VDC, 15W (39Ω)
- Material: Viton® (FKM) seals
- Also available: Solenoid Valve Rebuild Kit (RPI Part #SDK052) and Spanner Wrench (RPI Part #RXT020)

Fits: Vent Valve Assembly

Models: NX





#### RPI Part #SDV096 OEM Part # (No OEM Part # Available)

#### SOLENOID VALVE

- 120VAC, 18W (41Ω)
- Material: Viton® (FKM) seals
- Also available: Solenoid Valve Rebuild Kit (RPI Part #SDK052) and Spanner Wrench (RPI Part #RXT020)

Fits: Vacuum Manifold Assembly and Oil Recovery Valve

Models: 100S




## RPI Part #SDK052 OEM Part # (No OEM Part # Available)

#### SOLENOID VALVE REBUILD KIT

- Viton® (FKM) seals
- Includes all parts as shown
- Included in: Solenoid Valve (RPI Part #SDV031), Solenoid Valve (RPI Part #SDV049) and Solenoid Valve (RPI Part #SDV096)
- Also available: Spanner Wrench (RPI Part #RXT020)

Fits: Solenoid Valve Models: NX, 100NX & 100S




#### RPI Part #SDF072 0EM Part # (No OEM Part # Available)

#### FAN

- 24VDC: 4300 RPM
- Size: 4.68" Sq. x 1.50" Thk.
- Material: Glass reinforced plastic, black finish; ball bearing construction
- Includes: heat shrink tubing
   Fits: Vaporizer/Condenser

Assembly Models: 100NX




#### RPI Part #SDF095 0EM Part # (No OEM Part # Available)

#### FΔN

- 24VDC; 2750 RPM
- Size: 4.72" Sq. x 1.25" Thk.
- Material: Glass reinforced plastic, black finish; ball bearing construction
- Includes: (1pc) Cable Tie (8" Black) (RPI Part #RPT279) and (4 pcs) Rivet (RPI Part #RXR064)

Fits: Vaporizer/Condenser

Assembly Models: NX



#### RPI Part #RXR064 0EM Part # (No OEM Part # Available)

#### RIVE

- · 6 per package
- Fits hole dia.: .163"; Panel Thk: .188"
- . Material: Nylon® Color: Black
- Size: .250" dia. head x .220" lg.
- Included in: Fan (RPI Part #SDF095)

Fits: Vaporizer/Condenser Fan

Models: NX



#### RPI Part #RPT279

OEM Part # (No OEM Part # Available)

#### CABLE TIE (8" BLACK)

- · 25 per package
- With #10 mounting hole

Fits: Vaporizer/Condenser Fan

Models: NX & 100NX



## RPI Part #RPT083 OEM Part # (No OEM Part # Available)

#### **CABLE TIE (4" WHITE)**

- 25 per package
- Material: Nylon®

Models: NX, 100NX & 100S





# MORE NEW PARTS TO FIT STERRAD NX, 100NX & 100S

#### **HYDROGEN PEROXIDE STERILIZATION SYSTEMS!**

#### Now In Stock and Ready to Ship:

- Fans
- . Injector Assembly Rebuild Kit
- Injector Assembly Tool Kit
- Insulation Blankets
- Insulation Block Kit
- Screen Filter

- Solenoid Valves
- Solenoid Valve Rebuild Kit
- Specialty Tools
- Vacuum Test Gauge Kit
- Vaporizer/Condenser Service Kits
- and more!

See inside for all of the new parts RPI offers to fit Sterrad NX, 100NX and 100S!





ADDRESS SERVICES REQUESTED

PRESORTED STANDARD U.S. POSTAGE PAID VAN NUYS, CA PERMIT NO. 2399

625 Cochran Street • Simi Valley, California • 93065



## replacement parts industries, inc.

#### RPI Part #SDB094


OEM Part # (No OEM Part # Available)

#### INSULATION BLANKET

- Size: 5.0" Wide x 11.00" Long x 3/32" Thick
- Material: Flame retardant Aramide®
- Includes: (2) Velcro strips (hook only)

Fits: Vaporizer/Condenser Assembly

Models: NX



#### RPI Part #SDB097

OEM Part # (No OEM Part # Available)

#### INSULATION BLANKET

- Size: 5.0" Wide x 24.75" Long x 3/32" Thick
- Material: Flame retardant Aramid®
- Includes: (2) Velcro strips (hook only) Fits: Vaporizer/Condenser Assembly

Models: 100NX

