



Parts to fit
Stryker® 900 Series
Cast Cutters!



For Synthetic & Plaster Casts!

RPI Part #STB056
 OEM Part #940-23

CAST CUTTER BLADE (PiN)

- For both plaster and synthetic casts
- Size: 2.36" oval blade
- Material: Plasma ion nitride treated stainless steel
- Color: Matte grey
- Also available: Cast Cutter Blade (RPI Part #STB043), Cast Cutter Blade (TiN) (RPI Part #STB045), and Blade Guard (RPI Part #STG042)

Fits: Head Assembly
Models: 940 & 941



For Plaster Casts!

RPI Part #STB043
 OEM Part #940-25

CAST CUTTER BLADE

- For plaster casts
- Size: 2.36" oval blade
- Material: Stainless steel
- Color: Silver
- Also available: Cast Cutter Blade (TiN) (RPI Part #STB045), Cast Cutter Blade (PiN) (RPI Part #STB056), and Blade Guard (RPI Part #STG042)

Fits: Head Assembly
Models: 940 & 941



For Synthetic & Plaster Casts!

RPI Part #STB045
 OEM Part #940-35

CAST CUTTER BLADE (TiN)

- For both plaster and synthetic casts
- Size: 2.36" oval blade
- Material: Titanium nitride coated stainless steel
- Color: Gold
- Also available: Cast Cutter Blade (RPI Part #STB043), Cast Cutter Blade (PiN) (RPI Part #STB056), and Blade Guard (RPI Part #STG042)

Fits: Head Assembly
Models: 940 & 941



RPI Part #STG042
 OEM Part #940-1-134

BLADE GUARD

- Fits: 2.36" oval blade
- Material: Clear blue plastic
- Includes: Mounting hardware
- Also available: Cast Cutter Blade (RPI Part #STB043), Cast Cutter Blade (TiN) (RPI Part #STB045), and Cast Cutter Blade (PiN) (RPI Part #STB056)

Fits: Head Assembly
Models: 940 & 941



RPI Part #STF044
 OEM Part #986-500

CARTRIDGE FILTER

- 11 micron
- Size: 7.00" dia. x 11.38" lg.
- Includes: Filter with O-Ring

Fits: Vacuum Canister
Models: 986 Plaster Vacuum

For more RPI parts to fit Stryker 800 and 900 Series Cast Cutters as well as 800 Series Autopsy Saws, be sure to see the next page!

Complete List of RPI Parts to fit Stryker®

800 & 900 Series Cast Cutters and 800 Series Autopsy Saws!

RPI Part #	OEM Part #	Description	800 & 900 Series Cast Cutters						986 Plaster Vacuum	800 Series Autopsy Saws	
			840	841	848	851	940	941	986	810	811
PHC015	875-033007	Retaining Clip - 5/Pkg					X	X			
RPC039	<i>(No OEM Part # Available)</i>	Brush Cap - 6/Pkg	X	X	X	X				X	X
RPC667	58-18	Brush Cap - 6/Pkg					X	X			
RPC671	<i>(No OEM Part # Available)</i>	Cable (4 Conductor)					X	X			
RPF665	<i>(No OEM Part # Available)</i>	Axial Lead Fuse (2-1/2A) - 5/Pkg					X				
RPL680	<i>(No OEM Part # Available)</i>	Bearing Grease - 5/Pkg					X	X			
RPO335	45-29	O-Ring - 12/Pkg					X	X			
RPS639	<i>(No OEM Part # Available)</i>	Silicone Gasket Maker 598					X	X			
RPS672	<i>(No OEM Part # Available)</i>	Cord Strap (White 6") - 3/Pkg					X	X			
RPT172	<i>(No OEM Part # Available)</i>	Strain Relief Bushing Pliers					X	X			
RPT700	<i>(No OEM Part # Available)</i>	T-15 Torx Wrench					X	X			
STA026	870-1-8	Armature Assembly	X		X					X	
STA028	<i>(No OEM Part # Available)</i>	Armature Assembly					X				
STA040	12-24, 21-129, 45-29, 50-33, 940-1-115, 940-1-118 & 940-1-125	Output Drive Assembly					X	X			
STB002	870-1-13, 940-1-113 (Brush only), 38-49 (Spring only)	Brush - 2/Pkg	X	X	X	X	X	X		X	X
STB003	81-117	Front Bearing	X	X	X	X				X	X
STB004	81-118/81-147	Needle Bearing	X	X	X	X	X	X		X	X
STB005	81-138	Rear Bearing	X	X	X	X				X	X
STB015	840-25-11	Cast Cutter Blade (2-1/2")	X	X	X	X					
STB031	<i>(No OEM Part # Available)</i>	Front Bearing					X	X			
STB032	<i>(No OEM Part # Available)</i>	Rear Bearing					X	X			
STB038	<i>(No OEM Part # Available)</i>	Strain Relief Bushing - 5/Pkg					X	X			
STB039	<i>(No OEM Part # Available)</i>	Power Supply PCB					X				
STB043	940-25	Cast Cutter Blade					X	X			
STB045	940-35	Cast Cutter Blade (TiN)					X	X			
STB056	940-23	Cast Cutter (PiN)					X	X			
STC009	840-2-7	Cam	X	X	X	X				X	X
STC011	840-1-13	Collar	X	X	X	X				X	X
STC029	<i>(No OEM Part # Available)</i>	Cam					X	X			
STD012	840-2-4	Felt Oil Distributor	X	X	X	X				X	X
STF027	840-3-32	Field Assembly	X		X					X	
STF036	<i>(No OEM Part # Available)</i>	Thermal Fuse (104°C)					X	X			
STF044	986-500	Cartridge Filter							X		
STG013	840-1-10	Gasket	X	X	X	X				X	X
STG042	940-1-134	Blade Guard					X	X			
STK030	28-111 (Retaining Ring only), 11-275 (Washer only), 81-118/ 81-147 (Needle Bearing only)	Roller Bearing kit					X	X			
STK033	<i>(No OEM Part # Available)</i>	Armature Washer Kit					X	X			
STK037	<i>(No OEM Part # Available)</i>	Power Cable Kit					X	X			
STL016	840-1-5 (800 Series only)	Shaft Link	X	X	X	X	X	X		X	X
STL035	<i>(No OEM Part # Available)</i>	Power Switch Label - 3/Pkg					X	X			
STN017	840-1-14	Shaft Nut	X	X	X	X				X	X
STS018	840-3-21	Switch Assembly	X	X	X	X				X	X
STS021	47-5	Oil Seal	X	X	X	X				X	X
STS023	870-1-3	Shaft	X	X	X	X					
STS034	<i>(No OEM Part # Available)</i>	Power Switch Label - 3/Pkg					X	X			
STT010	<i>(No OEM Part # Available)</i>	Cam Tool	X	X	X	X	X	X		X	X
STW025	13-4	Internal Tooth Washer	X	X	X	X				X	X