



... Proud to be a Healthcare Components Group company ...



HEALTHCARE
Components Group

STRYKER


New Parts to fit

Stryker® 900 Series

Plaster Vacuums & Cast Cutters

Did You Know?

To remove the case screws for the Stryker 986 CastVac™, simply remove the perforated metal plate to access the Nylok nuts. Note: Removing the perforated metal plate will rip/tear the Motor Gasket (RPI Part #STG049).



replacement parts industries, inc. A Healthcare Components Group Company

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit the equipment noted. All product names used in this document are trademarks or registered trademarks of their respective holders. 2131ST Web 08/25

SECTION
9

PAGE
308

UPDATED
08/25



Parts to fit Stryker® 900 Series Plaster Vacuums



RPI Part #STH052
OEM Part #986-1-30

VACUUM HOSE

- 104" lg., 1.0" OD, White
- Includes: (2 pcs) Vacuum Cuff (RPI Part #STC051) and (6 pcs) Cord Strap (White 6") (RPI Part #RPS672)

Fits: Chassis

Models: 986, 987 & 996



RPI Part STK048
OEM Part # (No OEM Part # Available)

MOTOR KIT

- 120VAC, thru-flow discharge, single speed, double ball bearings
- Includes: Motor Gasket (RPI Part #STG049), (1 pkg) Grommet (3/16" ID x 5/16" Groove) (RPI Part #RXG107), and mounting hardware

Fits: Chassis

Models: 986 & 987



x2

RPI Part #STC051
OEM Part # (No OEM Part # Available)

VACUUM CUFF

- 2 per package
- Size: 1" dia., white TPE
- Includes: Lock button and spring
- Included in: Vacuum Hose (RPI Part #STH052)

Fits: Vacuum Hose

Models: 986, 987 & 996



RPI Part #STG049
OEM Part # (No OEM Part # Available)

MOTOR GASKET

- Size: 4.00" ID x 5.00" OD x .63" thk., PSA on top and bottom with release liner
- Material: Polyurethane closed cell foam
- Included in: Motor Kit (RPI Part #STK048)

Fits: Motor

Models: 986, 987 & 996



RPI Part #STB046
OEM Part # (No OEM Part # Available)

CONTROL PC BOARD

- Size: 4.00" x 3.00"
- Includes: Installed 1/4" quick connect terminals and #6-32 x .50" threaded standoffs

Fits: Chassis

Models: 986 & 987



RPI Part #STF044
OEM Part #986-500

CARTRIDGE FILTER


- 11 micron
- Size: 7.00" dia. x 11.38" lg.
- Includes: Filter with O-Ring

Fits: Vacuum Canister

Models: 986 Plaster Vacuum



Parts to fit Stryker® 900 Series Plaster Vacuums

**RPI Part #RXR105**


OEM Part # (No OEM Part # Available)

RECEPTACLE (HOSPITAL GRADE)

- Single hospital grade receptacle; Color: Grey
- Includes: Sealing gasket with PSA and mounting hardware as shown

Fits: Chassis

Models: 986, 987 & 996

**RPI Part #STL055**


OEM Part # (No OEM Part # Available)

LATCH

- 2 per package
- Low profile one piece molded polypropylene; Color: Black
- Includes: Latch, keeper, and mounting hardware

Fits: Vacuum Canister Housing and Lid

Models: 986, 987 & 996

**RPI Part #STK050**

OEM Part # (No OEM Part # Available)

POWER CORD KIT

- 12 ft. lg., 16/3 SJT, hospital grade 5/15, gray
- Includes: (1 pc) Strain Relief Bushing (RPI Part #RPB164), (2 pcs) QC Terminal (1/4 Female) (RPI Part #RPT472), and (1 pc) #6 Ring Terminal (RPI Part #RPT005)
- Also available: Strain Relief Bushing Pliers (RPI Part #RPT172)

Fits: Chassis

Models: 986

**RPI Part #STS047**

OEM Part # (No OEM Part # Available)

SWITCH

- Rocker switch (On-Off-On)
- 20A @ 125-277VAC; DPDT
- 1/4" tabs
- Panel Cutout: .830" x 1.450", .04" - .25" thk. panel

Fits: Vacuum Canister

Models: 986, 987 & 996

**RPI Part #RPB164**

OEM Part # (No OEM Part # Available)

STRAIN RELIEF BUSHING

- 8 per package
- Type: SR-6N3-4
- Accommodates: .325" - .360" dia. round cable.
- Mounting Hole: 5/8" dia. x 1/8" thk. panel max.
- Included in: Power Cord Kit (RPI Part #STK050)
- Also available: Strain Relief Bushing Pliers (RPI Part #RPT172)

Fits: Vacuum Housing

Models: 986, 987 & 996

x8

RPI Part #RXG107

OEM Part # (No OEM Part # Available)

GROMMET**(3/16" ID X 5/16" GROOVE)**

- 4 per package
- Material: Rubber, 55 Duro; Color: Black
- Included in: Motor Kit (RPI Part #STK048)

Fits: Motor Housing

Models: 986, 987 & 996



x4

RPI Part #RXG108

OEM Part # (No OEM Part # Available)

GROMMET**(1/4" ID X 11/32" GROOVE)**

- 4 per package
- Material: Rubber, 55 Duro; Color: Black

Fits: Vacuum Housing

Models: 986, 987 & 996



x4



Parts to fit Stryker® 900 Series Plaster Vacuums

**RPI Part #STS054****OEM Part #** (No OEM Part # Available)**LID SEAL**

- Size: 8.12" ID x 8.62" OD x 1/8" thk.
- Material: Grey closed cell vinyl foam with PSA backing
- Also available: Housing Seal (RPI Part #STS053)

Fits: Chassis**Models:** 986, 987 & 996**RPI Part #STS053****OEM Part #** (No OEM Part # Available)**HOUSING SEAL**

- Size: 1/2" wide x 28.0" lg. x 1/8" thk.
- Material: Black closed cell foam with PSA backing
- Includes: Installation template
- Also available: Lid Seal (RPI Part #STS054)

Fits: Chassis**Models:** 986, 987 & 996

Parts to fit Stryker® 900 Series Cast Cutters

**RPI Part #STA057****OEM Part #** (No OEM Part # Available)**POWER SUPPLY ASSEMBLY**

- Includes: Power Supply PCB (RPI Part #STB039), Power Cable Kit (RPI Part #STK037) and (6 pcs) Cord Strap (White 6") (RPI Part #RPS672)

Fits: Cast Saw Body**Models:** 940*For Synthetic & Plaster Casts!***RPI Part #STB056****OEM Part #940-23****CAST CUTTER BLADE (PiN)**

- For both plaster and synthetic casts
- Size: 2.36" oval blade
- Material: Plasma ion nitride treated stainless steel
- Color: Matte grey
- Also available: Blade Guard (RPI Part #STG042)

Fits: Head Assembly**Models:** 940 & 941*For Plaster Casts!***RPI Part #STB043****OEM Part #940-25****CAST CUTTER BLADE**

- For plaster casts
- Size: 2.36" oval blade
- Material: Stainless steel
- Color: Silver
- Also available: Blade Guard (RPI Part #STG042)

Fits: Head Assembly**Models:** 940 & 941*For Synthetic & Plaster Casts!***RPI Part #STB045****OEM Part #940-35****CAST CUTTER BLADE (TiN)**

- For both plaster and synthetic casts
- Size: 2.36" oval blade
- Material: Titanium nitride coated stainless steel
- Color: Gold
- Also available: Blade Guard (RPI Part #STG042)

Fits: Head Assembly**Models:** 940 & 941