PI Has The Best Tools To Make Your Job Easier!





Diagnostic Smart Kit® - For Servicing Tuttnauer® Automatic Autoclaves RPI Part #TUK108 (OEM Part #TEST-9)

• Calibrate/verify displayed pressure values and Temperature Sensor values.

• Check and monitor all test point values from the AJUNC or ANL boards.

• Includes: - Carrying Case (RPI Part #RPC476)

- Max Register Thermometer (RPI Part #RPT113)

- Pressure Gauge (Test) (RPI Part #TUG110)

- Ribbon Cable (RPI Part #TUC117)

- Simulator (PT-100) (RPI Part #TUT114)

- Test Point Board (RPI Part #TUB109)

- Trim Pot Adjustment Tool (RPI Part #RPT460)


- Wire Harness (No. 1) (RPI Part #TUH111)

- Wire Harness (No. 2) (RPI Part #TUH112)

- Wire Harness (No. 3) (RPI Part #TUH113) - Troubleshooting information and usage booklet

Models: EZ9, EZ10, EZ10k, 1730E/EK, 2340E/EA/EKA, 2540E/EA/EHS/EKA & 3870EHS





RPI Part #MIK074

· Check and adjust pressure and temperature

• Includes: - Cable Tie (High Temp) (RPI Part #RPT480)

- Carrying Case (RPI Part #RPC476)

- Max Register Thermometer (RPI Part #RPT113)

- Test Pressure Gauge (RPI Part #RPG461)

- Trim Pot Adjustment Tool (RPI Part #RPT460)

- Troubleshooting information and usage booklet

All parts in Kit are also Sold Separately!

Models:

M9-001 thru -019 with Serial # prefixes CZ, DA, DB, DX, DY, FD & OM;



M9D-001 thru -019 with Serial # prefixes FF & FG; M11-001 thru -019 with Serial # prefixes ES, ET, FP, FR & GB; M11D-001 thru -019 with Serial # prefixes GC & GD.

Field Service Smart Kit® - For Servicing SciCan® Cassette Autoclaves RPI Part #SCK028

• RPI exclusive service and repair kit!

• Includes: - Alignment Tool (RPI Part #SCT052)

- Calibration Solution (1 pouch) (RPI Part #RPS821)

- Carrying Case (RPI Part #RPC476)

- Cassette Seal Removal Tool (RPI Part #RPT372)

- Extension Test Cable (RPI Part #SCT026)

- Hex Balldriver (5/64") (RPI Part #RPT297)

- Hex Balldriver (9/64") (RPI Part #RPT836)

- Liquid Soap (RPI Part #RPS287)

- Pump Tester Bottle Kit (RPI Part #SCK024)

- Spanner Nut (RPI Part #RPT364)

- Start Switch (RPI Part #SCS021)

- Threadlocker 545 (RPI Part #RPA369)

- Water Quality Tester (RPI Part #RPT820)

- Troubleshooting information and usage booklet

Models: 2000 & 5000





replacement parts industries, inc.

More Tools To Make Your Job Easier!.

General Shop Aid Tools



RPI Part #RPT172

STRAIN RELIEF BUSHING PLIERS

• These hand pliers are adjustable for all bushings that assemble into a mounting hole up to 7/8" diameter!



RPI Part #RPT477

CABLE TIE TOOL

• This lightweight pocket sized tool is perfect for applying a wide variety of ties making for fast application and a clean cut-off!



RPI Part #RPT482

TERMINAL CRIMP TOOL

- Wire stripper #10 to #22 AWG
- · Crimps #10 to #22 insulated and non-insulated terminals
- This multi-function tool also includes pliers at tip, a wire loop former, wire cutter and can cut #4-40 to #10-32 screws!



RPI Part #RPT502

FUSE PULLER

 This handy tool is ideal for tight spaces and removes or installs glass, ceramic, or blade style fuses!



RPI Part #RPT460

TRIM POT ADJUSTMENT TOOL


- · Plastic with stainless steel blades
- Not only does this tool fit most trim pots, but it also functions as a small screwdriver and includes a shirt pocket clip!
- Also included in: Diagnostic Smart Kit[®] (RPI Part #TUK108) and Field Service Calibration Smart Kit® (RPI Part #MIK074)



RPI Part #RPT297

HEX BALLDRIVER (5/64")

- 5/64" Ball-point hex driver
- ~ 7" long steel tool with plastic handle
- Also included in: Field Service Smart Kit® (RPI Part #SCK028)



RPI Part #RPT836

HEX BALLDRIVER (9/64")

- 9/64" Ball-point hex driver
- ~ 10.50" long steel tool with plastic handle
- Also included in: Field Service Smart Kit® (RPI Part #SCK028)



RPI Part #RPT151

POZIDRIVE® SCREWDRIVER #1

• ~6-1/2" lg.



RPI Part #RPT152


POZIDRIVE® SCREWDRIVER #2

• ~8" lq.

RPI Part #RPT539

L-HEX WRENCH SET (METRIC)


- 9 piece set
- 1.5mm through 10mm
- Ball-point and chamfered ends




RPI Part #RPT540 L-HEX WRENCH SET

(INCH) • 13 piece set

- .050" through 3/8"
- Ball-point and chamfered ends




Dental Equipment Tools



RPI Part #RPT294

1/8" SLEEVE TOOL


- Use to install 1/8" sleeve clamps and uni-clamps on 1/8" OD tubing
- The tool provides 90° and 45° angles for ease of use



RPI Part #RPT295

1/4" SLEEVE TOOL


- Use to install 1/4" sleeve clamps and uni-clamps on 1/4" OD tubing
- The tool provides 90° and 45° angles for ease of use



RPI Part #RPT296

UNIVERSAL SLEEVE TOOL


- Use to install 1/8" and 1/4" sleeve clamps and uni-clamps on 1/8" and 1/4" OD tubing
- The tool provides 90° angles on both ends for ease of use



RPI Part #RPK581

O-RING REMOVAL TOOLS

· Pointed ends for hard to reach areas and blunt edges for larger surfaces



RPI Part #RPK481

O-RING INSTALLATION TOOLS

• Three sizes to accommodate O-Rings with ID of .029" through .208"



Dental Equipment Tools



RPI Part #ADT242

STEM ADJUSTMENT TOOL

- 5 per package; Fits: .125" hex socket
- Can be used in place of A-dec Handpiece Control or Flow Adjustment Keys
- Includes: 0-Ring (RPI Part #RP0486)



RPI Part #PCT615

GAS SPRING TOOL

- Great for Gas and Mechanical Spring adjustment!
- Our Tool has special markings right on the tool itself to help you make tension adjustments to Gas Springs. These markings indicate the maximum distance you should insert and turn the Spring Tool on Pelton & Crane ceiling lights or unit lights



RPI Part #BET029

ADJUSTMENT TOOL (SPRING TENSION)

 Perfect for adjusting spring tension and head position on Belmont and Healthco X-Calibur models



RPI Part #PCT685

SOCKET RETRACTOR TOOL


• For installation and removal of the light bulb used in the Pelton & Crane LF units



RPI Part #RPT705

SIGHT GLASS TOOL

 For easy removal and installation of sight glasses on Copeland Head lubricated compressors




RPI Part #RPT264

VALVE CORE REMOVAL/SYRINGE TOOL

 Syringe repair tool and removal tool for use with Valve Core (RPI Part #ADV001)


Sterilizer Tools



RPI Part #RPT501

DOUBLE ENDED SPANNER WRENCH

- .875" Hex x ~7" lg.
- Convenient to use when repairing difficult to reach solenoid valves in Instrumentation Laboratories analyzers and Midmark®, Pelton & Crane, and Scican sterilizers
- Includes: Spanner Nuts (RPI Part #'s ILS014 and RPT364) also sold separately
- Wow! Reversible and replaceable end attachments allow easy replacement of damaged pins!



RPI Part #RPT271

VALVE SEAT TOOL

 Specifically designed for installation and removal of valve seat and valve seat body on Pelton & Crane sterilizers!
 Models: Validator 8 (AB) & Validator 10 (AA)

RPI Part #RPT092

VALVE SEAT WRENCH

 Use to remove and install air valve seat and body

Models: Validator 8 (AB) & Validator 10 (AA)

RPI Part #RPT550 METAL TUBING BENDER


- Handles 3/16", 1/4", 5/16", and 3/8" OD metal tubing
- Easily create left, right, or offset bends up to 90°



RPI Part #RPT548

METAL TUBING CUTTER

- ~2" lg.
- 1/8" 5/8" OD tubing capacity
- Compact design makes this perfect to use where clearance around tubing is limited!
- Also available: replacement Tubing Cutter Wheel (RPI Part # RPT549)



RPI Part #RPT113

MAX REGISTER THERMOMETER

- Mercury
- Accuracy Rating: ±1 Division (±2°F)
- Dual Scale: 150°F to 350°F (65°C to 177°C)
- Perfect for temperature verification in sterilizers



RPI Part #RPK791 CLEANING KIT • Exclusively from RPI, this kit has the supplies you need to keep those sterilizers clean - from the outside casing to the small valves inside the device! • Includes all parts as shown - all parts also sold separately! RPB792 RPB793 RPB797 RPB797



replacement parts industries, inc.

20338 Corisco Street • Chatsworth, California • 91311

PRESORTED STANDARD U.S. POSTAGE PAID VAN NUYS, CA **PERMIT NO. 2399**


ADDRESS SERVICES REQUESTED



Imagine how RPI's Plastic and Metal Tubing Snakes can make your job easier!


There is a multitude of uses for RPI's Plastic and Metal Tubing Snakes (RPI Part #RPT620 & RPT629) that easily push and pull wire, cable, tubing, and fiber optic bundles behind walls, above ceilings, inside operatories, through crawl spaces, and under floors. Both tools are offered with a variety of attachments that allow you to customize them for each job!







PLASTIC & METAL TUBING SNAKES









PLASTIC TUBING SNAKE RPI Part #RPT620

- 15' of 5/32" dia. non-conductive fiberglass rod with PVC holder
- Includes: Bullet Tip (RPI Part #RPA627), Flexible Extension (RPI Part #RPA625), and Accessory Pouch (RPI Part #RPP634)
- One end has fixed 1/4" dia. brass tip with 3/16" dia. hole; Opposite end has a 1/4-20 male thread adapter for adding attachments

METAL TUBING SNAKE RPI Part #RPT629

- 12' of 1/8" dia. flexible cable
- Includes: Bullet Tip (RPI Part #RPA633), and Accessory Pouch (RPI Part #RPP634)
- One end has a fixed 3/16" dia. brass tip with 1/8" dia. hole; Opposite end has a #10-32 male thread adapter for adding attachments
- Adapter (RPI Part #RPA631) is required for use with the Pulling Eye Ring, Hook, and Magnet Attachments

Plastic Snake Attachments -



FIBER-OPTIC PULLER (5/16" dia.) RPI Part #RPA623



CONNECTOR (PLASTIC SNAKE) RPI Part #RPA621

Attachments for Both Snakes -

Other Attachments Also Available!



HOOK ATTACHMENT RPI Part #RPA624



PULLING EYE RING RPI Part #RPA626



MAGNET ATTACHMENT RPI Part #RPA622

Metal Snake Attachments



FIBER-OPTIC PULLER (1/4" dia.) RPI Part #RPA630



CONNECTOR (METAL SNAKE) RPI Part #RPA632



ADAPTER RPI Part #RPA631