


CONDENSATION COILS ... RPI HAS THE ONE YOU NEED!

Over the years, Tuttnauer has changed the configuration of their reservoirs by reversing the coil.

RPI has you covered – we offer all three styles!

RPI Part #TUC040
OEM Part #CU836101


Fits EZ9, EZ10,
1730E/EK & 1730M/MK



RPI Part #TUC041
OEM Part #CT836101 (Prior to 1993)

Fits 2340E/EK/M/MK
& 2540E/EK/M/MK
& Valueklave (1730MKV)

*Coil joins reservoir at the left rear and
vents towards the front of the machine*



RPI Part #TUC063
OEM Part #CT836101 (After 1993)

Fits 2340E/EK/M/MK
& 2540E/EK/M/MK

*Coil joins reservoir at the left front and
vents towards the rear of the machine*

