
New Parts For Utility Room Water Controls

RPI Part #WCE003 OEM Part #53035; PFM85240; 77000300: WFE-750

FILTER ELEMENT (5 MICRON)

- 2 per package
- Size: 2.50" OD x 9.88" lq.
- · Includes all parts as shown and PM Stickers
- · Recommended annual replacement for the Filter Element

Fits: 3/4" Filter Housing Assembly Models: See Models below



RPI Part #RX0065 0EM Part #53106: MRP70950: 77007010

0-RING

- 2 per package
- Size: 3.859" ID x 4.137" OD x .139" C/S
- Material: Nitrile (Buna-N), 70 Duro Fits: 3/4" Filter Housing Assembly

Models: See Models below



RPI Part #WCT007 0EM Part #MMS71173; 7700069

SPANNER WRENCH

• Wrench for 4.83" OD bowl Fits: 4.83" ID Filter Housing Assembly

Models: See Models below




Models: Air Techniques: 53020 & 53020-1; Apollo/Midmark: ASAHWFB34, ASAHWFB34B, ASAHWFB34C, ASAWCFB34, ASAWCFB34B, ASAWCFB34C, ASAWCFB34B, ASAWFB34B & ASAWFB34C; Matrx Medical: FBR-34R; Midmark: 77002146; Tech West: WB-750, WB-750-RC, WBF-750 & WBF-750-RC

RPI Part #WCE004 **0EM Part #64520006**

FILTER ELEMENT (5 MICRON)

- 2 per package
- Size: 2.50" OD x 9.88" lq.
- Includes all parts as shown and PM Stickers
- · Recommended annual replacement for the Filter Element

Fits: 3/4" Filter Housing Assembly Models: See Models below



RPI Part #RX0066

OEM Part # (No OEM Number Available)

0-RING

- 2 per package
- Size: 3.484" ID x 3.762" OD x .139" C/S
- Material: Nitrile (Buna-N), 70 Duro

Fits: 3/4" Filter Housing Assembly Models: See Models below



RPI Part #WCT008

OEM Part # (No OEM Number Available)

SPANNER WRENCH

• Wrench for 3.86" OD bowl Fits: 3.86" ID Filter Housing Assembly

Models: See Models below



Models: DentalEZ/Custom Air/Ramvac: WC110, WC210 & WC230


RPI Part #WCE010 OEM Part # 53153; PFM85241; 3 77000473; 64520012; WPE-100

FILTER ELEMENT (5 MICRON)

- 2 per package
- Size: 4.50" OD x 9.75" lg.
- · Includes all parts as shown and PM Stickers
- · Recommended annual replacement for the Filter Element

Fits: 1", 1-1/4", and 1-1/2" Filter Housing Assembly

Models: See Models below




RPI Part #RX0067 OEM Part #53172-1; MRP70951; 65977200; 64553032

0-RING

- 2 per package
- Size: 5.510" ID x 5.850" OD x .170" C/S
- Material: Nitrile (Buna-N), 70 Duro

Fits: 1", 1-1/4", and 1-1/2" Filter Housing Assembly

Models: See Models below



NOT SHOWN TO SIZE

RPI Part #WCT012 **0EM Part #77000063**

SPANNER WRENCH

• Wrench for 6.05" OD bowl Fits: 6.05" ID Filter Housing Assembly

Models: See Models below



Models: Air Techniques: 53170 & 53171; Apollo/Midmark®: ASAWFB100, ASAWFB100B, ASAWFB100C, ASAWFB125, ASAWFB125C, ASAWFB150 & ASAWFB125C; DentalEZ: WC-300, WC-310, WC-400 & WC-410; Matrx Medical: FBP-100R, FBP-125R & FBP-150R; Midmark®: 77002147, 77002148 & 77002149; Tech West: WBF-100 & WBF-100-RC;



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

SECTION

4

PAGE

1900

UPDATED

06/20

ENTAL WATER CONTROLS