






# New line of parts to fit Welch Allyn Wall Transformers!




**RPI Part #WAB001**  
**OEM Part #747284-501**  
**POWER SUPPLY BOARD (3.5V)**  
 • 3.5VDC  
 • Includes all parts as shown  
 • Also available: Position Control Knob (RPI Part #BUK043)  
**Fits:** Chassis  
**Model:** 74710



**RPI Part #WAK003**  
**OEM Part #747300-501** (Complete Handpiece Assembly)  
**COILED CORD KIT**  
 • Size: 2 Conductor #20 AWG x ~3 ft. lg.  
 • Includes all parts as shown  
 • Now you can replace just the Coiled Cord instead of the whole Handpiece Assembly - It's the RPI Advantage!  
 • Also available: Strain Relief Bushing Pliers (RPI Part #RPT172)  
**Fits:** Handpiece  
**Model:** 74710



**RPI Part #WAB004**  
**OEM Part #747145**  
**WALL BRACKET**  
 • Size: 4.16" lg. x 3.0" wd. x .06" thk.  
 • Material: Aluminum  
 • Includes: Mounting hardware  
**Fits:** Mounting Plate  
**Model:** 74710




**RPI Part #WAB007**  
**OEM Part #767035-501**  
**WALL BRACKET**  
 • Size: 6.0" lg. x 2.2" wd. x .06" thk.  
 • Material: Aluminum  
 • Includes: Mounting hardware  
**Fits:** Mounting Plate  
**Model:** 767



**RPI Part #BUK043**  
**OEM Part #747116-1**




**POSITION CONTROL KNOB**  
 • Material: Black plastic knob, aluminum face plate  
 • Size: .92" OD x .64" lg.  
 • Shaft Size: 1/4" dia.  
 • Includes: (2) installed set screws  
**Fits:** Control Panel  
**Model:** 74710



**RPI Part #WAS002**  
**OEM Part #747014-2**




**MICROSWITCH**  
 • Momentary  
 • 4A @ 250VAC; SPDT  
 • .058" dia. Pin Terminals  
**Fits:** Handle Mounting Bracket  
**Model:** 74710



**RPI Part #TUS009**  
**OEM Part #483311-2**




**POWER SWITCH**  
 • Rocker switch (On-Off)  
 • 20A @ 125/250VAC; DPST  
 • Color: Black body with illuminated green lens;  
 Legend: I/O  
 • .25" Tabs  
**Fits:** Main Housing  
**Model:** 74710



**RPI Part #WAS008**  
**OEM Part #236707-1**



**POWER SWITCH**  
 • Rocker Switch (ON - OFF)  
 • 16A @ 250VAC; DPST  
 • Color: Black body with illuminated green lens  
 • .25" Tabs  
**Fits:** Main Housing  
**Model:** 767



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@riparts.com](mailto:order@riparts.com) • WEBSITE [www.riparts.com](http://www.riparts.com)

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Welch Allyn equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Show 09/14




**replacement parts industries, inc.**

20338 Corisco Street • Chatsworth, California • 91311

PRESORTED STANDARD  
U.S. POSTAGE PAID  
VAN NUYS, CA  
PERMIT NO. 2399

ADDRESS SERVICES REQUESTED

## RPI Offers New Line of Parts to fit Welch Allyn Wall Transformers!



**RPI Part #CAK120**  
**OEM Part #747306-501** (Power Cord only)/711314 (Strain Relief only)

### POWER CORD KIT

- Length: 8 ft.
- Color: Gray
- 18/3 SJT, 60°C; 10A @ 125VAC
- Includes all parts as shown
- Also available: Strain Relief Bushing Pliers (RPI Part #RPT172)


**Fits:** Chassis  
**Model:** 74710

**RPI Part #RPC557**  
**OEM Part #761069**

### HOSPITAL GRADE POWER CORD

- Length: 8 ft.; Color: Gray
- 16/3 SJT, 60°C; 13A @ 125VAC
- Plug: NEMA 5-15P
- Connector: EN60320-C13

**Fits:** Power Input Module  
**Model:** 767




**RPI Part #RPF508**  
**OEM Part #236706-3111**

### FUSE (1/4A, 250V) - FAST ACTING

- 5 per package
- Type: Glass
- Size: 5mm dia. x 20mm lg.
- Included in: Fuse Kit (RPI Part #RPK520)

**Fits:** Power Supply Board  
**Model:** 767



## Lamps for Exam Heads


**RPI Part #LMP109**  
**OEM Part #WA-03100-U**

### LAMP - 6/pkg

- 3.5V, 2.52W
- Shape: T1.5; Xenon

**Fits:** Examination Head

**Models:** 20000 Otoscope, 20200 Otoscope, 20202 Otoscope, 21700 Otoscope, 25000 Otoscope, 25200 Otoscope, 25020 Otoscope, 73500 Handle Adapter, 28100 Tongue Blade Holder, 26530 Illuminator, 27000 Illuminator, 27050 Illuminator, 41100 Illuminator & 43300 Illuminator




**RPI Part #LMP110**  
**OEM Part #WA-03000-U**

### LAMP - 6/pkg

- 3.5V, 2.52W
- Shape: T1.5; Xenon
- Single contact 8mm flange

**Fits:** Examination Head

**Models:** 11600 Ophthalmoscope, 11605 Ophthalmoscope, 11610 Ophthalmoscope, 11710 Ophthalmoscope, 12400 Strabismoscope, 43700 Episcopes & 18000 Retinoscope




**RPI Part #LMP111**  
**OEM Part #WA-06500-U**

### LAMP - 6/pkg

- 3.5V, 2.73W
- Shape: TL1.5; Xenon

**Fits:** Examination Head

**Models:** 23810 Otoscope & 23820 Otoscope



**RPI Part #LMP112**  
**OEM Part #WA-03800-U**

### LAMP - 6/pkg

- 3.5V, 3.325W
- Shape: T1.5; Xenon

**Fits:** Examination Head

**Models:** 11800 Ophthalmoscope, 11810 Ophthalmoscope & 11820 Ophthalmoscope

