



## PATIENT CABLES & LEADWIRES

### RPI IS YOUR SOURCE FOR PATIENT CABLES & LEADWIRES!

#### PATIENT CABLES

- For ECG, EEG, monitoring and defibrillator equipment
- 3, 5, 7 and 10 leadwire sets available
- Snap or Pin style fixed lead patient ends
- Dual/Shielded, Din/Unshielded or Fixed lead patient ends
- Leadwire Adaptor Kit also available to adapt fixed banana plugs or pins to tab or snap electrodes
- Strain Relief protective feature molded into cable to extend life and integrity of the cables

#### ACQUISITION MODULE CABLES

- Cables to fit GE AM4, AM5 and CAM 14 Acquisition Modules
- 17 ft. durable coiled cords
- Includes molded strain reliefs and connector shields
- With over molded ferrite beads for superior noise reduction

#### HOLTER CABLES


- 3 Channel feature
- 5 and 7 Din/Unshielded patient ends
- Snap or Pinch style leadwire sets also available
- Strain Relief protective feature molded into cable to extend life and integrity of the cables

#### TELEMETRY CABLES

- 2, 3, 4 and 5 fixed lead patient ends
- Snap or Pinch style patient ends
- Strain Relief protective feature molded into cable to extend life and integrity of the cables

#### LEADWIRES

- Sets and Individual leadwires available in standard colors for easy identification
- Snap or Pinch style patient ends
- Dual/Shielded and Din/Unshielded leadwire cable ends
- Reuseable or disposable Leadwire sets available
- Sets and Individual leadwires for Veterinary use also available
- Strain Relief protective feature molded into leadwires to extend wires' life and integrity



FOR MORE INFORMATION ABOUT THE PARTS LISTED IN THIS FLYER, SEE RPI'S WEBSITE: [www.rpiparts.com](http://www.rpiparts.com)

**rpi** replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@rpiparts.com](mailto:order@rpiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)

The above parts are manufactured by Replacement Parts Industries, Inc. to fit various equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. REV C 01/16

# INDEX

## Acquisition, Patient, Holter and Telemetry Cables Page

Acquisition Module Cables to fit GE .....	4
Din/Unshielded Patient Cables .....	2
Dual/Shielded Patient Cables .....	3
Dual/Shielded Patient Cables to fit GE Multi-Link® .....	3
Fixed Patient Cables .....	1
Holter Cables .....	4
Telemetry Cables .....	4

## Patient Cable Leadwires and Leadwire Sets

Bonded Leadwire Sets to fit GE Multi-Link® .....	5
Combined Leadwire sets to fit GE Multi-Link® .....	5
Din/Unshielded Individual Leadwires .....	7
Din/Unshielded Leadwire Sets .....	6-7
Dual/Shielded Individual Leadwires .....	8
Dual/Shielded Leadwire Sets .....	8
Leadwire Sets to fit GE Multi-Link® .....	6

## Disposable Patient Cable Leadwire Sets

Din/Unshielded Leadwire Sets .....	9
Dual/Shielded Leadwire Sets .....	10

## Leadwire Accessories

Leadwire Adaptor Kit .....	1
Leadwire Combiners .....	10
Leadwire Organizers .....	10
Leadwire Shielded Accessory Kits .....	10

## Leadwires and Leadwire Sets for Veterinary Use


Din/Unshielded Individual Leadwires .....	11
Din/Unshielded Leadwire Sets .....	11
Dual/Shielded Individual Leadwires .....	11
Dual/Shielded Leadwire Sets .....	11

## OEM Model Cross Reference

Patient, Holter & Telemetry Cable OEM Cross Reference .....	13-16
---	-------

# FIXED PATIENT CABLES


## RPI Part #KCB002



### PATIENT CABLE - 3 LEAD FIXED SNAP

- Number of Leads: 3 • Length: 11 ft.


## RPI Part #KCB008



### PATIENT CABLE - 3 LEAD FIXED SNAP

- Number of Leads: 3 • Length: 11 ft.


## RPI Part #KCB036



### PATIENT CABLE - 3 LEAD FIXED SNAP

- OEM Part #11110-00029/3006218-02
- Number of Leads: 3 • Length: 10 ft.


## RPI Part #KEB012



### PATIENT CABLE - 5 LEAD FIXED SNAP

- Number of Leads: 5 • Length: 13 ft.


## RPI Part #KEB014



### PATIENT CABLE - 5 LEAD FIXED PIN

- Number of Leads: 5 • Length: 8 ft.
- Also available: Leadwire Adaptor Kit (RPI PART #LAK001)


## RPI Part #KEB030



### PATIENT CABLE - 5 LEAD FIXED BANANA

- Number of Leads: 5 • Length: 13 ft.
- Also available: Leadwire Adaptor Kit (RPI PART #LAK001)


## RPI Part #KKB002



### PATIENT CABLE - 10 LEAD FIXED BANANA

- Number of Leads: 10 • Length: 10 ft.
- Also available: Leadwire Adaptor Kit (RPI PART #LAK001)


## RPI Part #KKB003



### PATIENT CABLE - 10 LEAD FIXED BANANA

- Number of Leads: 10 • Length: 10 ft.
- Also available: Leadwire Adaptor Kit (RPI PART #LAK001)


## RPI Part #KKB006



### PATIENT CABLE - 10 LEAD FIXED BANANA

- Number of Leads: 10 • Length: 10 ft.
- Also available: Leadwire Adaptor Kit (RPI PART #LAK001)


## RPI Part #KKB007



### PATIENT CABLE - 10 LEAD FIXED BANANA

- Number of Leads: 10 • Length: 10 ft.
- Also available: Leadwire Adaptor Kit (RPI PART #LAK001)


## RPI Part #KKB008



### PATIENT CABLE - 10 LEAD FIXED BANANA

- Number of Leads: 10 • Length: 10 ft.
- Also available: Leadwire Adaptor Kit (RPI PART #LAK001)


## RPI Part #KKB010



### PATIENT CABLE - 10 LEAD FIXED PINCH

- OEM Part #60-00180-01/60-00184-01
- Number of Leads: 10 • Length: 18 ft.


## RPI Part #KKB011



### PATIENT CABLE - 10 LEAD FIXED SNAP

- OEM Part #60-00181-01/60-00185-01
- Number of Leads: 10 • Length: 18 ft.

## RPI Part #KKB012



### PATIENT CABLE - 10 LEAD FIXED BANANA

- OEM Part #223 418 09
- Number of Leads: 10 • Length: 9 ft.
- Also available: Leadwire Adaptor Kit (RPI PART #LAK001)

## RPI Part #LAK001

### LEADWIRE ADAPTOR KIT

- Includes 10 adaptors individually labeled with leadwire identifiers
- Adapts banana plugs or pins to tab or snap electrodes
- Widest contact area for superior trace quality
- Flat, wide bodied design helps prevent twisting
- Durable construction - built to last




Includes 10 adaptors individually labeled with leadwire identifiers





# DUAL/SHIELDED PATIENT CABLES


## RPI Part #KCC003



### PATIENT CABLE - 3 LEAD DUAL

- Number of Leads: 3 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM014 & LKM017), Disposable Leadwire Sets (RPI Part #'s LDM305 & LDM312), and Veterinary Leadwire Set (RPI Part #LVM705)


## RPI Part #KCC009



### PATIENT CABLE - 3 LEAD DUAL

- Number of Leads: 3 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM014 & LKM017), Disposable Leadwire Sets (RPI Part #'s LDM305 & LDM312), and Veterinary Leadwire Set (RPI Part #LVM705)


## RPI Part #KCC012



### PATIENT CABLE - 3 LEAD DUAL

- Number of Leads: 3 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM014 & LKM017), Disposable Leadwire Sets (RPI Part #'s LDM305 & LDM312), and Veterinary Leadwire Set (RPI Part #LVM705)


## RPI Part #KCC021



### PATIENT CABLE - 3 LEAD DUAL

- Number of Leads: 3 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM014 & LKM017), Disposable Leadwire Sets (RPI Part #'s LDM305 & LDM312), and Veterinary Leadwire Set (RPI Part #LVM705)


## RPI Part #KCC034



### PATIENT CABLE - 3 LEAD DUAL OEM Part #545302

- Number of Leads: 3 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM014 & LKM017), Disposable Leadwire Sets (RPI Part #'s LDM305 & LDM312), and Veterinary Leadwire Set (RPI Part #LVM705)


## RPI Part #KEC003



### PATIENT CABLE - 5 LEAD DUAL

- Number of Leads: 5 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM032 & LKM035), Disposable Leadwire Sets (RPI Part #'s LDM307 & LDM314), and Veterinary Leadwire Set (RPI Part #LVM707)


## RPI Part #KEC007



### PATIENT CABLE - 5 LEAD DUAL

- Number of Leads: 5 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM032 & LKM035), Disposable Leadwire Sets (RPI Part #'s LDM307 & LDM314), and Veterinary Leadwire Set (RPI Part #LVM707)


## RPI Part #KEC010



### PATIENT CABLE - 5 LEAD DUAL

- Number of Leads: 5 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM032 & LKM035), Disposable Leadwire Sets (RPI Part #'s LDM307 & LDM314), and Veterinary Leadwire Set (RPI Part #LVM707)


## RPI Part #KEC016



### PATIENT CABLE - 5 LEAD DUAL

- Number of Leads: 5 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM032 & LKM035), Disposable Leadwire Sets (RPI Part #'s LDM307 & LDM314), and Veterinary Leadwire Set (RPI Part #LVM707)

## RPI Part #KEC022




### PATIENT CABLE - 5 LEAD DUAL OEM Part #545303

- Number of Leads: 5 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM032 & LKM035), Disposable Leadwire Sets (RPI Part #'s LDM307 & LDM314), and Veterinary Leadwire Set (RPI Part #LVM707)

# DUAL/SHIELDED PATIENT CABLES TO FIT GE MULTI-LINK®


## RPI Part #KCD033



### PATIENT CABLE - 3 LEAD MULTI-LINK OEM Part #1554AA0

- Number of Leads: 3 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM052 & LKM053)

## RPI Part #KED021



### PATIENT CABLE - 5 LEAD MULTI-LINK OEM Part #2017003-001/412931-001

- Number of Leads: 5 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM052, LKM053, LKM054, LKM055, LKM056 & LKM057)


# ACQUISITION MODULE CABLES TO FIT GE

## RPI Part #GEC020

**ACQUISITION MODULE CABLE**  
**OEM Part #700044-202/700044-204**

- Length: ~17 ft. overall length

**Fits:** AM4 Acquisition Module and AM5 Acquisition Module




## RPI Part #GEC021

**ACQUISITION MODULE CABLE**  
**OEM Part #2016560-002/2016560-003/700657-002/700657-003**


- Length: ~17 ft. overall length

**Fits:** CAM 14 Acquisition Module




## HOLTER CABLES

### RPI Part #KEA004



**HOLTER CABLE - 5 LEAD 3 CHANNEL**

- Number of Leads: 5 • Length: 30"
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM026, LKM029 & LKM030), Disposable Leadwire Sets (RPI Part #'s LDM303 & LDM310), and Veterinary Leadwire Set (RPI Part #LVM703)




### RPI Part #KGA001




**HOLTER CABLE - 7 LEAD 3 CHANNEL**

- Number of Leads: 7 • Length: 40"
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM038 & LKM041), Disposable Leadwire Sets (RPI Part #'s LDM304 & LDM311), and Veterinary Leadwire Set (RPI Part #LVM704)



### RPI Part #KEB020



**HOLTER CABLE - 5 LEAD FIXED SNAP**  
**OEM Part #3-370-0001**


- Number of Leads: 5 • Length: 43"



## TELEMETRY CABLES

Accessories  
Also Available  
See Page 10


### RPI Part #KBB001



**TELEMETRY CABLE - 2 LEAD SNAP**

- Number of Leads: 2 • Length: 32"


### RPI Part #KCB028



**TELEMETRY CABLE - 3 LEAD SNAP**

- Number of Leads: 3 • Length: 27"
- Also available: 3 Lead Organizer (RPI Part #LAR004)

### RPI Part #KCB029




**TELEMETRY CABLE - 3 LEAD SNAP**

- Number of Leads: 3 • Length: 30" • Includes: Adaptor Block
- Also available: 3 Lead Organizer (RPI Part #LAR004)

## TELEMETRY CABLES (continued)


### RPI Part #KCB031



**TELEMETRY CABLE - 3 LEAD SNAP**

- Number of Leads: 3 • Length: 30"
- Also available: 3 Lead Organizer (RPI Part #LAR004)


### RPI Part #KCB032



**TELEMETRY CABLE - 3 LEAD PINCH**

- Number of Leads: 3 • Length: 30"
- Also available: 3 Lead Organizer (RPI Part #LAR004)


### RPI Part #KDB002



**TELEMETRY CABLE - 4 LEAD SNAP**

- Number of Leads: 4 • Length: 27"
- Also available: 5 Lead Organizer (RPI Part #LAR005)


### RPI Part #KEB015



**TELEMETRY CABLE - 5 LEAD SNAP**

- Number of Leads: 5 • Length: 24"
- Also available: 5 Lead Organizer (RPI Part #LAR005)


### RPI Part #KEB018



**TELEMETRY CABLE - 5 LEAD SNAP**

- Number of Leads: 5 • Length: 30"
- Also available: 5 Lead Organizer (RPI Part #LAR005)

### RPI Part #KEB019




**TELEMETRY CABLE - 5 LEAD PINCH**

- Number of Leads: 5 • Length: 30"
- Also available: 5 Lead Organizer (RPI Part #LAR005)

## BONDED LEADWIRE SETS TO FIT GE MULTI-LINK®

### RPI Part #LKM052




#### 3 LEADWIRE SET - 30" MULTI-LINK/SNAP

##### OEM Part #411203-001

- Number of Leads: 3 bonded leadwires • Length: 30"
- Includes: 3 Lead Organizer (RPI Part #LAR004)
- Also available: Bonded Leadwire Sets (RPI Part #'s LKM053, LKM054 & LKM055) and Combined Leadwire Sets (RPI Part #'s LKM056 & LKM057)

**Fits:** GE Multi-Link® Cable, Patient Cable - 3 Lead Multi-Link® (RPI Part #KCD033) and Patient Cable - 5 Lead Multi-Link® (RPI Part #KED021)

### RPI Part #LKM053




#### 3 LEADWIRE SET - 30" MULTI-LINK/PINCH

##### OEM Part #412682-001

- Number of Leads: 3 bonded leadwires • Length: 30"
- Includes: 3 Lead Organizer (RPI Part #LAR004)
- Also available: Bonded Leadwire Sets (RPI Part #'s LKM052, LKM054 & LKM055) and Combined Leadwire Sets (RPI Part #'s LKM056 & LKM057)

**Fits:** GE Multi-Link® Cable, Patient Cable - 3 Lead Multi-Link® (RPI Part #KCD033) and Patient Cable - 5 Lead Multi-Link® (RPI Part #KED021)

### RPI Part #LKM054




#### 5 LEADWIRE SET - 30" MULTI-LINK/SNAP

##### OEM Part #411202-001

- Number of Leads: 5 bonded leadwires • Length: 30"
- Includes: 5 Lead Organizer (RPI Part #LAR005)
- Also available: Bonded Leadwire Sets (RPI Part #'s LKM052, LKM053 & LKM055) and Combined Leadwire Sets (RPI Part #'s LKM056 & LKM057)

**Fits:** GE Multi-Link® Cable and Patient Cable - 5 Lead Multi-Link® (RPI Part #KED021)

### RPI Part #LKM055



#### 5 LEADWIRE SET - 30" MULTI-LINK/PINCH


##### OEM Part #412681-001

- Number of Leads: 5 bonded leadwires • Length: 30"
- Includes: 5 Lead Organizer (RPI Part #LAR005)
- Also available: Bonded Leadwire Sets (RPI Part #'s LKM052, LKM053 & LKM054) and Combined Leadwire Sets (RPI Part #'s LKM056 & LKM057)

**Fits:** GE Multi-Link® Cable and Patient Cable - 5 Lead Multi-Link® (RPI Part #KED021)

## COMBINED LEADWIRE SETS TO FIT GE MULTI-LINK®

### RPI Part #LKM056




#### 5 LEADWIRE SET - 30" MULTI-LINK/SNAP

##### OEM Part #411200-001

- Number of Leads: 5 individual leadwires with combiner • Length: 30"
- Includes: 5 Lead Organizer (RPI Part #LAR005) and combiner
- Also available: Bonded Leadwire Sets (RPI Part #'s LKM052, LKM053, LKM054 & LKM055) and Combined Leadwire Set (RPI Part #LKM057)

**Fits:** GE Multi-Link® Cable and Patient Cable - 5 Lead Multi-Link® (RPI Part #KED021)

### RPI Part #LKM057




#### 5 LEADWIRE SET - 30" MULTI-LINK/PINCH

##### OEM Part #414556-001

- Number of Leads: 5 individual leadwires with combiner • Length: 30"
- Includes: 5 Lead Organizer (RPI Part #LAR005) and combiner
- Also available: Bonded Leadwire Sets (RPI Part #'s LKM052, LKM053, LKM054 & LKM055) and Combined Leadwire Set (RPI Part #LKM056)

**Fits:** GE Multi-Link® Cable and Patient Cable - 5 Lead Multi-Link® (RPI Part #KED021)

### RPI Part #LKM058




#### 6 LEADWIRE SET - 30" MULTI-LINK/SNAP

##### OEM Part #421930-001

- Number of Leads: 6 individual leadwires with combiner • Length: 30"
- Includes: (2) 3 Lead Organizer (RPI Part #LAR004) and 6 Lead Combiner
- Also available: Combined Leadwire Set (RPI Part #LKM059)

### RPI Part #LKM059



#### 6 LEADWIRE SET - 30" MULTI-LINK/PINCH

##### OEM Part #421932-001

- Number of Leads: 6 individual leadwires with combiner • Length: 30"
- Includes: (2) 3 Lead Organizer (RPI Part #LAR004) and 6 Lead Combiner
- Also available: Combined Leadwire Set (RPI Part #LKM058)

# LEADWIRE SETS TO FIT GE MULTI-LINK®


## RPI Part #LKM049

### 10 LEADWIRE SET - CAM 14/BANANA

#### OEM Part #420101-002

- Number of Leads: 10
- Length: 27": Red (V1), Yellow (V2), Green (V3), Blue (V4), Orange (V5), Violet (V6); and, 37": Green (RL), White (RA), Black (LA) & Red (LL)
- Includes: (2) 5 Lead Organizer (RPI Part #LAR005)
- Also available: 10 Leadwire Set - CAM 14/Pinch (RPI Part #LKM050) & 10 Leadwire Set - CAM 14/Snap (RPI Part #LKM051)

Fits: CAM 14 Module




## RPI Part #LKM050

### 10 LEADWIRE SET - CAM 14/PINCH

#### OEM Part #420101-002

- Number of Leads: 10
- Length: 27": Red (V1), Yellow (V2), Green (V3), Blue (V4), Orange (V5), Violet (V6); and, 37": Green (RL), White (RA), Black (LA) & Red (LL)
- Includes: (2) 5 Lead Organizer (RPI Part #LAR005)
- Also available: 10 Leadwire Set - CAM 14/Banana (RPI Part #LKM049) & 10 Leadwire Set - CAM 14/Snap (RPI Part #LKM051)

Fits: CAM 14 Module




## RPI Part #LKM051

### 10 LEADWIRE SET - CAM 14/SNAP

#### OEM Part #420101-002

- Number of Leads: 10
- Length: 27": Red (V1), Yellow (V2), Green (V3), Blue (V4), Orange (V5), Violet (V6); and, 37": Green (RL), White (RA), Black (LA) & Red (LL)
- Includes: (2) 5 Lead Organizer (RPI Part #LAR005)
- Also available: 10 Leadwire Set - CAM 14/Banana (RPI Part #LKM049) & 10 Leadwire Set - CAM 14/Pinch (RPI Part #LKM050)

Fits: CAM 14 Module




## DIN/UNSHIELDED LEADWIRE SETS

Accessories  
Also Available  
See Page 10

## RPI Part #LKM002

### 3 LEADWIRE SET - 24" DIN/PINCH

- Number of Leads: 3 • Length: 24"
- Also available: Disposable Leadwire Set (RPI Part #LDM308) & 3 Lead Organizer (RPI Part #LAR004)



## RPI Part #LKM005

### 3 LEADWIRE SET - 24" DIN/SNAP

- Number of Leads: 3 • Length: 24"
- Also available: Disposable Leadwire Set (RPI Part #LDM301) & 3 Lead Organizer (RPI Part #LAR004)



## RPI Part #LKM008

### 3 LEADWIRE SET - 24" DIN/PINCH

- Number of Leads: 3 • Length: 24"
- Also available: Disposable Leadwire Set (RPI Part #LDM309) & 3 Lead Organizer (RPI Part #LAR004)



## RPI Part #LKM011

### 3 LEADWIRE SET - 24" DIN/SNAP

- Number of Leads: 3 • Length: 24"
- Also available: Disposable Leadwire Set (RPI Part #LDM302) & 3 Lead Organizer (RPI Part #LAR004)






## DIN/UNSHIELDED LEADWIRE SETS (continued)

Accessories  
Also Available  
See Page 10


### RPI Part #LKM026



#### 5 LEADWIRE SET - 24" DIN/PINCH

- Number of Leads: 5 • Length: 24"
- Also available: Disposable Leadwire Set (RPI Part #LDM310) & 5 Lead Organizer (RPI Part #LAR005)


### RPI Part #LKM029



#### 5 LEADWIRE SET - 24" DIN/SNAP

- Number of Leads: 5 • Length: 24"
- Also available: 24" Disposable Leadwire Set (RPI Part #LDM303) & 5 Lead Organizer (RPI Part #LAR005)


### RPI Part #LKM030



#### 5 LEADWIRE SET - 40" DIN/SNAP

- Number of Leads: 5 • Length: 40"
- Also available: 24" Disposable Leadwire Set (RPI Part #LDM303) & 5 Lead Organizer (RPI Part #LAR005)


### RPI Part #LKM038



#### 7 LEADWIRE SET - 24" DIN/PINCH

- Number of Leads: 7 • Length: 24"
- Also available: Disposable Leadwire Set (RPI Part #LDM311), 3 Lead Organizer (RPI Part #LAR004) & 5 Lead Organizer (RPI Part #LAR005)


### RPI Part #LKM041



#### 7 LEADWIRE SET - 24" DIN/SNAP


- Number of Leads: 7 • Length: 24"
- Also available: Disposable Leadwire Set (RPI Part #LDM304), 3 Lead Organizer (RPI Part #LAR004) & 5 Lead Organizer (RPI Part #LAR005)

## DIN/UNSHIELDED INDIVIDUAL LEADWIRES



- Equipment End: Din/Unshielded, .060" Safety Socket
- Patient End: Pinch
- Length: 24"

DESCRIPTION	LEADWIRE COLOR
<b>RPI Part #LPB002</b> LEADWIRE BLACK 24" DIN/PINCH	Black
<b>RPI Part #LPE002</b> LEADWIRE ORANGE 24" DIN/PINCH	Orange
<b>RPI Part #LPG002</b> LEADWIRE GREEN 24" DIN/PINCH	Green
<b>RPI Part #LPN002</b> LEADWIRE BROWN 24" DIN/PINCH	Brown
<b>RPI Part #LPR002</b> LEADWIRE RED 24" DIN/PINCH	Red
<b>RPI Part #LPU002</b> LEADWIRE BLUE 24" DIN/PINCH	Blue
<b>RPI Part #LPW002</b> LEADWIRE WHITE 24" DIN/PINCH	White



- Equipment End: Din/Unshielded, .060" Safety Socket
- Patient End: Snap
- Length: 24"


DESCRIPTION	LEADWIRE COLOR
<b>RPI Part #LSB002</b> LEADWIRE BLACK 24" DIN/SNAP	Black
<b>RPI Part #LSE002</b> LEADWIRE ORANGE 24" DIN/SNAP	Orange
<b>RPI Part #LSG002</b> LEADWIRE GREEN 24" DIN/SNAP	Green
<b>RPI Part #LSN002</b> LEADWIRE BROWN 24" DIN/SNAP	Brown
<b>RPI Part #LSR002</b> LEADWIRE RED 24" DIN/SNAP	Red
<b>RPI Part #LSU002</b> LEADWIRE BLUE 24" DIN/SNAP	Blue
<b>RPI Part #LSW002</b> LEADWIRE WHITE 24" DIN/SNAP	White

Depending on equipment model, these leadwires are also available as sets – see *DIN/UNSHIELDED LEADWIRE SETS* on pages 6-7.

# DUAL/SHIELDED LEADWIRE SETS

Accessories  
Also Available  
See Page 10


## RPI Part #LKM014



### 3 LEADWIRE SET - 40" DUAL/PINCH

- Number of Leads: 3 • Length: 40"
- Includes: 3 Lead Shielded Accessory Kit (RPI Part #LAK006)
- Also available: Disposable Leadwire Set (RPI Part #LDM312)


## RPI Part #LKM017



### 3 LEADWIRE SET - 40" DUAL/SNAP

- Number of Leads: 3 • Length: 40"
- Includes: 3 Lead Shielded Accessory Kit (RPI Part #LAK006)
- Also available: Disposable Leadwire Set (RPI Part #LDM305)


## RPI Part #LKM020



### 3 LEADWIRE SET - 40" DUAL/PINCH

- Number of Leads: 3 • Length: 40"
- Includes: 3 Lead Shielded Accessory Kit (RPI Part #LAK006)
- Also available: Disposable Leadwire Set (RPI Part #LDM313)


## RPI Part #LKM023



### 3 LEADWIRE SET - 40" DUAL/SNAP

- Number of Leads: 3 • Length: 40"
- Includes: 3 Lead Shielded Accessory Kit (RPI Part #LAK006)
- Also available: Disposable Leadwire Set (RPI Part #LDM306)


## RPI Part #LKM032



### 5 LEADWIRE SET - 40" DUAL/PINCH

- Number of Leads: 5 • Length: 40"
- Includes: 5 Lead Shielded Accessory Kit (RPI Part #LAK007)
- Also available: Disposable Leadwire Set (RPI Part #LDM314)


## RPI Part #LKM035



### 5 LEADWIRE SET - 40" DUAL/SNAP


- Number of Leads: 5 • Length: 40"
- Includes: 5 Lead Shielded Accessory Kit (RPI Part #LAK007)
- Also available: Disposable Leadwire Set (RPI Part #LDM307)

# DUAL/SHIELDED INDIVIDUAL LEADWIRES



- Equipment End: Dual/Shielded
- Patient End: Pinch
- Length: 40"

DESCRIPTION	LEADWIRE COLOR
<b>RPI Part #LPB005</b> LEADWIRE BLACK 40" DUAL/PINCH	Black
<b>RPI Part #LPG005</b> LEADWIRE GREEN 40" DUAL/PINCH	Green
<b>RPI Part #LPN005</b> LEADWIRE BROWN 40" DUAL/PINCH	Brown
<b>RPI Part #LPR005</b> LEADWIRE RED 40" DUAL/PINCH	Red
<b>RPI Part #LPW005</b> LEADWIRE WHITE 40" DUAL/PINCH	White



- Equipment End: Dual/Shielded
- Patient End: Snap
- Length: 40"


DESCRIPTION	LEADWIRE COLOR
<b>RPI Part #LSB005</b> LEADWIRE BLACK 40" DUAL/SNAP	Black
<b>RPI Part #LSG005</b> LEADWIRE GREEN 40" DUAL/SNAP	Green
<b>RPI Part #LSN005</b> LEADWIRE BROWN 40" DUAL/SNAP	Brown
<b>RPI Part #LSR005</b> LEADWIRE RED 40" DUAL/SNAP	Red
<b>RPI Part #LSW005</b> LEADWIRE WHITE 40" DUAL/SNAP	White

Depending on equipment model, these leadwires are also available as sets – see *DUAL/SHIELDED LEADWIRE SETS* at the top of this page.

# DISPOSABLE: DIN/UNSHIELDED LEADWIRE SETS

Accessories  
Also Available  
See Page 10


## RPI Part #LDM301



### 3 LEADWIRE SET DISPOSABLE 24" DIN/SNAP

- Number of Leads: 3 • Length: 24" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM005) & 3 Lead Organizer (RPI Part #LAR004)


## RPI Part #LDM302



### 3 LEADWIRE SET DISPOSABLE 24" DIN/SNAP

- Number of Leads: 3 • Length: 24" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM011) & 3 Lead Organizer (RPI Part #LAR004)


## RPI Part #LDM303



### 5 LEADWIRE SET DISPOSABLE 24" DIN/SNAP

- Number of Leads: 5 • Length: 24" • 10 sets per package
- Also available: Reusable 24" Leadwire Set (RPI Part #LKM029), Reusable 40" Leadwire Set (RPI Part #LKM030) & 5 Lead Organizer (RPI Part #LAR005)


## RPI Part #LDM304



### 7 LEADWIRE SET DISPOSABLE 24" DIN/SNAP

- Number of Leads: 7 • Length: 24" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM041)


## RPI Part #LDM308



### 3 LEADWIRE SET DISPOSABLE 24" DIN/PINCH

- Number of Leads: 3 • Length: 24" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM002) & 3 Lead Organizer (RPI Part #LAR004)


## RPI Part #LDM309



### 3 LEADWIRE SET DISPOSABLE 24" DIN/PINCH

- Number of Leads: 3 • Length: 24" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM008) & 3 Lead Organizer (RPI Part #LAR004)


## RPI Part #LDM310



### 5 LEADWIRE SET DISPOSABLE 24" DIN/PINCH

- Number of Leads: 5 • Length: 24" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM026) & 5 Lead Organizer (RPI Part #LAR005)

## RPI Part #LDM311




### 7 LEADWIRE SET DISPOSABLE 24" DIN/PINCH

- Number of Leads: 7 • Length: 24" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM038)

# DISPOSABLE: DUAL/SHIELDED LEADWIRE SETS

Accessories  
Also Available  
See Below


## RPI Part #LDM305



### 3 LEADWIRE SET DISPOSABLE 30" DUAL/SNAP

- Number of Leads: 3 • Length: 30" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM017) & 3 Lead Shielded Accessory Kit (RPI Part #LAK006)


## RPI Part #LDM306



### 3 LEADWIRE SET DISPOSABLE 30" DUAL/SNAP

- Number of Leads: 3 • Length: 30" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM023) & 3 Lead Shielded Accessory Kit (RPI Part #LAK006)


## RPI Part #LDM307



### 5 LEADWIRE SET DISPOSABLE 30" DUAL/SNAP

- Number of Leads: 5 • Length: 30" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM035) & 5 Lead Shielded Accessory Kit (RPI Part #LAK007)


## RPI Part #LDM312



### 3 LEADWIRE SET DISPOSABLE 30" DUAL/PINCH

- Number of Leads: 3 • Length: 30" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM014) & 3 Lead Shielded Accessory Kit (RPI Part #LAK006)


## RPI Part #LDM313



### 3 LEADWIRE SET DISPOSABLE 30" DUAL/PINCH

- Number of Leads: 3 • Length: 30" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM020) & 3 Lead Shielded Accessory Kit (RPI Part #LAK006)

## RPI Part #LDM314



### 5 LEADWIRE SET DISPOSABLE 30" DUAL/PINCH

- Number of Leads: 5 • Length: 30" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM032) & 5 Lead Shielded Accessory Kit (RPI Part #LAK007)

## ACCESSORIES

DESCRIPTION	ACCESSORY
<b>RPI Part #LAC002</b> <b>3 LEAD COMBINER</b> <ul style="list-style-type: none"> <li>• Number of Leads: 3</li> <li>• Use with Dual/Shielded Leadwires only</li> </ul> <b>Also included in:</b> <ul style="list-style-type: none"> <li>• 3 Lead Shielded Accessory Kit (RPI Part #LAK006)</li> </ul>	
<b>RPI Part #LAR004</b> <b>3 LEAD ORGANIZER</b> <ul style="list-style-type: none"> <li>• Number of Leads: 3</li> </ul> <b>Also included in:</b> <ul style="list-style-type: none"> <li>• 3 Lead Shielded Accessory Kit (RPI Part #LAK006)</li> </ul>	
<b>RPI Part #LAK006</b> <b>3 LEAD SHIELDED ACCESSORY KIT</b> <ul style="list-style-type: none"> <li>• Number of Leads: 3</li> <li>• Shielded Accessory Kit</li> <li>• Use with Dual/Shielded Leadwires only</li> <li>• Includes all parts as shown</li> </ul>	


DESCRIPTION	ACCESSORY
<b>RPI Part #LAC003</b> <b>5 LEAD COMBINER</b> <ul style="list-style-type: none"> <li>• Number of Leads: 5</li> <li>• Use with Dual/Shielded Leadwires only</li> </ul> <b>Also included in:</b> <ul style="list-style-type: none"> <li>• 5 Lead Shielded Accessory Kit (RPI Part #LAK007)</li> </ul>	
<b>RPI Part #LAR005</b> <b>5 LEAD ORGANIZER</b> <ul style="list-style-type: none"> <li>• Number of Leads: 5</li> </ul> <b>Also included in:</b> <ul style="list-style-type: none"> <li>• 5 Lead Shielded Accessory Kit (RPI Part #LAK007)</li> </ul>	
<b>RPI Part #LAK007</b> <b>5 LEAD SHIELDED ACCESSORY KIT</b> <ul style="list-style-type: none"> <li>• Number of Leads: 5</li> <li>• Shielded Accessory Kit</li> <li>• Use with Dual/Shielded Leadwires only</li> <li>• Includes all parts as shown</li> </ul>	



## VETERINARY: DIN/UNSHIELDED LEADWIRE SETS

Accessories  
Also Available  
See Page 10


### RPI Part #LVM701



#### 3 LEADWIRE SET - 24" DIN/ALLIGATOR

- Reuseable
- Number of Leads: 3 • Length: 24"
- Recommended for veterinary use only
- Also available: 3 Lead Organizer (RPI Part #LAR004)


### RPI Part #LVM702



#### 3 LEADWIRE SET - 24" DIN/ALLIGATOR

- Reuseable
- Number of Leads: 3 • Length: 24"
- Recommended for veterinary use only
- Also available: 3 Lead Organizer (RPI Part #LAR004)


### RPI Part #LVM703



#### 5 LEADWIRE SET - 24" DIN/ALLIGATOR

- Reuseable
- Number of Leads: 5 • Length: 24"
- Recommended for veterinary use only
- Also available: 5 Lead Organizer (RPI Part #LAR005)

### RPI Part #LVM704




#### 7 LEADWIRE SET - 24" DIN/ALLIGATOR

- Reuseable
- Number of Leads: 7 • Length: 24"
- Recommended for veterinary use only
- Also available: 3 Lead Organizer (RPI Part #LAR004)  
& 5 Lead Organizer (RPI Part #LAR005)

## VETERINARY: DUAL/SHIELDED LEADWIRE SETS

Accessories  
Also Available  
See Page 10


### RPI Part #LVM705



#### 3 LEADWIRE SET - 30" DUAL/ALLIGATOR

- Reuseable
- Number of Leads: 3 • Length: 30"
- Recommended for veterinary use only
- Also available: 3 Lead Shielded Accessory Kit (RPI Part #LAK006)

### RPI Part #LVM707



#### 5 LEADWIRE SET - 30" DUAL/ALLIGATOR

- Reuseable
- Number of Leads: 5 • Length: 30"
- Recommended for veterinary use only
- Also available: 5 Lead Shielded Accessory Kit (RPI Part #LAK007)

### RPI Part #LVM706




#### 3 LEADWIRE SET - 30" DUAL/ALLIGATOR








- Reuseable
- Number of Leads: 3 • Length: 30"
- Recommended for veterinary use only
- Also available: 3 Lead Shielded Accessory Kit (RPI Part #LAK006)

## VETERINARY: INDIVIDUAL LEADWIRES


## DIN/UNSHIELDED INDIVIDUAL








- Equipment End: Din/Unshielded, .060 Safety Socket
- Patient End: Alligator Clip
- Length: 24"

DESCRIPTION	LEADWIRE COLOR
<b>RPI Part #LGB002</b> LEADWIRE BLACK 24" DIN/ALLIGATOR	Black 
<b>RPI Part #LGE002</b> LEADWIRE ORANGE 24" DIN/ALLIGATOR	Orange 
<b>RPI Part #LGG002</b> LEADWIRE GREEN 24" DIN/ALLIGATOR	Green 
<b>RPI Part #LGN002</b> LEADWIRE BROWN 24" DIN/ALLIGATOR	Brown 
<b>RPI Part #LGR002</b> LEADWIRE RED 24" DIN/ALLIGATOR	Red 
<b>RPI Part #LGU002</b> LEADWIRE BLUE 24" DIN/ALLIGATOR	Blue 
<b>RPI Part #LGW002</b> LEADWIRE WHITE 24" DIN/ALLIGATOR	White 

## DUAL/SHIELDED INDIVIDUAL



- Equipment End: Dual/Shielded
- Patient End: Alligator Clip
- Length: 30"

DESCRIPTION	LEADWIRE COLOR
<b>RPI Part #LGB004</b> LEADWIRE BLACK 30" DUAL/ALLIGATOR	Black 
<b>RPI Part #LGG004</b> LEADWIRE GREEN 30" DUAL/ALLIGATOR	Green 
<b>RPI Part #LGN004</b> LEADWIRE BROWN 30" DUAL/ALLIGATOR	Brown 
<b>RPI Part #LGR004</b> LEADWIRE RED 30" DUAL/ALLIGATOR	Red 
<b>RPI Part #LGW004</b> LEADWIRE WHITE 30" DUAL/ALLIGATOR	White 

Depending on equipment model, these leadwires are also available as sets – see *VETERINARY: DIN/UNSHIELDED LEADWIRE SETS* and *VETERINARY: DUAL/SHIELDED LEADWIRE SETS* on page 11.

## NOTES

[illegible]

# PATIENT, HOLTER & TELEMETRY CABLE OEM CROSS REFERENCE

PATIENT CABLE TYPE Found on	FIXED Page 1	DIN/ UNSHIELDED Page 2	DUAL/ SHIELDED Page 3
<b>PATIENT, HOLTER &amp; TELEMETRY CABLES</b>			
<b>ABBOTT</b> AAEMM, AAXM, ABBV, ABBY SERIES, ABEM, ABTM & EX-50 • Leadwire sets also available, see pages 5-12.		KCA004	KEC010
<b>ABBOTT</b> EKR50 & EXPI • Leadwire sets also available, see pages 5-12.			KEC010
<b>ABBOTT</b>  LIFECARE 5100 & 5200 • Leadwire sets also available, see pages 5-12.	KCB008 <sup>1</sup> KEB012 <sup>2</sup>	KCA007	KCC009
<b>ACCUCORE/PIC/VITALCOM</b> TELEMETRY SYSTEM 3 LEAD ^ KCB028 (Telemetry cable), see page 4.	KCB028^		
<b>ACCUCORE/PIC/VITALCOM</b> TELEMETRY SYSTEM 4 LEAD ^ KDB002 (Telemetry cable), see page 4.	KDB002^		
<b>ACCUCORE/PIC/VITALCOM</b> TELEMETRY SYSTEM 5 LEAD ^ KEB015 (Telemetry cable), see page 4.	KEB015^		
<b>ADVANCED TECH LABS</b> SEE ATL			
<b>AGILENT</b> SEE HEWLETT PACKARD/AGILENT/ZYMED/PHILIPS MEDICAL			
<b>A.I.R</b> 3000 HOLTER (3 Channel 7 LEAD) ^ KGA001 (Holter cable), see page 4. • Leadwire sets also available, see pages 5-12.		KGA001^	
<b>AIR SHIELDS/HILL ROM/VICKER/NARCO</b> HRRM-71, HRRM-72, PORTAFIB I, PORTAFIB II, & SYSTEM 6 NEONATAL • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
<b>AIR SHIELDS/HILL ROM/VICKER/NARCO</b>  PORTAFIB III & SYSTEM 5 NEW • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003 <sup>1</sup> KEC010 <sup>2</sup>
<b>AMERICAN OPTICAL</b> 1125, 1700 & PULSAR 4 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
<b>AMERICAN OPTICAL</b> 5260, 5330, 26125 & 33167 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	
<b>ATL</b> 400; 800; 800 PLUS; VERA MARK; and ULTRAMARK 4, 4 PLUS, 5, 7 & 9 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
<b>BECTON DICKINSON</b> SEE ELECTRODYNE/BECTON DICKINSON/DATAMEDIX/PPG			
<b>BIOCHEM INT'L</b> 3101, 3404, 6100, 9100, ADVISOR & AUTOCORR • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
<b>BIONET</b> CARDIOCARE 2000 & CARDIOTOUGH 3000	KKB002		
<b>BIRTCHEER</b> 300, 335 & 337	KEB014 <sup>2</sup>		
<b>BIRTCHEER</b> 365 CARDIOTRACER, 375 CARDIOTRACER & 375-11	KKB002		
<b>BRAEMER</b> DL125 & DL1700 ^ KGA001 (Holter cable), see page 4. • Leadwire sets also available, see pages 5-12.		KGA001^	
<b>BURDICK/SIEMENS</b> ATRIA 3000; ECLIPSE 4, 8, 400, 800, 850, 850 PLUS; ECLIPSE LE II; PLUS	KKB008		
<b>BURDICK/SIEMENS</b> EK-1, EK-2, EK-3, EK-4, EK-5 & EK-5A	KEB014 <sup>2</sup>		

<sup>1</sup> 3 Leadwire Patient Cable

<sup>2</sup> 5 Leadwire Patient Cable

<sup>3</sup> Old RPI Part #BUC009

PATIENT CABLE TYPE Found on	FIXED Page 1	DIN/ UNSHIELDED Page 2	DUAL/ SHIELDED Page 3
<b>PATIENT, HOLTER &amp; TELEMETRY CABLES (continued)</b>			
<b>BURDICK/SIEMENS</b> E350, E350 I, E550, E560, EK-10, ECLIPSE LE, ECLIPSE LEI, ELITE & ELITE II	KKB003		
<b>BURDICK/SIEMENS</b> MA521 • Leadwire sets also available, see pages 5-12.			KCC034
<b>BURDICK/SIEMENS</b> MEDIC 3, MEDIC 4 & MEDIC 5 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003 <sup>1</sup> KEC010 <sup>2</sup>
<b>BURDICK/SIEMENS</b> SIRECUST 630 & SIRECUST DU • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
<b>CAMBRIDGE</b> VS-1, VS-2 & VS-3	KEB014 <sup>2</sup>		
<b>COLIN MED. INST.</b> BP 308, BP 308E, BP 408 & CMB 3000 • Leadwire sets also available, see pages 5-12.		KCA004	
<b>COROMETRICS</b> 115, 116, 118 & 151 ^ KCD033 (Patient Cable - 3 Lead Multi-Link <sup>®</sup> ), see page 3 • Leadwire sets also available, see pages 5-12.			KCD033^
<b>COROMETRICS</b> 505, 506, 515 & 515A • Leadwire sets also available, see pages 5-12.		KCA013	
<b>COROMETRICS</b> 556 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	
<b>CRITICARE SYS.</b> 1100 & 1508	KEB012		
<b>CRITICARE SYS.</b> 504 US • Leadwire sets also available, see pages 5-12.		KCA004	
<b>CRITICARE SYS.</b> 602-A USP & POET 2 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003 <sup>1</sup> KEC010 <sup>2</sup>
<b>CRITIKON INC.</b> 8700 SERIES, 9700 SERIES, DINAMAP, DINAMAP PLUS, DINAMAP PRO & DINAMAP MPS • Leadwire sets also available, see pages 5-12.	KCB008	KCA007	KCC009
<b>DATAMEDIX</b> SEE ELECTRODYNE/BECTON DICKINSON/DATAMEDIX/PPG			
<b>DATASCOPE</b> 2000 (RESP), 2001 (RESP), 2002 (RESP), 2100 (RESP), 2200 (RESP), EL (RESP) & PASSPORT (RESP) • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
<b>DATASCOPE</b> 2000 (ECG), 2001 (ECG), 2002 (ECG), 2100 (ECG), 2200 (ECG), 3000, BALLOON PUMP, MD-3, MD-3A, MD-4, PASSPORT (ECG), EL (ECG), PASSPORT XG, PASSPORT IL, PASSPORT XL, SPECTRUM, SYS 82, SYS 83, SYS 90, SYS 90T, SYS 97 & TRIO • Leadwire sets also available, see pages 5-12.		KCA016	
<b>DATASCOPE</b> TELEMETRY 3 LEAD ^ KCB028 (Telemetry cable), see page 4.	KCB028^		
<b>DATASCOPE</b> TELEMETRY 4 LEAD ^ KCB002 (Telemetry cable), see page 4.	KDB002^		
<b>DATASCOPE</b> TELEMETRY 5 LEAD ^ KEB015 (Telemetry cable), see page 4.	KEB015^		
<b>DATEx</b> 245P, AESPIRE, AS/3, AS/3 COMPACT, CARDIOCAP I, CARDIOCAP II, CS/3, S/5, S/5 COMPACT, S/5 LIGHT & SATELLITE PLUS • Leadwire sets also available, see pages 5-12.			KCC034 <sup>1</sup> KEC010 <sup>2</sup>

(continued)

# PATIENT, HOLTER & TELEMETRY CABLE OEM CROSS REFERENCE *(continued)*

PATIENT CABLE TYPE Found on	FIXED Page 1	DIN/ UNSHIELDED Page 2	DUAL/ SHIELDED Page 3
<b>PATIENT, HOLTER &amp; TELEMETRY CABLES</b> <i>(continued)</i>			
<b>DMI</b> 3 CH RECORDER ^ KGA001 (Holter cable), see page 4. • Leadwire sets also available, see pages 5-12.		KGA001^	
<b>ELECTRODYNE/BECTON DICKINSON/ DATAMEDIX/PPG</b> MD3A • Leadwire sets also available, see pages 5-12.			KEC010
<b>ELECTRODYNE/BECTON DICKINSON/ DATAMEDIX/PPG</b> GUARDIAN 1001 & GUARDIAN N1001 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003 <sup>①</sup> KEC010 <sup>②</sup>
<b>GE MEDICAL</b> SEE MARQUETTE/GE MEDICAL			
<b>HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL</b> <b>DEFIBRILLATOR (12 PIN):</b> M3500B & M5500B • Leadwire sets also available, see pages 5-12.		KCA010 <sup>①</sup> KEA005 <sup>②</sup>	KCC012 <sup>①</sup> KEC007 <sup>②</sup>
<b>HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL</b> <b>DEFIBRILLATOR (8 PIN):</b> M3500B & M5500B • Leadwire sets also available, see pages 5-12.		KCA022 <sup>①④⑤</sup> KCA025 <sup>①④</sup> KEA017 <sup>②⑥</sup>	
<b>HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL</b> <b>C3 BEDSIDE MONITOR:</b> 862474, 862475, 862478, 862479 & M1205A; <b>CODEMASTER XL &amp; DEFIBRILLATOR (12 PIN):</b> M1722A & M1722B; <b>DEFIBRILLATOR:</b> M3535A & M3536A; <b>ECG MODULE:</b> M1001B, M1002A & M1002B; <b>ECG MONITOR:</b> 78352C (12 PIN), M1001A, M1165A, M1166A & M1167A; <b>INTELL/VUE:</b> M8001A MP20, M8002A MP30, M8003A MP40, M8004A MP50, M8005A MP60, M8007A MP70 & M8010A MP90; <b>MERLIN CMS:</b> M1175A; <b>MMS PATIENT MONITOR:</b> M3001A; <b>NEONATAL MONITOR:</b> 78834C; <b>TELEMETRY MONITOR:</b> M2601A; <b>VERIDIA CMS:</b> M1177A; <b>VERIDIA PATIENT MONITOR:</b> M3000A, M3046A & M4375A • Leadwire sets also available, see pages 5-12.		KCA010 <sup>①</sup> KEA005 <sup>②</sup>	KCC012 <sup>①</sup> KEC007 <sup>②</sup>
<b>HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL</b> <b>CODEMASTER XL &amp; DEFIBRILLATOR (8 PIN):</b> M1723A & M1723B; <b>DEFIBRILLATOR:</b> 43100A (8 PIN), 43110A (8 PIN), 43110M, 43110MC, 43120A (8 PIN), 43200A, 78619A, 78620A, 78660A, 78660B, 78661A, 78661B, 78670A, 78671A; <b>ECG MODULE:</b> 78211A (8 PIN); <b>ULTRASOUND:</b> 77000A, 77020A, 77020AC, 77020AR, 77030A, 77035A, 77040A, M2400A, M2406A, M2475A (8 PIN), M2475B, SONOS 500, SONOS 4500, SONOS 5500, SONOS 7506 & 77025A SONOS • Leadwire sets also available, see pages 5-12.		KCA022 <sup>①④⑤</sup> KCA025 <sup>①④</sup> KEA017 <sup>②⑥</sup>	
<b>HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL</b> <b>CODEMASTER XL &amp; DEFIBRILLATOR (8 PIN):</b> M1722A & M1722B; <b>DEFIBRILLATOR (8 PIN):</b> M3500B + M5500B; <b>ECG MODULE:</b> 78208A & 78210A; <b>ECG MONITOR:</b> 778330A, 78331A, 78332A, 78333A, 78341A, 78342A, 78345A, 78346A, 78351A, 78352A (8 PIN), 78353A, 78353B, 78354A, 78532A, 78532B, 78534A, 78534B & 78534C • Leadwire sets also available, see pages 5-12.		KCA022 <sup>①④⑤</sup> KCA025 <sup>①④</sup> KEA017 <sup>②⑥</sup>	
<b>HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL</b> <b>DEFIBRILLATOR (6 PIN):</b> 43100A, 43110A, 43120A & M3921A; <b>PORTABLE PATIENT MONITOR:</b> M3921A • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003

① 3 Leadwire Patient Cable      ② 5 Leadwire Patient Cable      ③ Old RPI Part #BUC009  
④ Old RPI Part #HPC008      ⑤ Pin colors are White, Green, & Black

PATIENT CABLE TYPE Found on	FIXED Page 1	DIN/ UNSHIELDED Page 2	DUAL/ SHIELDED Page 3
<b>PATIENT, HOLTER &amp; TELEMETRY CABLES</b> <i>(continued)</i>			
<b>HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL</b> <b>ECG MONITOR:</b> 1500A, 1500B, 1504A, 1506A, 1511A & 1511B		KEB030 <sup>③</sup>	
<b>HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL</b> <b>10 LEAD FIXED CABLE:</b> M1770A PAGEWRITER 300Pi, M1771A PAGEWRITER 200 & M1772A PAGEWRITER 100		KKB006	
<b>HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL</b> <b>HOLTER MONITOR:</b> Tri-Track 1100-005 & Tri-Track 1100-008 ^ KGA001 (Holter cable), see page 4.		KGA001^	
<b>HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL</b> <b>78100A SAFETY TELEMETRY</b> ^ KCB029 (Telemetry cable), see page 4.		KCB029 <sup>②</sup> ^	
<b>HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL</b>			KCB031 <sup>①</sup> ^ KCB032 <sup>①</sup> ^ KEB018 <sup>②</sup> ^ KEB019 <sup>②</sup> ^
<b>M2601A TELEMETRY TRANSMITTER</b> ^ KCB031 (Telemetry cable), see page 4. ^ KCB032 (Telemetry cable), see page 4. ^ KEB018 (Telemetry cable), see page 4. ^ KEB019 (Telemetry cable), see page 4.			
<b>HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL</b> <b>5 LEAD TELEMETRY:</b> DIGITRACK 2020 ^ KEB015 (Telemetry cable), see page 4.		KEB015^	
<b>HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL</b> <b>4 LEAD TELEMETRY</b> ^ KCB002 (Telemetry cable), see page 4.		KDB002^	
<b>HILL ROM</b> SEE AIR SHIELDS/HILL ROM/VICKER/NARCO			
<b>INSTROMEDIX</b> KING OF HEARTS ^ KBB001 (Telemetry cable), see page 4.		KBB001^	
<b>INVIVO RESEARCH</b> MILLENNIA, 3500 & 3550 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003 <sup>①</sup> KEC010 <sup>②</sup>
<b>INVIVO RESEARCH</b> OMNI-TRACK 300, M6, M12, NVS & TVS • Leadwire sets also available, see pages 5-12.		KCA004	KEC010
<b>IREX</b> SYSTEM III • Leadwire sets also available, see pages 5-12.		KCA013	
<b>IREX</b> 880 • Leadwire sets also available, see pages 5-12.	KCB008	KCA007	KCC009
<b>IVY BIOMEDICAL</b> 101 (Before 1989), 402, 700 & 3100 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
<b>KERNS &amp; WHITE</b> KEC-56 • Leadwire sets also available, see pages 5-12.		KCA013	
<b>LASER SYSTEMS</b> VISTA1 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
<b>LUMISCOPE</b> 5 LEAD ECG CABLE • Leadwire sets also available, see pages 5-12.		KEA017 <sup>②</sup>	
<b>MARQUETTE/GE MEDICAL</b> 3 CH HOLTER CABLE ^ KGA001 (Holter cable), see page 4. • Leadwire sets also available, see pages 5-12.		KGA001^	

① Old RPI Part #HPC004      ⑤ Old RPI Part #HPC006      ③ Old RPI Part #HPC005  
② Pin colors are White, Red, & Black

*(continued)*



# PATIENT, HOLTER & TELEMETRY CABLE OEM CROSS REFERENCE *(continued)*

PATIENT CABLE TYPE Found on	FIXED Page 1	DIN/ UNSHIELDED Page 2	DUAL/ SHIELDED Page 3
<b>PATIENT, HOLTER &amp; TELEMETRY CABLES</b> <i>(continued)</i>			
<b>MARQUETTE/GE MEDICAL</b> 900 DEFIB SERIES, 1200 DEFIB SERIES, & 1500 DEFIB SERIES • Leadwire sets also available, see pages 5-12.	KCB008 <sup>①</sup> KEB012 <sup>②</sup>	KCA007	KCC009
<b>MARQUETTE/GE MEDICAL</b> 1200 ST, HELLIGE MICROSMART, MAC 400/500/600/800/1000/1100/1200/1600	KKB012		
<b>MARQUETTE/GE MEDICAL</b> APEXPRO TELEMETRY UNITS • Leadwire sets also available, see pages 5-12.			
<b>MARQUETTE/GE MEDICAL</b> CASE, CASE 8000, MAC 5000 & MAC 5000 ST ^ GEC021 (Acquisition Module cable), see page 4. • Leadwire sets also available, see pages 5-12.	GEC021^		
<b>MARQUETTE/GE MEDICAL</b> CASE 12, CASE 15, CASE 16, CENTRA, MAC 6, MAC 8, MAC 12, MAC PC, MAC VU & MAX-PERSONAL ^ GEC020 (Acquisition Module cable), see page 4.	GEC020^		
<b>MARQUETTE/GE MEDICAL</b>  DASH 2000 PRO, DASH 3000 PRO, DASH 4000 PRO, DINAMAP PRO 100, EAGLE SYSTEMS, MAC-LAB, SOLAR 3000, SOLAR 8000, SOLAR 9500, TRAM 100, TRAM 200, TRAM 300, TRAM 400, TRAM 500, TRAM 600, TRAM 800, TRAM AR & UNITY SYSTEMS ^ KED021 (Patient Cable - 5 Lead Multi-Link®), see page 3. • Leadwire sets also available, see pages 5-12.		KCA019 <sup>①</sup> KEA001 <sup>②</sup>	KCC021 <sup>①</sup> KEC003 <sup>②</sup> KED021 <sup>②</sup> ^
<b>MARQUETTE/GE MEDICAL</b> HELLIGE EK-51	KKB007		
<b>MATRIX MEDICAL</b> LIFE DEFENSE PLUS (DEFIB) • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
<b>MDE</b> SEE MEDICAL DATA ELEC.			
<b>MEDICAL DATA ELEC. (MDE)</b> 400673	KCB008		
<b>MEDICAL DATA ELEC. (MDE)</b> ESCORT 100, ESCORT 101, ESCORT 102, ESCORT 200, ESCORT 300 & PRISM (WITHOUT 90° CONNECTOR) • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
<b>MEDICAL DATA ELEC. (MDE)</b> ESCORT W/ HP OPTION • Leadwire sets also available, see pages 5-12.		KCA010 <sup>①</sup> KEA005 <sup>②</sup>	KCC012 <sup>①</sup> KEC007 <sup>②</sup>
<b>MEDICAL DATA ELEC. (MDE)</b> 3 LEAD TELEMETRY ^ KCB028 (Telemetry cable), see page 4.	KCB028^		
<b>MEDICAL DATA ELEC. (MDE)</b> 5 LEAD TELEMETRY ^ KEB015 (Telemetry cable), see page 4.	KEB015^		
<b>MEDICAL RESEARCH LABS</b> 260SL, 450SL, 521A, PRE-AMP 450, PORTAPAK 80, PORTAPAK 90 & PORTAPAK 1000 ISOLATED • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
<b>MEDTRONIC</b> SEE PHYSIO CONTROL/MEDTRONIC			
<b>MENNEN MEDICAL</b> TELEMETRY 4 LEAD ^ KDB002 (Telemetry cable), see page 4.	KDB002^		
<b>MENNEN MEDICAL</b> TELEMETRY 5 LEAD ^ KEB015 (Telemetry cable), see page 4.	KEB015^		
<b>MIDMARK</b> IQ MARK ^ KEB020 (Holter cable), see page 4.	KEB020^		

① 3 Leadwire Patient Cable  
② Old RPI Part #HPC008  
③ 5 Leadwire Patient Cable  
④ Pin colors are White, Green, & Black  
⑤ Old RPI Part #BUC009

PATIENT CABLE TYPE Found on	FIXED Page 1	DIN/ UNSHIELDED Page 2	DUAL/ SHIELDED Page 3
<b>PATIENT, HOLTER &amp; TELEMETRY CABLES</b> <i>(continued)</i>			
<b>MINDRAY</b> BENEVIEW T5, BENEVIEW T6, BENEVIEW T8 & MPM MODULE • Leadwire sets also available, see pages 5-12.		KEA005	KEC007
<b>MORTARA INC.</b> ELI 250	KKB007		
<b>MTM</b> 2 CH HOLTER CABLE ^ KEA004 (Holter cable), see page 5. • Leadwire sets also available, see pages 5-12.		KEA004^	
<b>MTM</b> 3 CH HOLTER CABLE ^ KGA001 (Holter cable), see page 5. • Leadwire sets also available, see pages 5-12.		KGA001^	
<b>NARCO</b> SEE AIR SHIELDS/HILL ROM/VICKER/NARCO			
<b>NELLCOR/PURTAN BENNET</b>  791018, MBP-4000, N-200, N-250 & N-1000 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003 <sup>①</sup> KEC010 <sup>②</sup>
<b>NIHON KOHDEN</b> ECG: 5151A, 5403A, 6151A, 6306A, 6353A, 6543A, 6551A, 7514A, 8110A, 8240A, 8330A, 8350A, 8370A & 9340A	KKB002		
<b>NIHON KOHDEN</b> BSM-8301A, BSM-8302, BSM-8500A, BSM-8502, OEC-5102A (Serial#20156), OEC-6105A (Serial#>22026), OEC-7102A, OEC-8108, TEC-5200A, TEC-7100A, TEC-7200A, TEC-7300A & 8000 SERIES (DEFIB) • Leadwire sets also available, see pages 5-12.		KCA004	KEC010
<b>NORTH AMERICAN DRAEGER</b> VITALERT 1000 & VITALERT 2000 • Leadwire sets also available, see pages 5-12.		KCA004	KEC010
<b>NOVAMETRIX MED</b> LIFE DEFENSE 1 • Leadwire sets also available, see pages 5-12.		KCA004	KEC010
<b>PACE TECH</b> MINIMAX 4000CL, MINIPACK 300, MINIPACK 3000, MINIPACK 3100, VITALMAX 2100, VITALMAX 2200, VITALMAX 4000, VITALMAX 4000CL, VITALMAX 4100 & VITALMAX 800 PLUS • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
<b>PARKS ELECTRONICS</b> 2 CH HOLTER CABLE ^ KEA004 (Holter cable), see page 4. • Leadwire sets also available, see pages 5-12.		KEA004^	
<b>PC SCIENTIFIC</b>  2 CH HOLTER CABLE 3 CH HOLTER CABLE ^ KEA004 (Holter cable), see page 4. ^^ KGA001 (Holter cable), see page 4. • Leadwire sets also available, see pages 5-12.		KEA004^ KGA001^^	
<b>PHILIPS MEDICAL</b> SEE HEWLETT PACKARD/AGILENT/ZYMED/PHILIPS MEDICAL			
<b>PHYSIO CONTROL/MEDTRONIC</b> LIFEPAK 6	KEB012		
<b>PHYSIO CONTROL/MEDTRONIC</b> LIFEPAK 12 & LIFEPAK 20 • Leadwire sets also available, see pages 5-12.	KCB036	KCA035	
<b>PHYSIO CONTROL/MEDTRONIC</b>  LIFEPAK 5, LIFEPAK 7, LIFEPAK 8, LIFEPAK 9, LIFEPAK 10, LIFEPAK 510, VSM-1, VSM-2, VSM-2ESF & VSM-3 W/O RESISTORS • Leadwire sets also available, see pages 5-12.	KCB008 <sup>①</sup> KEB012 <sup>②</sup>	KCA007	KCC009

① Old RPI Part #HPC004  
② Pin colors are White, Red, & Black  
③ Old RPI Part #HPC006  
④ Old RPI Part #HPC005  
*(continued)*

# PATIENT, HOLTER & TELEMETRY CABLE OEM CROSS REFERENCE *(continued)*

PATIENT CABLE TYPE Found on	FIXED Page 1	DIN/ UNSHIELDED Page 2	DUAL/ SHIELDED Page 3
<b>PATIENT, HOLTER &amp; TELEMETRY CABLES</b> <i>(continued)</i>			
<b>PIC</b> SEE ACCUCORE/PIC/VITALCOM			
<b>PPG</b> SEE ELECTRODYNE/BECTON DICKINSON/DATAMEDIX/PPG			
<b>PROTOCOL SYS./WELCH ALLYN</b>  PROPAQ 206 EL (HP OPTION) • Leadwire sets also available, see pages 5-12.			KCC012 <sup>①</sup> KCA010 <sup>②</sup>
<b>PROTOCOL SYS./WELCH ALLYN</b> 6200, ATLAS, ENCORE, PROPAQ 102, PROPAQ 103, PROPAQ 104, PROPAQ 106, PROPAQ CS & ULTRASMART • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
<b>PURTAN BENNET</b> SEE NELLCOR/PURTAN BENNET			
<b>QUINTON</b>  Q-STRESS ^ KKB010 (Patient Cable - 10 Lead Fixed Pinch), see page 1. ^^ KKB011 (Patient Cable - 10 Lead Fixed Snap), see page 1.	KKB010^ KKB011^^		
<b>QUINTON</b> Q-TEL 3LD TELEMETRY ^ KCB028 (Telemetry cable), see page 4.	KCB028^		
<b>QUINTON</b> Q-TEL 4LD TELEMETRY ^ KDB002 (Telemetry cable), see page 4.	KDB002^		
<b>QUINTON</b> 365 ^ KEB015 (Telemetry cable), see page 4.	KEB015^		
<b>ROZINN ELECTRONICS</b> 151 3 CH 7 LD ^ KGA001 (Holter cable), see page 4. • Leadwire sets also available, see pages 5-12.		KGA001^	
<b>ROZINN ELECTRONICS</b> 151 3 CH 5 LD ^ KEA004 (Holter cable), see page 4. • Leadwire sets also available, see pages 5-12.		KEA004^	
<b>SCHILLER</b> 3 CH HOLTER ^ KGA001 (Holter cable), see page 4. • Leadwire sets also available, see pages 5-12.		KGA001^	
<b>SCHILLER</b> AT-1, AT-2, AT-2 PLUS & AT-10	KKB007		
<b>SIEMENS</b> SEE BURDICK/SIEMENS			
<b>SIMONSEN &amp; WELLS</b> 3 LD TELEMETRY ^ KCB028 (Telemetry cable), see page 4.	KCB028^		
<b>SIMONSEN &amp; WELLS</b> 4 LD TELEMETRY ^ KDB002 (Telemetry cable), see page 4.	KDB002^		
<b>SIMONSEN &amp; WELLS</b> 8000 (New Generation), Athena, Diascope 1 & Diascope 2 • Leadwire sets also available, see pages 5-12.			KCC034
<b>SPACE LABS</b> ULTRAVIEW MONITORS (90496) • Leadwire sets also available, see pages 5-12.			KEC016
<b>SPACE LABS</b> 900 RESP W/O RESISTORS	KEB012		
<b>SPACE LABS</b> 400 Series, 500 Series, 900 Series, 701S, 703S, 720S, 4045, 90408, 90418, 90470, 90600 Series, 90700 SERIES, TEK-413, TEK-414, TEK-511, TEK-512, TEK-514 & TEK-521 • Leadwire sets also available, see pages 5-12.		KCA004	KEC010

PATIENT CABLE TYPE Found on	FIXED Page 1	DIN/ UNSHIELDED Page 2	DUAL/ SHIELDED Page 3
<b>PATIENT, HOLTER &amp; TELEMETRY CABLES</b> <i>(continued)</i>			
<b>TEMTECH</b> LIFEGUARD 9, M401, M501 & M503	KEB012		
<b>VICKER</b> SEE AIR SHIELDS/HILL ROM/VICKER/NARCO			
<b>VITALCOM</b> SEE ACCUCORE/PIC/VITALCOM			
<b>WELCH ALLYN</b> SEE PROTOCOL SYS./WELCH ALLYN			
<b>ZOLL</b> D-90, M SERIES, PD-1200, PD-1400, PD-1600 & PD-2000 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	
<b>ZYMED</b> SEE HEWLETT PACKARD/AGILENT/ZYMED/PHILIPS MEDICAL			

① 3 Leadwire Patient Cable  
② Old RPI Part #HPC008

③ 5 Leadwire Patient Cable  
④ Pin colors are White, Green, & Black

⑤ Old RPI Part #BUC009

⑥ Old RPI Part #HPC004  
⑦ Pin colors are White, Red, & Black

⑧ Old RPI Part #HPC006

⑨ Old RPI Part #HPC005

## NOTES

[illegible]

ADDRESS SERVICES REQUESTED

## ***RPI is Your Source for Patient Cables and Leadwires!***



### **PATIENT CABLES AND LEADWIRES TO FIT THESE OEMS AND MORE!**

<b>ABBOTT</b>	<b>GE MEDICAL</b>	<b>NOVAMETRIX MED</b>
<b>ACCUCORE</b>	<b>HEWLETT PACKARD</b>	<b>PAGE TECH</b>
<b>ADVANCED TECH LABS (ATL)</b>	<b>HILL ROM</b>	<b>PARKS ELECTRONICS</b>
<b>AGILENT</b>	<b>INSTROMEDIX</b>	<b>PC SCIENTIFIC</b>
<b>A.I.R.</b>	<b>INVIVO RESEARCH</b>	<b>PHILIPS MEDICAL</b>
<b>AIR SHIELDS</b>	<b>IREX</b>	<b>PHYSIO CONTROL</b>
<b>AMERICAN OPTICAL</b>	<b>IVY BIOMEDICAL</b>	<b>PIC</b>
<b>APPLIED CARDIAC SYSTEMS</b>	<b>KERNS &amp; WHITE</b>	<b>PPG</b>
<b>BECTON DICKINSON</b>	<b>LASER SYSTEMS</b>	<b>PROTOCOL SYS.</b>
<b>BIOCHEM INT'L</b>	<b>LUMISCOPE</b>	<b>PURTAN BENNET</b>
<b>BIONET</b>	<b>MARQUETTE</b>	<b>QUINTON</b>
<b>BIRTCHER</b>	<b>MATRIX MEDICAL</b>	<b>ROZINN ELECTRONICS</b>
<b>BRAEMER</b>	<b>MEDICAL DATA ELEC. (MDE)</b>	<b>SCHILLER</b>
<b>BURDICK</b>	<b>MEDICAL RESEARCH LABS</b>	<b>SIEMENS</b>
<b>COLIN MED. INST.</b>	<b>MEDTRONIC</b>	<b>SIMONSEN &amp; WELLS</b>
<b>COROMETRICS</b>	<b>MENNEN MEDICAL</b>	<b>SPACE LABS</b>
<b>CRITICARE SYS.</b>	<b>MINDRAY</b>	<b>TEMTECH</b>
<b>CRITIKON INC.</b>	<b>MORTARA INC.</b>	<b>VICKER</b>
<b>DATAMEDIX</b>	<b>MTM</b>	<b>VITALCOM</b>
<b>DATASCOPE</b>	<b>NARCO</b>	<b>WELCH ALLYN</b>
<b>DATRIX</b>	<b>NELLCOR</b>	<b>ZOLL</b>
<b>DMI</b>	<b>NIHON KOHDEN</b>	<b>ZYMED and lots more ...</b>
<b>ELECTRODYNE</b>	<b>NORTH AMERICAN DRAEGER</b>	<b>see our website for updates!</b>

*Note: RPI parts may fit or be utilized with equipment not listed above. Please check the equipment manual to ensure an exact fit and function.*

**CALL (800) 221-9723 or (818) 882-8611 • FAX (818) 882-7028 • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)  
E-MAIL [moreinfo@rpiparts.com](mailto:moreinfo@rpiparts.com) • EMAIL [order@rpiparts.com](mailto:order@rpiparts.com) • EMAIL [techhelp@rpiparts.com](mailto:techhelp@rpiparts.com)**