



PATIENT CABLES & LEADWIRES

RPI IS YOUR SOURCE FOR PATIENT CABLES & LEADWIRES!

PATIENT CABLES

- For ECG, EEG, monitoring and defibrillator equipment
- 3, 5, 7 and 10 leadwire sets available
- Snap or Pin style fixed lead patient ends
- Dual/Shielded, Din/Unshielded or Fixed lead patient ends
- Leadwire Adaptor Kit also available to adapt fixed banana plugs or pins to tab or snap electrodes
- Strain Relief protective feature molded into cable to extend life and integrity of the cables

ACQUISITION MODULE CABLES

- Cables to fit GE AM4, AM5 and CAM 14 Acquisition Modules
- 17 ft. durable coiled cords
- Includes molded strain reliefs and connector shields
- With over molded ferrite beads for superior noise reduction

HOLTER CABLES


- 3 Channel feature
- 5 and 7 Din/Unshielded patient ends
- Snap or Pinch style leadwire sets also available
- Strain Relief protective feature molded into cable to extend life and integrity of the cables

TELEMETRY CABLES

- 2, 3, 4 and 5 fixed lead patient ends
- Snap or Pinch style patient ends
- Strain Relief protective feature molded into cable to extend life and integrity of the cables

LEADWIRES

- Sets and Individual leadwires available in standard colors for easy identification
- Snap or Pinch style patient ends
- Dual/Shielded and Din/Unshielded leadwire cable ends
- Reuseable or disposable Leadwire sets available
- Sets and Individual leadwires for Veterinary use also available
- Strain Relief protective feature molded into leadwires to extend wires' life and integrity



FOR MORE INFORMATION ABOUT THE PARTS LISTED IN THIS FLYER, SEE RPI'S WEBSITE: www.rpiparts.com

rpi replacement parts industries, inc.

CALL (800) 221-9723 • E-MAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit various equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. REV D (10/17)

INDEX

Acquisition, Patient, Holter and Telemetry Cables Page

Acquisition Module Cables to fit GE	4
Din/Unshielded Patient Cables	2
Dual/Shielded Patient Cables	3
Dual/Shielded Patient Cables to fit GE Multi-Link®	3
Fixed Patient Cables	1
Holter Cables	4
Telemetry Cables	4

Patient Cable Leadwires and Leadwire Sets

Bonded Leadwire Sets to fit GE Multi-Link®	5
Combined Leadwire sets to fit GE Multi-Link®	5
Din/Unshielded Individual Leadwires	7
Din/Unshielded Leadwire Sets	6-7
Dual/Shielded Individual Leadwires	8
Dual/Shielded Leadwire Sets	8
Leadwire Sets to fit GE Multi-Link®	6

Disposable Patient Cable Leadwire Sets

Din/Unshielded Leadwire Sets	9
Dual/Shielded Leadwire Sets	10

Leadwire Accessories

Leadwire Adaptor Kit	1
Leadwire Combiners	10
Leadwire Organizers	10
Leadwire Shielded Accessory Kits	10

Leadwires and Leadwire Sets for Veterinary Use


Din/Unshielded Individual Leadwires	11
Din/Unshielded Leadwire Sets	11
Dual/Shielded Individual Leadwires	11
Dual/Shielded Leadwire Sets	11

OEM Model Cross Reference

Patient, Holter & Telemetry Cable OEM Cross Reference	13-16
---	-------

FIXED PATIENT CABLES


RPI Part #KCB002



PATIENT CABLE - 3 LEAD FIXED SNAP

- Number of Leads: 3 • Length: 11 ft.


RPI Part #KCB008



PATIENT CABLE - 3 LEAD FIXED SNAP

- Number of Leads: 3 • Length: 11 ft.

RPI Part #KCB036




PATIENT CABLE - 3 LEAD FIXED SNAP

OEM Part #11110-00029/3006218-02

- Number of Leads: 3 • Length: 10 ft.


RPI Part #KEB012



PATIENT CABLE - 5 LEAD FIXED SNAP

- Number of Leads: 5 • Length: 13 ft.


RPI Part #KEB014



PATIENT CABLE - 5 LEAD FIXED PIN

- Number of Leads: 5 • Length: 8 ft.
- Also available: Leadwire Adaptor Kit (RPI PART #LAK001)


RPI Part #KEB030



PATIENT CABLE - 5 LEAD FIXED BANANA

- Number of Leads: 5 • Length: 13 ft.
- Also available: Leadwire Adaptor Kit (RPI PART #LAK001)


RPI Part #KKB002



PATIENT CABLE - 10 LEAD FIXED BANANA

- Number of Leads: 10 • Length: 10 ft.
- Also available: Leadwire Adaptor Kit (RPI PART #LAK001)


RPI Part #KKB003



PATIENT CABLE - 10 LEAD FIXED BANANA

- Number of Leads: 10 • Length: 10 ft.
- Also available: Leadwire Adaptor Kit (RPI PART #LAK001)


RPI Part #KKB006



PATIENT CABLE - 10 LEAD FIXED BANANA

- Number of Leads: 10 • Length: 10 ft.
- Also available: Leadwire Adaptor Kit (RPI PART #LAK001)


RPI Part #KKB007



PATIENT CABLE - 10 LEAD FIXED BANANA

- Number of Leads: 10 • Length: 10 ft.
- Also available: Leadwire Adaptor Kit (RPI PART #LAK001)


RPI Part #KKB008



PATIENT CABLE - 10 LEAD FIXED BANANA

- Number of Leads: 10 • Length: 10 ft.
- Also available: Leadwire Adaptor Kit (RPI PART #LAK001)

RPI Part #KKB010




PATIENT CABLE - 10 LEAD FIXED PINCH

OEM Part #60-00180-01/60-00184-01

- Number of Leads: 10 • Length: 18 ft.

RPI Part #KKB011




PATIENT CABLE - 10 LEAD FIXED SNAP

OEM Part #60-00181-01/60-00185-01

- Number of Leads: 10 • Length: 18 ft.

RPI Part #KKB012



PATIENT CABLE - 10 LEAD FIXED BANANA

OEM Part #223 418 09

- Number of Leads: 10 • Length: 9 ft.
- Also available: Leadwire Adaptor Kit (RPI PART #LAK001)

RPI Part #LAK001

LEADWIRE ADAPTOR KIT


- Includes 10 adaptors individually labeled with leadwire identifiers
- Adapts banana plugs or pins to tab or snap electrodes
- Widest contact area for superior trace quality
- Flat, wide bodied design helps prevent twisting
- Durable construction - built to last



Includes 10 adaptors individually labeled with leadwire identifiers

DUAL/SHIELDED PATIENT CABLES


RPI Part #KCC003



PATIENT CABLE - 3 LEAD DUAL

- Number of Leads: 3 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM014 & LKM017), Disposable Leadwire Sets (RPI Part #'s LDM305 & LDM312), and Veterinary Leadwire Set (RPI Part #LVM705)


RPI Part #KCC009



PATIENT CABLE - 3 LEAD DUAL

- Number of Leads: 3 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM014 & LKM017), Disposable Leadwire Sets (RPI Part #'s LDM305 & LDM312), and Veterinary Leadwire Set (RPI Part #LVM705)


RPI Part #KCC012



PATIENT CABLE - 3 LEAD DUAL

- Number of Leads: 3 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM014 & LKM017), Disposable Leadwire Sets (RPI Part #'s LDM305 & LDM312), and Veterinary Leadwire Set (RPI Part #LVM705)


RPI Part #KCC021



PATIENT CABLE - 3 LEAD DUAL

- Number of Leads: 3 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM014 & LKM017), Disposable Leadwire Sets (RPI Part #'s LDM305 & LDM312), and Veterinary Leadwire Set (RPI Part #LVM705)


RPI Part #KCC034



PATIENT CABLE - 3 LEAD DUAL OEM Part #545302

- Number of Leads: 3 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM014 & LKM017), Disposable Leadwire Sets (RPI Part #'s LDM305 & LDM312), and Veterinary Leadwire Set (RPI Part #LVM705)


RPI Part #KEC003



PATIENT CABLE - 5 LEAD DUAL

- Number of Leads: 5 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM032 & LKM035), Disposable Leadwire Sets (RPI Part #'s LDM307 & LDM314), and Veterinary Leadwire Set (RPI Part #LVM707)


RPI Part #KEC007



PATIENT CABLE - 5 LEAD DUAL

- Number of Leads: 5 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM032 & LKM035), Disposable Leadwire Sets (RPI Part #'s LDM307 & LDM314), and Veterinary Leadwire Set (RPI Part #LVM707)


RPI Part #KEC010



PATIENT CABLE - 5 LEAD DUAL

- Number of Leads: 5 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM032 & LKM035), Disposable Leadwire Sets (RPI Part #'s LDM307 & LDM314), and Veterinary Leadwire Set (RPI Part #LVM707)


RPI Part #KEC016



PATIENT CABLE - 5 LEAD DUAL

- Number of Leads: 5 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM032 & LKM035), Disposable Leadwire Sets (RPI Part #'s LDM307 & LDM314), and Veterinary Leadwire Set (RPI Part #LVM707)

RPI Part #KEC022




PATIENT CABLE - 5 LEAD DUAL OEM Part #545303

- Number of Leads: 5 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM032 & LKM035), Disposable Leadwire Sets (RPI Part #'s LDM307 & LDM314), and Veterinary Leadwire Set (RPI Part #LVM707)

DUAL/SHIELDED PATIENT CABLES TO FIT GE MULTI-LINK®


RPI Part #KCD033



PATIENT CABLE - 3 LEAD MULTI-LINK OEM Part #1554AA0

- Number of Leads: 3 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM052 & LKM053)

RPI Part #KED021



PATIENT CABLE - 5 LEAD MULTI-LINK OEM Part #2017003-001/412931-001

- Number of Leads: 5 • Length: 10 ft.
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM052, LKM053, LKM054, LKM055, LKM056 & LKM057)

ACQUISITION MODULE CABLES TO FIT GE

RPI Part #GEC020

ACQUISITION MODULE CABLE
OEM Part #700044-202/700044-204

- Length: ~17 ft. overall length

Fits: AM4 Acquisition Module and AM5 Acquisition Module



RPI Part #GEC021

ACQUISITION MODULE CABLE
OEM Part #2016560-002/2016560-003/700657-002/700657-003


- Length: ~17 ft. overall length

Fits: CAM 14 Acquisition Module




HOLTER CABLES

RPI Part #KEA004



HOLTER CABLE - 5 LEAD 3 CHANNEL

- Number of Leads: 5 • Length: 30"
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM026, LKM029 & LKM030), Disposable Leadwire Sets (RPI Part #'s LDM303 & LDM310), and Veterinary Leadwire Set (RPI Part #LVM703)




RPI Part #KGA001




HOLTER CABLE - 7 LEAD 3 CHANNEL

- Number of Leads: 7 • Length: 40"
- Suggested Leadwire Sets: Reusable Leadwire Sets (RPI Part #'s LKM038 & LKM041), Disposable Leadwire Sets (RPI Part #'s LDM304 & LDM311), and Veterinary Leadwire Set (RPI Part #LVM704)



RPI Part #KEB020



HOLTER CABLE - 5 LEAD FIXED SNAP
OEM Part #3-370-0001


- Number of Leads: 5 • Length: 43"



TELEMETRY CABLES


Accessories
Also Available
See Page 10

RPI Part #KBB001




TELEMETRY CABLE - 2 LEAD SNAP

- Number of Leads: 2 • Length: 32"




RPI Part #KCB028




TELEMETRY CABLE - 3 LEAD SNAP

- Number of Leads: 3 • Length: 27"
- Also available: 3 Lead Organizer (RPI Part #LAR004)



RPI Part #KCB029




TELEMETRY CABLE - 3 LEAD SNAP

- Number of Leads: 3 • Length: 30" • Includes: Adaptor Block
- Also available: 3 Lead Organizer (RPI Part #LAR004)



TELEMETRY CABLES (continued)

RPI Part #KCB031




TELEMETRY CABLE - 3 LEAD SNAP

- Number of Leads: 3 • Length: 30"
- Also available: 3 Lead Organizer (RPI Part #LAR004)



RPI Part #KCB032




TELEMETRY CABLE - 3 LEAD PINCH

- Number of Leads: 3 • Length: 30"
- Also available: 3 Lead Organizer (RPI Part #LAR004)



RPI Part #KDB002




TELEMETRY CABLE - 4 LEAD SNAP

- Number of Leads: 4 • Length: 27"
- Also available: 5 Lead Organizer (RPI Part #LAR005)



RPI Part #KEB015




TELEMETRY CABLE - 5 LEAD SNAP

- Number of Leads: 5 • Length: 24"
- Also available: 5 Lead Organizer (RPI Part #LAR005)



RPI Part #KEB018




TELEMETRY CABLE - 5 LEAD SNAP

- Number of Leads: 5 • Length: 30"
- Also available: 5 Lead Organizer (RPI Part #LAR005)



RPI Part #KEB019




TELEMETRY CABLE - 5 LEAD PINCH

- Number of Leads: 5 • Length: 30"
- Also available: 5 Lead Organizer (RPI Part #LAR005)



BONDED LEADWIRE SETS TO FIT GE MULTI-LINK®

RPI Part #LKM052




3 LEADWIRE SET - 30" MULTI-LINK/SNAP

OEM Part #411203-001

- Number of Leads: 3 bonded leadwires • Length: 30"
- Includes: 3 Lead Organizer (RPI Part #LAR004)
- Also available: Bonded Leadwire Sets (RPI Part #'s LKM053, LKM054 & LKM055) and Combined Leadwire Sets (RPI Part #'s LKM056 & LKM057)

Fits: GE Multi-Link® Cable, Patient Cable - 3 Lead Multi-Link® (RPI Part #KCD033) and Patient Cable - 5 Lead Multi-Link® (RPI Part #KED021)

RPI Part #LKM053




3 LEADWIRE SET - 30" MULTI-LINK/PINCH

OEM Part #412682-001

- Number of Leads: 3 bonded leadwires • Length: 30"
- Includes: 3 Lead Organizer (RPI Part #LAR004)
- Also available: Bonded Leadwire Sets (RPI Part #'s LKM052, LKM054 & LKM055) and Combined Leadwire Sets (RPI Part #'s LKM056 & LKM057)

Fits: GE Multi-Link® Cable, Patient Cable - 3 Lead Multi-Link® (RPI Part #KCD033) and Patient Cable - 5 Lead Multi-Link® (RPI Part #KED021)

RPI Part #LKM054




5 LEADWIRE SET - 30" MULTI-LINK/SNAP

OEM Part #411202-001

- Number of Leads: 5 bonded leadwires • Length: 30"
- Includes: 5 Lead Organizer (RPI Part #LAR005)
- Also available: Bonded Leadwire Sets (RPI Part #'s LKM052, LKM053 & LKM055) and Combined Leadwire Sets (RPI Part #'s LKM056 & LKM057)

Fits: GE Multi-Link® Cable and Patient Cable - 5 Lead Multi-Link® (RPI Part #KED021)

RPI Part #LKM055



5 LEADWIRE SET - 30" MULTI-LINK/PINCH


OEM Part #412681-001

- Number of Leads: 5 bonded leadwires • Length: 30"
- Includes: 5 Lead Organizer (RPI Part #LAR005)
- Also available: Bonded Leadwire Sets (RPI Part #'s LKM052, LKM053 & LKM054) and Combined Leadwire Sets (RPI Part #'s LKM056 & LKM057)

Fits: GE Multi-Link® Cable and Patient Cable - 5 Lead Multi-Link® (RPI Part #KED021)

COMBINED LEADWIRE SETS TO FIT GE MULTI-LINK®

RPI Part #LKM056




5 LEADWIRE SET - 30" MULTI-LINK/SNAP

OEM Part #411200-001

- Number of Leads: 5 individual leadwires with combiner • Length: 30"
- Includes: 5 Lead Organizer (RPI Part #LAR005) and combiner
- Also available: Bonded Leadwire Sets (RPI Part #'s LKM052, LKM053, LKM054 & LKM055) and Combined Leadwire Set (RPI Part #LKM057)

Fits: GE Multi-Link® Cable and Patient Cable - 5 Lead Multi-Link® (RPI Part #KED021)

RPI Part #LKM057




5 LEADWIRE SET - 30" MULTI-LINK/PINCH

OEM Part #414556-001

- Number of Leads: 5 individual leadwires with combiner • Length: 30"
- Includes: 5 Lead Organizer (RPI Part #LAR005) and combiner
- Also available: Bonded Leadwire Sets (RPI Part #'s LKM052, LKM053, LKM054 & LKM055) and Combined Leadwire Set (RPI Part #LKM056)

Fits: GE Multi-Link® Cable and Patient Cable - 5 Lead Multi-Link® (RPI Part #KED021)

RPI Part #LKM058




6 LEADWIRE SET - 30" MULTI-LINK/SNAP

OEM Part #421930-001

- Number of Leads: 6 individual leadwires with combiner • Length: 30"
- Includes: (2) 3 Lead Organizer (RPI Part #LAR004) and 6 Lead Combiner
- Also available: Combined Leadwire Set (RPI Part #LKM059)

RPI Part #LKM059



6 LEADWIRE SET - 30" MULTI-LINK/PINCH

OEM Part #421932-001

- Number of Leads: 6 individual leadwires with combiner • Length: 30"
- Includes: (2) 3 Lead Organizer (RPI Part #LAR004) and 6 Lead Combiner
- Also available: Combined Leadwire Set (RPI Part #LKM058)

LEADWIRE SETS TO FIT GE MULTI-LINK®


RPI Part #LKM049

10 LEADWIRE SET - CAM 14/BANANA

OEM Part #420101-002

- Number of Leads: 10
- Length: 27": Red (V1), Yellow (V2), Green (V3), Blue (V4), Orange (V5), Violet (V6); and, 37": Green (RL), White (RA), Black (LA) & Red (LL)
- Includes: (2) 5 Lead Organizer (RPI Part #LAR005)
- Also available: 10 Leadwire Set - CAM 14/Pinch (RPI Part #LKM050) & 10 Leadwire Set - CAM 14/Snap (RPI Part #LKM051)

Fits: CAM 14 Module




RPI Part #LKM050

10 LEADWIRE SET - CAM 14/PINCH

OEM Part #420101-002

- Number of Leads: 10
- Length: 27": Red (V1), Yellow (V2), Green (V3), Blue (V4), Orange (V5), Violet (V6); and, 37": Green (RL), White (RA), Black (LA) & Red (LL)
- Includes: (2) 5 Lead Organizer (RPI Part #LAR005)
- Also available: 10 Leadwire Set - CAM 14/Banana (RPI Part #LKM049) & 10 Leadwire Set - CAM 14/Snap (RPI Part #LKM051)

Fits: CAM 14 Module




RPI Part #LKM051

10 LEADWIRE SET - CAM 14/SNAP

OEM Part #420101-002

- Number of Leads: 10
- Length: 27": Red (V1), Yellow (V2), Green (V3), Blue (V4), Orange (V5), Violet (V6); and, 37": Green (RL), White (RA), Black (LA) & Red (LL)
- Includes: (2) 5 Lead Organizer (RPI Part #LAR005)
- Also available: 10 Leadwire Set - CAM 14/Banana (RPI Part #LKM049) & 10 Leadwire Set - CAM 14/Pinch (RPI Part #LKM050)

Fits: CAM 14 Module




DIN/UNSHIELDED LEADWIRE SETS

Accessories
Also Available
See Page 10

RPI Part #LKM002

3 LEADWIRE SET - 24" DIN/PINCH


- Number of Leads: 3 • Length: 24"
- Also available: Disposable Leadwire Set (RPI Part #LDM308) & 3 Lead Organizer (RPI Part #LAR004)



RPI Part #LKM005

3 LEADWIRE SET - 24" DIN/SNAP

- Number of Leads: 3 • Length: 24"
- Also available: Disposable Leadwire Set (RPI Part #LDM301) & 3 Lead Organizer (RPI Part #LAR004)



RPI Part #LKM008

3 LEADWIRE SET - 24" DIN/PINCH

- Number of Leads: 3 • Length: 24"
- Also available: Disposable Leadwire Set (RPI Part #LDM309) & 3 Lead Organizer (RPI Part #LAR004)



RPI Part #LKM011

3 LEADWIRE SET - 24" DIN/SNAP


- Number of Leads: 3 • Length: 24"
- Also available: Disposable Leadwire Set (RPI Part #LDM302) & 3 Lead Organizer (RPI Part #LAR004)



DIN/UNSHIELDED LEADWIRE SETS (continued)

Accessories
Also Available
See Page 10


RPI Part #LKM026



5 LEADWIRE SET - 24" DIN/PINCH

- Number of Leads: 5 • Length: 24"
- Also available: Disposable Leadwire Set (RPI Part #LDM310) & 5 Lead Organizer (RPI Part #LAR005)


RPI Part #LKM029



5 LEADWIRE SET - 24" DIN/SNAP

- Number of Leads: 5 • Length: 24"
- Also available: 24" Disposable Leadwire Set (RPI Part #LDM303) & 5 Lead Organizer (RPI Part #LAR005)


RPI Part #LKM030



5 LEADWIRE SET - 40" DIN/SNAP

- Number of Leads: 5 • Length: 40"
- Also available: 24" Disposable Leadwire Set (RPI Part #LDM303) & 5 Lead Organizer (RPI Part #LAR005)


RPI Part #LKM038



7 LEADWIRE SET - 24" DIN/PINCH

- Number of Leads: 7 • Length: 24"
- Also available: Disposable Leadwire Set (RPI Part #LDM311), 3 Lead Organizer (RPI Part #LAR004) & 5 Lead Organizer (RPI Part #LAR005)


RPI Part #LKM041



7 LEADWIRE SET - 24" DIN/SNAP


- Number of Leads: 7 • Length: 24"
- Also available: Disposable Leadwire Set (RPI Part #LDM304), 3 Lead Organizer (RPI Part #LAR004) & 5 Lead Organizer (RPI Part #LAR005)

DIN/UNSHIELDED INDIVIDUAL LEADWIRES



- Equipment End: Din/Unshielded, .060" Safety Socket
- Patient End: Pinch
- Length: 24"

DESCRIPTION	LEADWIRE COLOR
RPI Part #LPB002 LEADWIRE BLACK 24" DIN/PINCH	Black
RPI Part #LPE002 LEADWIRE ORANGE 24" DIN/PINCH	Orange
RPI Part #LPG002 LEADWIRE GREEN 24" DIN/PINCH	Green
RPI Part #LPN002 LEADWIRE BROWN 24" DIN/PINCH	Brown
RPI Part #LPR002 LEADWIRE RED 24" DIN/PINCH	Red
RPI Part #LPU002 LEADWIRE BLUE 24" DIN/PINCH	Blue
RPI Part #LPW002 LEADWIRE WHITE 24" DIN/PINCH	White



- Equipment End: Din/Unshielded, .060" Safety Socket
- Patient End: Snap
- Length: 24"


DESCRIPTION	LEADWIRE COLOR
RPI Part #LSB002 LEADWIRE BLACK 24" DIN/SNAP	Black
RPI Part #LSE002 LEADWIRE ORANGE 24" DIN/SNAP	Orange
RPI Part #LSG002 LEADWIRE GREEN 24" DIN/SNAP	Green
RPI Part #LSN002 LEADWIRE BROWN 24" DIN/SNAP	Brown
RPI Part #LSR002 LEADWIRE RED 24" DIN/SNAP	Red
RPI Part #LSU002 LEADWIRE BLUE 24" DIN/SNAP	Blue
RPI Part #LSW002 LEADWIRE WHITE 24" DIN/SNAP	White

Depending on equipment model, these leadwires are also available as sets – see *DIN/UNSHIELDED LEADWIRE SETS* on pages 6-7.

DUAL/SHIELDED LEADWIRE SETS

Accessories
Also Available
See Page 10


RPI Part #LKM014



3 LEADWIRE SET - 40" DUAL/PINCH

- Number of Leads: 3 • Length: 40"
- Includes: 3 Lead Shielded Accessory Kit (RPI Part #LAK006)
- Also available: Disposable Leadwire Set (RPI Part #LDM312)


RPI Part #LKM017



3 LEADWIRE SET - 40" DUAL/SNAP

- Number of Leads: 3 • Length: 40"
- Includes: 3 Lead Shielded Accessory Kit (RPI Part #LAK006)
- Also available: Disposable Leadwire Set (RPI Part #LDM305)


RPI Part #LKM020



3 LEADWIRE SET - 40" DUAL/PINCH

- Number of Leads: 3 • Length: 40"
- Includes: 3 Lead Shielded Accessory Kit (RPI Part #LAK006)
- Also available: Disposable Leadwire Set (RPI Part #LDM313)


RPI Part #LKM023



3 LEADWIRE SET - 40" DUAL/SNAP

- Number of Leads: 3 • Length: 40"
- Includes: 3 Lead Shielded Accessory Kit (RPI Part #LAK006)
- Also available: Disposable Leadwire Set (RPI Part #LDM306)


RPI Part #LKM032



5 LEADWIRE SET - 40" DUAL/PINCH

- Number of Leads: 5 • Length: 40"
- Includes: 5 Lead Shielded Accessory Kit (RPI Part #LAK007)
- Also available: Disposable Leadwire Set (RPI Part #LDM314)


RPI Part #LKM035



5 LEADWIRE SET - 40" DUAL/SNAP


- Number of Leads: 5 • Length: 40"
- Includes: 5 Lead Shielded Accessory Kit (RPI Part #LAK007)
- Also available: Disposable Leadwire Set (RPI Part #LDM307)

DUAL/SHIELDED INDIVIDUAL LEADWIRES



- Equipment End: Dual/Shielded
- Patient End: Pinch
- Length: 40"

DESCRIPTION	LEADWIRE COLOR
RPI Part #LPB005 LEADWIRE BLACK 40" DUAL/PINCH	Black
RPI Part #LPG005 LEADWIRE GREEN 40" DUAL/PINCH	Green
RPI Part #LPN005 LEADWIRE BROWN 40" DUAL/PINCH	Brown
RPI Part #LPR005 LEADWIRE RED 40" DUAL/PINCH	Red
RPI Part #LPW005 LEADWIRE WHITE 40" DUAL/PINCH	White



- Equipment End: Dual/Shielded
- Patient End: Snap
- Length: 40"


DESCRIPTION	LEADWIRE COLOR
RPI Part #LSB005 LEADWIRE BLACK 40" DUAL/SNAP	Black
RPI Part #LSG005 LEADWIRE GREEN 40" DUAL/SNAP	Green
RPI Part #LSN005 LEADWIRE BROWN 40" DUAL/SNAP	Brown
RPI Part #LSR005 LEADWIRE RED 40" DUAL/SNAP	Red
RPI Part #LSW005 LEADWIRE WHITE 40" DUAL/SNAP	White

Depending on equipment model, these leadwires are also available as sets – see *DUAL/SHIELDED LEADWIRE SETS* at the top of this page.

DISPOSABLE: DIN/UNSHIELDED LEADWIRE SETS

Accessories
Also Available
See Page 10


RPI Part #LDM301



3 LEADWIRE SET DISPOSABLE 24" DIN/SNAP

- Number of Leads: 3 • Length: 24" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM005) & 3 Lead Organizer (RPI Part #LAR004)


RPI Part #LDM302



3 LEADWIRE SET DISPOSABLE 24" DIN/SNAP

- Number of Leads: 3 • Length: 24" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM011) & 3 Lead Organizer (RPI Part #LAR004)


RPI Part #LDM303



5 LEADWIRE SET DISPOSABLE 24" DIN/SNAP

- Number of Leads: 5 • Length: 24" • 10 sets per package
- Also available: Reusable 24" Leadwire Set (RPI Part #LKM029), Reusable 40" Leadwire Set (RPI Part #LKM030) & 5 Lead Organizer (RPI Part #LAR005)


RPI Part #LDM304



7 LEADWIRE SET DISPOSABLE 24" DIN/SNAP

- Number of Leads: 7 • Length: 24" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM041)


RPI Part #LDM308



3 LEADWIRE SET DISPOSABLE 24" DIN/PINCH

- Number of Leads: 3 • Length: 24" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM002) & 3 Lead Organizer (RPI Part #LAR004)


RPI Part #LDM309



3 LEADWIRE SET DISPOSABLE 24" DIN/PINCH

- Number of Leads: 3 • Length: 24" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM008) & 3 Lead Organizer (RPI Part #LAR004)


RPI Part #LDM310



5 LEADWIRE SET DISPOSABLE 24" DIN/PINCH

- Number of Leads: 5 • Length: 24" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM026) & 5 Lead Organizer (RPI Part #LAR005)

RPI Part #LDM311




7 LEADWIRE SET DISPOSABLE 24" DIN/PINCH

- Number of Leads: 7 • Length: 24" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM038)

DISPOSABLE: DUAL/SHIELDED LEADWIRE SETS

Accessories
Also Available
See Below


RPI Part #LDM305



3 LEADWIRE SET DISPOSABLE 30" DUAL/SNAP

- Number of Leads: 3 • Length: 30" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM017) & 3 Lead Shielded Accessory Kit (RPI Part #LAK006)


RPI Part #LDM306



3 LEADWIRE SET DISPOSABLE 30" DUAL/SNAP

- Number of Leads: 3 • Length: 30" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM023) & 3 Lead Shielded Accessory Kit (RPI Part #LAK006)


RPI Part #LDM307



5 LEADWIRE SET DISPOSABLE 30" DUAL/SNAP

- Number of Leads: 5 • Length: 30" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM035) & 5 Lead Shielded Accessory Kit (RPI Part #LAK007)


RPI Part #LDM312



3 LEADWIRE SET DISPOSABLE 30" DUAL/PINCH

- Number of Leads: 3 • Length: 30" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM014) & 3 Lead Shielded Accessory Kit (RPI Part #LAK006)


RPI Part #LDM313



3 LEADWIRE SET DISPOSABLE 30" DUAL/PINCH

- Number of Leads: 3 • Length: 30" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM020) & 3 Lead Shielded Accessory Kit (RPI Part #LAK006)

RPI Part #LDM314



5 LEADWIRE SET DISPOSABLE 30" DUAL/PINCH

- Number of Leads: 5 • Length: 30" • 10 sets per package
- Also available: Reusable Leadwire Set (RPI Part #LKM032) & 5 Lead Shielded Accessory Kit (RPI Part #LAK007)

ACCESSORIES


DESCRIPTION	ACCESSORY
RPI Part #LAC002 3 LEAD COMBINER <ul style="list-style-type: none"> • Number of Leads: 3 • Use with Dual/Shielded Leadwires only Also included in: <ul style="list-style-type: none"> • 3 Lead Shielded Accessory Kit (RPI Part #LAK006) 	
RPI Part #LAR004 3 LEAD ORGANIZER <ul style="list-style-type: none"> • Number of Leads: 3 Also included in: <ul style="list-style-type: none"> • 3 Lead Shielded Accessory Kit (RPI Part #LAK006) 	
RPI Part #LAK006 3 LEAD SHIELDED ACCESSORY KIT <ul style="list-style-type: none"> • Number of Leads: 3 • Shielded Accessory Kit • Use with Dual/Shielded Leadwires only • Includes all parts as shown 	

DESCRIPTION	ACCESSORY
RPI Part #LAC003 5 LEAD COMBINER <ul style="list-style-type: none"> • Number of Leads: 5 • Use with Dual/Shielded Leadwires only Also included in: <ul style="list-style-type: none"> • 5 Lead Shielded Accessory Kit (RPI Part #LAK007) 	
RPI Part #LAR005 5 LEAD ORGANIZER <ul style="list-style-type: none"> • Number of Leads: 5 Also included in: <ul style="list-style-type: none"> • 5 Lead Shielded Accessory Kit (RPI Part #LAK007) 	
RPI Part #LAK007 5 LEAD SHIELDED ACCESSORY KIT <ul style="list-style-type: none"> • Number of Leads: 5 • Shielded Accessory Kit • Use with Dual/Shielded Leadwires only • Includes all parts as shown 	

VETERINARY: DIN/UNSHIELDED LEADWIRE SETS

Accessories
Also Available
See Page 10


RPI Part #LVM701



3 LEADWIRE SET - 24" DIN/ALLIGATOR

- Reuseable
- Number of Leads: 3 • Length: 24"
- Recommended for veterinary use only
- Also available: 3 Lead Organizer (RPI Part #LAR004)


RPI Part #LVM703



5 LEADWIRE SET - 24" DIN/ALLIGATOR

- Reuseable
- Number of Leads: 5 • Length: 24"
- Recommended for veterinary use only
- Also available: 5 Lead Organizer (RPI Part #LAR005)


RPI Part #LVM702



3 LEADWIRE SET - 24" DIN/ALLIGATOR

- Reuseable
- Number of Leads: 3 • Length: 24"
- Recommended for veterinary use only
- Also available: 3 Lead Organizer (RPI Part #LAR004)

RPI Part #LVM704




7 LEADWIRE SET - 24" DIN/ALLIGATOR

- Reuseable
- Number of Leads: 7 • Length: 24"
- Recommended for veterinary use only
- Also available: 3 Lead Organizer (RPI Part #LAR004) & 5 Lead Organizer (RPI Part #LAR005)

VETERINARY: DUAL/SHIELDED LEADWIRE SETS

Accessories
Also Available
See Page 10

RPI Part #LVM705



3 LEADWIRE SET - 30" DUAL/ALLIGATOR

- Reuseable
- Number of Leads: 3 • Length: 30"
- Recommended for veterinary use only
- Also available: 3 Lead Shielded Accessory Kit (RPI Part #LAK006)


RPI Part #LVM706



3 LEADWIRE SET - 30" DUAL/ALLIGATOR

- Reuseable
- Number of Leads: 3 • Length: 30"
- Recommended for veterinary use only
- Also available: 3 Lead Shielded Accessory Kit (RPI Part #LAK006)

RPI Part #LVM707



5 LEADWIRE SET - 30" DUAL/ALLIGATOR








- Reuseable
- Number of Leads: 5 • Length: 30"
- Recommended for veterinary use only
- Also available: 5 Lead Shielded Accessory Kit (RPI Part #LAK007)

VETERINARY: INDIVIDUAL LEADWIRES


DIN/UNSHIELDED INDIVIDUAL








- Equipment End: Din/Unshielded, .060 Safety Socket
- Patient End: Alligator Clip
- Length: 24"

DESCRIPTION	LEADWIRE COLOR
RPI Part #LGB002 LEADWIRE BLACK 24" DIN/ALLIGATOR	Black 
RPI Part #LGE002 LEADWIRE ORANGE 24" DIN/ALLIGATOR	Orange 
RPI Part #LGG002 LEADWIRE GREEN 24" DIN/ALLIGATOR	Green 
RPI Part #LGN002 LEADWIRE BROWN 24" DIN/ALLIGATOR	Brown 
RPI Part #LGR002 LEADWIRE RED 24" DIN/ALLIGATOR	Red 
RPI Part #LGU002 LEADWIRE BLUE 24" DIN/ALLIGATOR	Blue 
RPI Part #LGW002 LEADWIRE WHITE 24" DIN/ALLIGATOR	White 

DUAL/SHIELDED INDIVIDUAL



- Equipment End: Dual/Shielded
- Patient End: Alligator Clip
- Length: 30"

DESCRIPTION	LEADWIRE COLOR
RPI Part #LGB004 LEADWIRE BLACK 30" DUAL/ALLIGATOR	Black 
RPI Part #LGG004 LEADWIRE GREEN 30" DUAL/ALLIGATOR	Green 
RPI Part #LGN004 LEADWIRE BROWN 30" DUAL/ALLIGATOR	Brown 
RPI Part #LGR004 LEADWIRE RED 30" DUAL/ALLIGATOR	Red 
RPI Part #LGW004 LEADWIRE WHITE 30" DUAL/ALLIGATOR	White 

Depending on equipment model, these leadwires are also available as sets – see *VETERINARY: DIN/UNSHIELDED LEADWIRE SETS* and *VETERINARY: DUAL/SHIELDED LEADWIRE SETS* on page 11.

NOTES

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

PATIENT, HOLTER & TELEMETRY CABLE OEM CROSS REFERENCE

PATIENT CABLE TYPE Found on	FIXED Page 1	DIN/ UNSHIELDED Page 2	DUAL/ SHIELDED Page 3
PATIENT, HOLTER & TELEMETRY CABLES			
ABBOTT AAEMM, AAXM, ABBV, ABBY SERIES, ABEM, ABTM & EX-50 • Leadwire sets also available, see pages 5-12.		KCA004	KEC010
ABBOTT EKR50 & EXPI • Leadwire sets also available, see pages 5-12.			KEC010
ABBOTT LIFECARE 5100 & 5200 • Leadwire sets also available, see pages 5-12.	KCB008 ^① KEB012 ^②	KCA007	KCC009
ACCUCORE/PIC/VITALCOM TELEMETRY SYSTEM 3 LEAD ^KCB028 (Telemetry cable), see page 4.	KCB028^		
ACCUCORE/PIC/VITALCOM TELEMETRY SYSTEM 4 LEAD ^KDB002 (Telemetry cable), see page 4.	KDB002^		
ACCUCORE/PIC/VITALCOM TELEMETRY SYSTEM 5 LEAD ^KEB015 (Telemetry cable), see page 4.	KEB015^		
ADVANCED TECH LABS SEE ATL			
AGILENT SEE HEWLETT PACKARD/AGILENT/ZYMED/PHILIPS MEDICAL			
A.I.R 3000 HOLTER (3 Channel 7 LEAD) ^ KGA001 (Holter cable), see page 4. • Leadwire sets also available, see pages 5-12.		KGA001^	
AIR SHIELDS/HILL ROM/VICKER/NARCO HRRM-71, HRRM-72, PORTAFIB I, PORTAFIB II, & SYSTEM 6 NEONATAL • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
AIR SHIELDS/HILL ROM/VICKER/NARCO PORTAFIB III & SYSTEM 5 NEW • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003 ^① KEC010 ^②
AMERICAN OPTICAL 1125, 1700 & PULSAR 4 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
AMERICAN OPTICAL 5260, 5330, 26125 & 33167 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	
ATL 400; 800; 800 PLUS; VERA MARK; and ULTRAMARK 4, 4 PLUS, 5, 7 & 9 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
BECTON DICKINSON SEE ELECTRODYNE/BECTON DICKINSON/DATAMEDIX/PPG			
BIOCHEM INT'L 3101, 3404, 6100, 9100, ADVISOR & AUTOCORR • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
BIONET CARDIOCARE 2000 & CARDIOTOUCH 3000	KKB002		
BIRTCHEER 300, 335 & 337	KEB014 ^②		
BIRTCHEER 365 CARDIOTRACER, 375 CARDIOTRACER & 375-11	KKB002		
BRAEMER DL125 & DL1700 ^ KGA001 (Holter cable), see page 4. • Leadwire sets also available, see pages 5-12.		KGA001^	
BURDICK/SIEMENS ATRIA 3000; ECLIPSE 4, 8, 400, 800, 850, 850 PLUS; ECLIPSE LE II; PLUS	KKB008		
BURDICK/SIEMENS EK-1, EK-2, EK-3, EK-4, EK-5 & EK-5A	KEB014 ^②		

① 3 LEADWIRE PATIENT CABLE

② 5 LEADWIRE PATIENT CABLE

③ Old RPI Part #BUC009

PATIENT CABLE TYPE Found on	FIXED Page 1	DIN/ UNSHIELDED Page 2	DUAL/ SHIELDED Page 3
PATIENT, HOLTER & TELEMETRY CABLES (continued)			
BURDICK/SIEMENS E350, E350 I, E550, E560, EK-10, ECLIPSE LE, ECLIPSE LEI, ELITE & ELITE II	KKB003		
BURDICK/SIEMENS MA521 • Leadwire sets also available, see pages 5-12.			KCC034
BURDICK/SIEMENS MEDIC 3, MEDIC 4 & MEDIC 5 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003 ^① KEC010 ^②
BURDICK/SIEMENS SIRECUST 630 & SIRECUST DU • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
CAMBRIDGE VS-1, VS-2 & VS-3	KEB014 ^②		
COLIN MED. INST. BP 308, BP 308E, BP 408 & CMB 3000 • Leadwire sets also available, see pages 5-12.		KCA004	
COROMETRICS 115, 116, 118 & 151 ^ KCD033 (Patient Cable - 3 Lead Multi-Link [®]), see page 3 • Leadwire sets also available, see pages 5-12.			KCD033^
COROMETRICS 505, 506, 515 & 515A • Leadwire sets also available, see pages 5-12.		KCA013	
COROMETRICS 556 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	
CRITICARE SYS. 1100 & 1508	KEB012		
CRITICARE SYS. 504 US • Leadwire sets also available, see pages 5-12.		KCA004	
CRITICARE SYS. 602-A USP & POET 2 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003 ^① KEC010 ^②
CRITIKON INC. 8700 SERIES, 9700 SERIES, DINAMAP, DINAMAP PLUS, DINAMAP PRO & DINAMAP MPS • Leadwire sets also available, see pages 5-12.	KCB008	KCA007	KCC009
DATAMEDIX SEE ELECTRODYNE/BECTON DICKINSON/DATAMEDIX/PPG			
DATASCOPE 2000 (RESP), 2001 (RESP), 2002 (RESP), 2100 (RESP), 2200 (RESP), EL (RESP) & PASSPORT (RESP) • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
DATASCOPE 2000 (ECG), 2001 (ECG), 2002 (ECG), 2100 (ECG), 2200 (ECG), 3000, BALLOON PUMP, MD-3, MD-3A, MD-4, PASSPORT (ECG), EL (ECG), PASSPORT XG, PASSPORT IL, PASSPORT XL, SPECTRUM, SYS 82, SYS 83, SYS 90, SYS 90T, SYS 97 & TRIO • Leadwire sets also available, see pages 5-12.		KCA016	
DATASCOPE TELEMETRY 3 LEAD ^ KCB028 (Telemetry cable), see page 4.	KCB028^		
DATASCOPE TELEMETRY 4 LEAD ^ KCB002 (Telemetry cable), see page 4.	KDB002^		
DATASCOPE TELEMETRY 5 LEAD ^ KEB015 (Telemetry cable), see page 4.	KEB015^		
DATEx 245P, AESPIRE, AS/3, AS/3 COMPACT, CARDIOCAP I, CARDIOCAP II, CS/3, S/5, S/5 COMPACT, S/5 LIGHT & SATELLITE PLUS • Leadwire sets also available, see pages 5-12.			KCC034 ^① KEC010 ^②

(continued)

PATIENT, HOLTER & TELEMETRY CABLE OEM CROSS REFERENCE (continued)

PATIENT CABLE TYPE Found on	FIXED Page 1	DIN/ UNSHIELDED Page 2	DUAL/ SHIELDED Page 3
PATIENT, HOLTER & TELEMETRY CABLES (continued)			
DMI 3 CH RECORDER ^ KGA001 (Holter cable), see page 4. • Leadwire sets also available, see pages 5-12.		KGA001^	
ELECTRODYNE/BECTON DICKINSON/ DATAMEDIX/PPG MD3A • Leadwire sets also available, see pages 5-12.			KEC010
ELECTRODYNE/BECTON DICKINSON/ DATAMEDIX/PPG GUARDIAN 1001 & GUARDIAN N1001 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003 ^① KEC010 ^②
GE MEDICAL SEE MARQUETTE/GE MEDICAL			
HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL DEFIBRILLATOR (12 PIN): M3500B & M5500B • Leadwire sets also available, see pages 5-12.		KCA010 ^① KEA005 ^②	KCC012 ^① KEC007 ^②
HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL DEFIBRILLATOR (8 PIN): M3500B & M5500B • Leadwire sets also available, see pages 5-12.		KCA022 ^{①④⑤} KCA025 ^{①⑤} KEA017 ^{②⑥}	
HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL C3 BEDSIDE MONITOR: 862474, 862475, 862478, 862479 & M1205A; CODEMASTER XL & DEFIBRILLATOR (12 PIN): M1722A & M1722B; DEFIBRILLATOR: M3535A & M3536A; ECG MODULE: M1001B, M1002A & M1002B; ECG MONITOR: 78352C (12 PIN), M1001A, M1165A, M1166A & M1167A; INTELL/VUE: M8001A MP20, M8002A MP30, M8003A MP40, M8004A MP50, M8005A MP60, M8007A MP70 & M8010A MP90; MERLIN CMS: M1175A; MMS PATIENT MONITOR: M3001A; NEONATAL MONITOR: 78834C; TELEMETRY MONITOR: M2601A; VERIDIA CMS: M1177A; VERIDIA PATIENT MONITOR: M3000A, M3046A & M4375A • Leadwire sets also available, see pages 5-12.		KCA010 ^① KEA005 ^②	KCC012 ^① KEC007 ^②
HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL CODEMASTER XL & DEFIBRILLATOR (8 PIN): M1723A & M1723B; DEFIBRILLATOR: 43100A (8 PIN), 43110A (8 PIN), 43110M, 43110MC, 43120A (8 PIN), 43200A, 78619A, 78620A, 78660A, 78660B, 78661A, 78661B, 78670A, 78671A; ECG MODULE: 78211A (8 PIN); ULTRASOUND: 77000A, 77020A, 77020AC, 77020AR, 77030A, 77035A, 77040A, M2400A, M2406A, M2475A (8 PIN), M2475B, SONOS 500, SONOS 4500, SONOS 5500, SONOS 7506 & 77025A SONOS • Leadwire sets also available, see pages 5-12.		KCA022 ^{①④⑤} KCA025 ^{①⑤} KEA017 ^{②⑥}	
HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL CODEMASTER XL & DEFIBRILLATOR (8 PIN): M1722A & M1722B; DEFIBRILLATOR (8 PIN): M3500B + M5500B; ECG MODULE: 78208A & 78210A; ECG MONITOR: 778330A, 78331A, 78332A, 78333A, 78341A, 78342A, 78345A, 78346A, 78351A, 78352A (8 PIN), 78353A, 78353B, 78354A, 78532A, 78532B, 78534A, 78534B & 78534C • Leadwire sets also available, see pages 5-12.		KCA022 ^{①④⑤} KCA025 ^{①⑤} KEA017 ^{②⑥}	
HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL DEFIBRILLATOR (6 PIN): 43100A, 43110A, 43120A & M3921A; PORTABLE PATIENT MONITOR: M3921A • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003

① 3 LEADWIRE PATIENT CABLE ② 5 LEADWIRE PATIENT CABLE ③ Old RPI Part #BUC009
④ Old RPI Part #HPC008 ⑤ PIN COLORS ARE WHITE, GREEN, & BLACK

PATIENT CABLE TYPE Found on	FIXED Page 1	DIN/ UNSHIELDED Page 2	DUAL/ SHIELDED Page 3
PATIENT, HOLTER & TELEMETRY CABLES (continued)			
HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL ECG MONITOR: 1500A, 1500B, 1504A, 1506A, 1511A & 1511B		KEB030 ^③	
HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL 10 LEAD FIXED CABLE: M1770A PAGEWRITER 300Pi, M1771A PAGEWRITER 200 & M1772A PAGEWRITER 100		KKB006	
HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL HOLTER MONITOR: Tri-Track 1100-005 & Tri-Track 1100-008 ^ KGA001 (Holter cable), see page 4.		KGA001^	
HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL 78100A SAFETY TELEMETRY ^ KCB029 (Telemetry cable), see page 4.		KCB029 ^② ^	
HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL M2601A TELEMETRY TRANSMITTER ^ KCB031 (Telemetry cable), see page 4. ^ KCB032 (Telemetry cable), see page 4. ^ KEB018 (Telemetry cable), see page 4. ^ KEB019 (Telemetry cable), see page 4.			KCB031 ^① ^ KCB032 ^① ^ KEB018 ^② ^ KEB019 ^② ^
HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL 5 LEAD TELEMETRY: DIGITRACK 2020 ^ KEB015 (Telemetry cable), see page 4.		KEB015^	
HEWLETT PACKARD/AGILENT/ZYMED/ PHILIPS MEDICAL 4 LEAD TELEMETRY ^ KCB002 (Telemetry cable), see page 4.		KDB002^	
HILL ROM SEE AIR SHIELDS/HILL ROM/VICKER/NARCO			
INSTROMEDIX KING OF HEARTS ^ KBB001 (Telemetry cable), see page 4.		KBB001^	
INVIVO RESEARCH MILLENNIA, 3500 & 3550 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003 ^① KEC010 ^②
INVIVO RESEARCH OMNI-TRACK 300, M6, M12, NVS & TVS • Leadwire sets also available, see pages 5-12.		KCA004	KEC010
IREX SYSTEM III • Leadwire sets also available, see pages 5-12.		KCA013	
IREX 880 • Leadwire sets also available, see pages 5-12.	KCB008	KCA007	KCC009
IVY BIOMEDICAL 101 (Before 1989), 402, 700 & 3100 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
KERNS & WHITE KEC-56 • Leadwire sets also available, see pages 5-12.		KCA013	
LASER SYSTEMS VISTA1 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
LUMISCOPE 5 LEAD ECG CABLE • Leadwire sets also available, see pages 5-12.		KEA017 ^②	
MARQUETTE/GE MEDICAL 3 CH HOLTER CABLE ^ KGA001 (Holter cable), see page 4. • Leadwire sets also available, see pages 5-12.		KGA001^	

① Old RPI Part #HPC004 ② Old RPI Part #HPC006 ③ Old RPI Part #HPC005
④ PIN COLORS ARE WHITE, RED, & BLACK (continued)

PATIENT, HOLTER & TELEMETRY CABLE OEM CROSS REFERENCE (continued)

PATIENT CABLE TYPE Found on	FIXED Page 1	DIN/ UNSHIELDED Page 2	DUAL/ SHIELDED Page 3
PATIENT, HOLTER & TELEMETRY CABLES (continued)			
MARQUETTE/GE MEDICAL 900 DEFIB SERIES, 1200 DEFIB SERIES, & 1500 DEFIB SERIES • Leadwire sets also available, see pages 5-12.	KCB008 ^① KEB012 ^②	KCA007	KCC009
MARQUETTE/GE MEDICAL 1200 ST, HELLOG MICROSMART, MAC 400/500/600/800/1000/1100/1200/1600	KKB012		
MARQUETTE/GE MEDICAL APEXPRO TELEMETRY UNITS • Leadwire sets also available, see pages 5-12.			
MARQUETTE/GE MEDICAL CASE, CASE 8000, MAC 5000 & MAC 5000 ST ^ GEC021 (Acquisition Module cable), see page 4. • Leadwire sets also available, see pages 5-12.	GEC021^		
MARQUETTE/GE MEDICAL CASE 12, CASE 15, CASE 16, CENTRA, MAC 6, MAC 8, MAC 12, MAC PC, MAC VU & MAX-PERSONAL ^ GEC020 (Acquisition Module cable), see page 4.	GEC020^		
MARQUETTE/GE MEDICAL DASH 2000 PRO, DASH 3000 PRO, DASH 4000 PRO, DINAMAP PRO 100, EAGLE SYSTEMS, MAC-LAB, SOLAR 3000, SOLAR 8000, SOLAR 9500, TRAM 100, TRAM 200, TRAM 300, TRAM 400, TRAM 500, TRAM 600, TRAM 800, TRAM AR & UNITY SYSTEMS ^ KED021 (Patient Cable - 5 Lead Multi-Link®), see page 3. • Leadwire sets also available, see pages 5-12.		KCA019 ^① KEA001 ^②	KCC021 ^① KEC003 ^② KED021 ^③
MARQUETTE/GE MEDICAL HELLIGE EK-51	KKB007		
MATRIX MEDICAL LIFE DEFENSE PLUS (DEFIB) • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
MDE SEE MEDICAL DATA ELEC.			
MEDICAL DATA ELEC. (MDE) 400673	KCB008		
MEDICAL DATA ELEC. (MDE) ESCORT 100, ESCORT 101, ESCORT 102, ESCORT 200, ESCORT 300 & PRISM (WITHOUT 90° CONNECTOR) • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
MEDICAL DATA ELEC. (MDE) ESCORT W/ HP OPTION • Leadwire sets also available, see pages 5-12.		KCA010 ^① KEA005 ^②	KCC012 ^① KEC007 ^②
MEDICAL DATA ELEC. (MDE) 3 LEAD TELEMETRY ^ KCB028 (Telemetry cable), see page 4.	KCB028^		
MEDICAL DATA ELEC. (MDE) 5 LEAD TELEMETRY ^ KEB015 (Telemetry cable), see page 4.	KEB015^		
MEDICAL RESEARCH LABS 260SL, 450SL, 521A, PRE-AMP 450, PORTAPAK 80, PORTAPAK 90 & PORTAPAK 1000 ISOLATED • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
MEDTRONIC SEE PHYSIO CONTROL/MEDTRONIC			
MENNEN MEDICAL TELEMETRY 4 LEAD ^ KDB002 (Telemetry cable), see page 4.	KDB002^		
MENNEN MEDICAL TELEMETRY 5 LEAD ^ KEB015 (Telemetry cable), see page 4.	KEB015^		
MIDMARK IQ MARK ^ KEB020 (Holter cable), see page 4.	KEB020^		

① 3 LEADWIRE PATIENT CABLE
② 5 LEADWIRE PATIENT CABLE
③ Old RPI Part #BUC009
④ Old RPI Part #HPC008
⑤ PIN COLORS ARE WHITE, GREEN, & BLACK

PATIENT CABLE TYPE Found on	FIXED Page 1	DIN/ UNSHIELDED Page 2	DUAL/ SHIELDED Page 3
PATIENT, HOLTER & TELEMETRY CABLES (continued)			
MINDRAY BENEVIEW T5, BENEVIEW T6, BENEVIEW T8 & MPM MODULE • Leadwire sets also available, see pages 5-12.		KEA005	KEC007
MORTARA INC. ELI 250	KKB007		
MTM 2 CH HOLTER CABLE ^ KEA004 (Holter cable), see page 5. • Leadwire sets also available, see pages 5-12.		KEA004^	
MTM 3 CH HOLTER CABLE ^ KGA001 (Holter cable), see page 5. • Leadwire sets also available, see pages 5-12.		KGA001^	
NARCO SEE AIR SHIELDS/HILL ROM/VICKER/NARCO			
NELLCOR/PURTAN BENNET 791018, MBP-4000, N-200, N-250 & N-1000 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003 ^① KEC010 ^②
NIHON KOHDEN ECG: 5151A, 5403A, 6151A, 6306A, 6353A, 6543A, 6551A, 7514A, 8110A, 8240A, 8330A, 8350A, 8370A & 9340A	KKB002		
NIHON KOHDEN BSM-8301A, BSM-8302, BSM-8500A, BSM-8502, OEC-5102A (Serial#20156), OEC-6105A (Serial#>22026), OEC-7102A, OEC-8108, TEC-5200A, TEC-7100A, TEC-7200A, TEC-7300A & 8000 SERIES (DEFIB) • Leadwire sets also available, see pages 5-12.		KCA004	KEC010
NORTH AMERICAN DRAEGER VITALERT 1000 & VITALERT 2000 • Leadwire sets also available, see pages 5-12.		KCA004	KEC010
NOVAMETRIX MED LIFE DEFENSE 1 • Leadwire sets also available, see pages 5-12.		KCA004	KEC010
PAGE TECH MINIMAX 4000CL, MINIPACK 300, MINIPACK 3000, MINIPACK 3100, VITALMAX 2100, VITALMAX 2200, VITALMAX 4000, VITALMAX 4000CL, VITALMAX 4100 & VITALMAX 800 PLUS • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
PARKS ELECTRONICS 2 CH HOLTER CABLE ^ KEA004 (Holter cable), see page 4. • Leadwire sets also available, see pages 5-12.		KEA004^	
PC SCIENTIFIC 2 CH HOLTER CABLE 3 CH HOLTER CABLE ^ KEA004 (Holter cable), see page 4. ^^ KGA001 (Holter cable), see page 4. • Leadwire sets also available, see pages 5-12.		KEA004^ KGA001^^	
PHILIPS MEDICAL SEE HEWLETT PACKARD/AGILENT/ZYMED/PHILIPS MEDICAL			
PHYSIO CONTROL/MEDTRONIC LIFEPAK 6	KEB012		
PHYSIO CONTROL/MEDTRONIC LIFEPAK 12 & LIFEPAK 20 • Leadwire sets also available, see pages 5-12.	KCB036	KCA035	
PHYSIO CONTROL/MEDTRONIC LIFEPAK 5, LIFEPAK 7, LIFEPAK 8, LIFEPAK 9, LIFEPAK 10, LIFEPAK 510, VSM-1, VSM-2, VSM-2ESF & VSM-3 W/O RESISTORS • Leadwire sets also available, see pages 5-12.	KCB008 ^① KEB012 ^②	KCA007	KCC009

① Old RPI Part #HPC004
② Old RPI Part #HPC006
③ PIN COLORS ARE WHITE, RED, & BLACK
④ Old RPI Part #HPC005

(continued)

PATIENT, HOLTER & TELEMETRY CABLE OEM CROSS REFERENCE *(continued)*

PATIENT CABLE TYPE Found on	FIXED Page 1	DIN/ UNSHIELDED Page 2	DUAL/ SHIELDED Page 3
PATIENT, HOLTER & TELEMETRY CABLES <i>(continued)</i>			
PIC SEE ACCUCORE/PIC/VITALCOM			
PPG SEE ELECTRODYNE/BECTON DICKINSON/DATAMEDIX/PPG			
PROTOCOL SYS./WELCH ALLYN PROPAQ 206 EL (HP OPTION) • Leadwire sets also available, see pages 5-12.			KCC012 ^① KCA010 ^②
PROTOCOL SYS./WELCH ALLYN 6200, ATLAS, ENCORE, PROPAQ 102, PROPAQ 103, PROPAQ 104, PROPAQ 106, PROPAQ CS & ULTRASMART • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	KCC003
PURTAN BENNET SEE NELLCOR/PURTAN BENNET			
QUINTON Q-STRESS ^ KKB010 (Patient Cable - 10 Lead Fixed Pinch), see page 1. ^^ KKB011 (Patient Cable - 10 Lead Fixed Snap), see page 1.	KKB010^ KKB011^^		
QUINTON Q-TEL 3LD TELEMETRY ^ KCB028 (Telemetry cable), see page 4.	KCB028^		
QUINTON Q-TEL 4LD TELEMETRY ^ KDB002 (Telemetry cable), see page 4.	KDB002^		
QUINTON 365 ^ KEB015 (Telemetry cable), see page 4.	KEB015^		
ROZINN ELECTRONICS 151 3 CH 7 LD ^ KGA001 (Holter cable), see page 4. • Leadwire sets also available, see pages 5-12.		KGA001^	
ROZINN ELECTRONICS 151 3 CH 5 LD ^ KEA004 (Holter cable), see page 4. • Leadwire sets also available, see pages 5-12.		KEA004^	
SCHILLER 3 CH HOLTER ^ KGA001 (Holter cable), see page 4. • Leadwire sets also available, see pages 5-12.		KGA001^	
SCHILLER AT-1, AT-2, AT-2 PLUS & AT-10	KKB007		
SIEMENS SEE BURDICK/SIEMENS			
SIMONSEN & WELLS 3 LD TELEMETRY ^ KCB028 (Telemetry cable), see page 4.	KCB028^		
SIMONSEN & WELLS 4 LD TELEMETRY ^ KDB002 (Telemetry cable), see page 4.	KDB002^		
SIMONSEN & WELLS 8000 (New Generation), Athena, Diascope 1 & Diascope 2 • Leadwire sets also available, see pages 5-12.			KCC034
SPACE LABS ULTRAVIEW MONITORS (90496) • Leadwire sets also available, see pages 5-12.			KEC016
SPACE LABS 900 RESP W/O RESISTORS	KEB012		
SPACE LABS 400 Series, 500 Series, 900 Series, 701S, 703S, 720S, 4045, 90408, 90418, 90470, 90600 Series, 90700 SERIES, TEK-413, TEK-414, TEK-511, TEK-512, TEK-514 & TEK-521 • Leadwire sets also available, see pages 5-12.		KCA004	KEC010

PATIENT CABLE TYPE Found on	FIXED Page 1	DIN/ UNSHIELDED Page 2	DUAL/ SHIELDED Page 3
PATIENT, HOLTER & TELEMETRY CABLES <i>(continued)</i>			
TEMTECH LIFEGUARD 9, M401, M501 & M503	KEB012		
VICKER SEE AIR SHIELDS/HILL ROM/VICKER/NARCO			
VITALCOM SEE ACCUCORE/PIC/VITALCOM			
WELCH ALLYN SEE PROTOCOL SYS./WELCH ALLYN			
ZOLL D-90, M SERIES, PD-1200, PD-1400, PD-1600 & PD-2000 • Leadwire sets also available, see pages 5-12.	KCB002	KCA001	
ZYMED SEE HEWLETT PACKARD/AGILENT/ZYMED/PHILIPS MEDICAL			

① 3 LEADWIRE PATIENT CABLE
② Old RPI Part #HPC008

③ 5 LEADWIRE PATIENT CABLE
④ PIN COLORS ARE WHITE, GREEN, & BLACK
⑤ Old RPI Part #BUC009

⑥ Old RPI Part #HPC004

⑦ Old RPI Part #HPC006
⑧ PIN COLORS ARE WHITE, RED, & BLACK
⑨ Old RPI Part #HPC005

NOTES

[illegible]

RPI is Your Source for Patient Cables and Leadwires!



PATIENT CABLES AND LEADWIRES TO FIT THESE OEMS AND MORE!

ABBOTT
ACCUCORE
ADVANCED TECH LABS (ATL)
AGILENT
A.I.R.
AIR SHIELDS
AMERICAN OPTICAL
APPLIED CARDIAC SYSTEMS
BECTON DICKINSON
BIOCHEM INT'L
BIONET
BIRTCHER
BRAEMER
BURDICK
COLIN MED. INST.
COROMETRICS
CRITICARE SYS.
CRITIKON INC.
DATAMEDIX
DATASCOPE
DATRIX
DMI
ELECTRODYNE
GE MEDICAL

HEWLETT PACKARD
HILL ROM
INSTROMEDIX
INVIVO RESEARCH
IREX
IVY BIOMEDICAL
KERNS & WHITE
LASER SYSTEMS
LUMISCOPE
MARQUETTE
MATRIX MEDICAL
MEDICAL DATA ELEC. (MDE)
MEDICAL RESEARCH LABS
MEDTRONIC
MENNEN MEDICAL
MIDMARK
MINDRAY
MORTARA INC.
MTM
NARCO
NELLCOR
NIHON KOHDEN
NORTH AMERICAN DRAEGER
NOVAMETRIX MED

PACE TECH
PARKS ELECTRONICS
PC SCIENTIFIC
PHILIPS MEDICAL
PHYSIO CONTROL
PIC
PPG
PROTOCOL SYS.
PURTAN BENNET
QUINTON
ROZINN ELECTRONICS
SCHILLER
SIEMENS
SIMONSEN & WELLS
SPACE LABS
TEMTECH
VICKER
VITALCOM
WELCH ALLYN
ZOLL
ZYMED

*and lots more ...
see our website for updates!*

Note: RPI parts may fit or be utilized with equipment not listed above. Please check the equipment manual to ensure an exact fit and function.

**CALL (800) 221-9723 • WEBSITE www.rpiparts.com
E-MAIL moreinfo@rpiparts.com • E-MAIL order@rpiparts.com • E-MAIL techhelp@rpiparts.com**