





RPI Part #CSP042
OEM Part #61301602498

SPARE PARTS PACKAGE
• Includes all parts as shown
Fits: Machine in general
Models: 3522


RPI Part #CSP042 Spare Parts Package to fit Castle/Getinge® Bulk Sterilizers



RPI Part #AMK174
OEM Part #56161
SOLENOID VALVE REPAIR KIT




(2x)
RPI Part #AMK177
OEM Part #61301601739
SOLENOID VALVE REPAIR KIT




RPI Part #AMK183
OEM Part #54190
SOLENOID VALVE REPAIR KIT




RPI Part #CSG002
OEM Part #61301605017
DOOR GASKET (16" x 16")




(2x)
RPI Part #CSK048
OEM Part #68164
SOLENOID VALVE REPAIR KIT




(3x)
RPI Part #CSK066
OEM Part #513369
CHECK VALVE REBUILD KIT



RPI Part #CSR028
OEM Part #86407
SOLID STATE RELAY



RPI Part #CSR029
OEM Part #86408
SOLID STATE RELAY



(2x)
RPI Part #CST059
OEM Part #61301602492
STEAM TRAP