


Suggested Cable Ties for use with RPI Hoses to fit Midmark•Ritter® Tables and Chairs!

Service Tip: Don't Forget the Cable Ties when Installing New Hoses!

RPI Part #RPT085

CABLE TIE (8" BLACK)


- 25 per package
- Without mounting hole
- 3-16 x 8" lg.
- Color: Black



RPI Part #RPT279

CABLE TIE (8" BLACK)

- 25 per package
- With # 10 mounting hole
- 3-16 x 8" lg.
- Color: Black



RPI Part #RPT375

CABLE TIE (12" BLACK)

- 25 per package
- Without mounting hole
- 3-16 x 12" lg.
- Color: Black

